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INTRODUCTION

*Please read this manually carefully prior to using your Compex stimulator.
You are strongly advised to read chapter I "Warnings" of this manual.*

CONTENTS

I. WARNINGS

1. Counter-indications	148
2. Safety measures	148

II. PRESENTATION

1. Reception of equipment and accessories	151
2. Guarantee	151
3. Maintenance	151
4. Storage and transportation conditions	152
5. Conditions of use	152
6. Elimination	152
7. Standards	152
8. Patents	152
9. Normalised symbols	153
10. Technical characteristics	153

III. HOW DOES ELECTROSTIMULATION WORK? 154

IV. USAGE GUIDELINES

1. Electrode positions	156
2. Stimulation positions	156
3. Adjusting stimulation energies	157
4. Progression in the levels	157

V. THE **m2** TECHNOLOGY

1. Practical rules of use	157
---------------------------	-----

VI. HOW IT WORKS

1. Description of the stimulator	158
2. Connections	159
3. Preliminary settings	159
4. Selecting a motivation	160
5. Selecting a programme	160
6. Personalising a programme	160
7. During the stimulation session	161
8. Consumption and recharging	163
9. Problems and solutions	164

VII. PROGRAMMES AND SPECIFIC APPLICATIONS

1. Your motivation: Sport	166
2. Your motivation: Health	170
3. Your motivation: Aesthetic	173

VIII. EMC TABLE (ELECTROMAGNETIC COMPATIBILITY)

175

147

I WARNINGS

1. Counter-indications

Important counter-indications

- Cardiac stimulator (pacemaker)
- Epilepsy
- Pregnancy (do not use on abdominal region)
- Serious arterial circulation problems in lower limbs
- Abdominal or inguinal hernia

Precautions when using the Compex

- After trauma or a recent operation (less than 6 months)
- Muscular atrophy
- Persistent pains
- Need for muscular rehabilitation

Osteosynthesis equipment

The presence of osteosynthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc.) is not a counter-indication for the use of Compex programmes. The electrical currents of the Compex are specially designed to have no harmful effect on osteosynthesis equipment.

2. Safety measures

What you should not do with the Compex and the *mi* system

- Do not use the Compex or the *mi*-SENSOR system in water or in a humid atmosphere (sauna, hydrotherapy, etc.).

Important!

- Do not use Compex stimulator programmes if you have sensitivity problems.
- Never use the Compex for prolonged periods without medical advice.
- Consult your doctor if you are in any doubt whatsoever.
- Read this manual carefully, in particular chapter VII, which provides information concerning the effects and indications of each stimulation programme.

person who is sitting or lying down.

In rare instances, people of a nervous disposition may experience a vasovagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, all that is required is to stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).

- Never allow movement resulting from muscular contraction during a stimulation session. You should always stimulate isometrically; this means that the extremities of the limb in which a muscle is being stimulated must be firmly fixed, so as to prevent the movement that results from contraction.
- Do not use the Compex or **mi-SENSOR** system if you are connected to a high-frequency surgical instrument as this could cause skin irritation or burns under the electrodes.
- Do not use the Compex or **mi-SENSOR** system within X metre (see EMC table) of short wave or microwave devices as this could alter the currents generated by the stimulator. If you are in any doubt as to the use of the Compex close to another medical device, seek advice from the manufacturer of the latter or from your doctor.
- Do not use the Compex stimulator or **mi-SENSOR** system in areas where unprotected devices are used to emit electromagnetic radiation. Portable communications equipment can interfere with electrical medical equipment.
- Only use stimulation cables supplied by Compex.
- Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first.
- Never connect stimulation cables to an external power supply. There is a risk of electric shock.
- Only use rechargeable battery units supplied by Compex.
- Never recharge the stimulator without first disconnecting the stimulation cables.
- Always use the charger provided by Compex to recharge the batteries.
- Never use the Compex or the charger if it is damaged (case, cables, etc.) or if the battery compartment is open. There is a risk of electric shock.
- Disconnect the charger immediately if the Compex "bleeps" continuously, if there is abnormal heating or smell, or if smoke comes from the charger or the Compex.
- Do not recharge the battery in a confined space (carrying case, etc.). There is a risk of fire or electric shock.
- Keep the Compex and its accessories out of reach of children.
- Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the Compex, the battery compartment and the charger.
- Sudden temperature changes can cause condensation to build up inside the stimulator. To prevent this, allow it to reach ambient temperature before use.
- Do not use the Compex while driving or operating machinery.
- Do not use the stimulator at altitudes of over 3,000 metres.

Where never to apply the electrodes

- In the vicinity of the head
- Counter-laterally, i.e. do not use two poles connected to the same channel on opposite sides of the body.
- On or near skin lesions of any kind (wounds, swelling, burns, irritation, eczema, etc.).

Precautions when using the

Mi-SENSOR system

- To access the Compex's **Mi** technology functions, it is essential to connect the **Mi-SENSOR** cable before switching on the device.
- Avoid connecting the **Mi-SENSOR** cable when the Compex is switched on.
- Do not disconnect the **Mi-SENSOR** cable while the sensor is being used.
- To work correctly, the **Mi-SENSOR** system must not be blocked or subjected to pressure.

Precautions when using the electrodes

- Only use electrodes supplied by Compex. Other electrodes may have electrical properties that are unsuitable for the Compex stimulator.
- Always turn off the stimulator before moving or removing any electrodes during a session.
- Do not place the electrodes in water.
- Do not apply solvents of any kind to the electrodes.
- For best results, wash and clean the skin of any oil and dry it before attaching the electrodes.
- Attach the electrodes in such a way that their entire surface is in contact with the skin.
- For obvious reasons of hygiene, each user must have his/her own electrode set. Do not use the same electrodes on different people.

- Never use a set of adhesive electrodes for more than 15 sessions as their bonding power deteriorates over time and optimal contact is very important for both user comfort and efficacy.
- Some people with very sensitive skin may experience redness under the electrodes after a session. Generally, this redness is completely harmless and disappears after 10 to 20 minutes. Never start another stimulation session in the same area, however, if the redness is still visible.

II PRESENTATION

1. Reception of equipment and accessories

Your kit contains:

514000 1 stimulator

68302x 1 charger

601131 1 set of black Snap connectors for electrodes, colour-coded (blue, green, yellow, red)

602076 2 bags of small electrodes (5x5 cm)

602077 2 bags of large electrodes (5x10 cm)

88572x 2 user and application manuals

88562x 2 quick start guides "Test your Compex in 5 minutes"

949000 1 belt clip

680029 1 travel pouch

2. Guarantee

Refer to the flyer attached.

3. Maintenance

Do not sterilise the stimulator.

Clean using a soft cloth and an alcohol-based, solvent-free cleaning product.

Use only a minimum amount of liquid when cleaning the Compex.

Do not attempt to repair the stimulator or any of its accessories. Never dismantle the Compex or the charger containing high-voltage parts because of risk of electric shock.

Compex Médical SA declines all responsibility for any damages or consequences resulting from unauthorised attempts to open, modify, or repair the stimulator. This may only be done by persons or repair services authorised by Compex Médical SA.

Your stimulator does not require calibration. Each Compex stimulator is always tested and validated prior to distribution. Its characteristics do not vary under normal conditions.

Nonetheless, as the Compex is a high-quality electrical instrument, its lifespan depends on the use that is made of it and the care and maintenance it receives during its lifetime. If your stimulator contains parts that seem worn or defective, please contact your nearest Compex Médical SA customer service centre regarding an upgrade.

Medical and health professionals are obliged to service the unit in compliance with relevant national laws and regulations. This involves testing performance and safety parameters at regular intervals.

4. Storage and transportation conditions

The Compex contains rechargeable batteries. This is why it must be stored and transported in accordance with the following conditions:

Storage and transportation temperature:
from -20°C to +45°C

Maximum relative humidity: 75%

Atmospheric pressure: 700 hPa to 1060 hPa

5. Conditions of use

Temperature: 0°C to 40°C

Relative humidity: 30% to 75%

Atmospheric pressure: 700 hPa to 1060 hPa

Do not use in areas where there is a risk of explosion.

6. Elimination

The main purpose of the 2002/96/EEC Directive is to prevent the creation of waste electrical and electronic equipment (WEEE) and to reduce the amount of waste to be disposed of by encouraging reuse, recycling, and other forms of recovery. The wheellie bin pictogram with the bar through it means that the equipment cannot be discarded with household refuse, but that it must be collected selectively. The equipment must be delivered to a suitable collection point for treatment. By doing so, you will be contributing to the safeguarding of natural resources and health.

Batteries must be disposed of in compliance with relevant national regulatory requirements.

7. Standards

The Compex complies with current medical standards.

To guarantee your safety, the Compex has been designed, manufactured, and distributed in compliance with the requirements of European Directive 93/42/EC on medical devices.

The Compex also complies with the IEC 60601-1 standard on general safety requirements for electro-medical devices, the IEC 60601-1-2 standard on electromagnetic compatibility, and the IEC 60601-2-10 standard on particular safety requirements for nerve and muscle stimulators.

Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation). **T**he Compex also complies with Directive 2002/96/EEC on waste electrical and electronic equipment (WEEE).

8. Patents

The Compex incorporates several innovations with patents pending.

9. Normalised symbols

Caution: Read User Manual or operating instructions (symbol no. 0434 IEC 60878).

The Complex is a class II device with internal electric power and type BF applied parts (symbol no. 5333 IEC 60878).

Waste electrical and electronic equipment (WEEE) marking according to EN 50419 standard.

The On/Off button is a multi-function button (symbol no. 5009 IEC 60878).

Functions

On/Off (two stable positions)

Waiting or on standby for a part of the unit

Stop (turns system off)

10. Technical characteristics

General information

941210 Battery unit

Nickel metal hybrid (NiMH) rechargeable (4.8V / ≥ 1200 mA/h).

68302x Chargers

Only battery chargers bearing the following information can be used to recharge the batteries supplied with Complex stimulators:

Europe

683020

Type TR503-02-E-133A03

Input 90-264V AC/47-63 Hz/0.5A max.

Output 9V/400 mA/6W

UK

683022

Type TR503-02-U-133A03

Input 90-264V AC/47-63 Hz/0.5A max.

Output 9V/400 mA/6W

601131 Black Snap-connection stimulation cables

Device connector: 6-pin

Electrode connector: female Snap

Length: 1500 mm

601160 Stimulation cable fitted with *mi*-SENSOR system (accessory sold separately)

Device connector: 6-pin

Electrode connector: female Snap

Length: 1500 mm

Protection index

IPXO (IEC 60529)

Neuro-stimulation

All electrical specifications are given for an impedance of 500-1000 ohms per channel.

Channels: Four independent and individually adjustable channels that are electrically isolated from each other and earthed.

Pulse shape: Constant rectangular current with pulse compensation to eliminate any direct current component to prevent residual polarisation at skin level.

Maximum pulse intensity: 100 mA.

Pulse intensity increments: manual adjustment of stimulation intensity from 0 to 694 (energy) in minimum increments of 0.5 mA.

Pulse width: 60 to 400 μ s.

Maximum electrical charge per pulse: 80 microcoulombs (2x40 μ C, compensated).

Standard pulse ramp-up time: 3 μ s (20%-80% of maximum current).

III HOW DOES ELECTROSTIMULATION WORK?

The principle of electrostimulation is to stimulate nerve fibres by means of electrical impulses transmitted by electrodes. The electrical pulses generated by Compex stimulators are high quality pulses - offering safety, comfort and efficiency - which can stimulate different types of nerve fibres:


1. the motor nerves, to stimulate a muscular response. The quantity and the benefits obtained depend on the stimulation parameters and this is known as electro-muscular stimulation (EMS).
2. certain types of sensitive nerve fibres to obtain analgesic or pain-relieving effects.

1. Motor nerve stimulation (EMS)

In voluntary activity, the order for muscular work comes from the brain, which sends a command to the nerve fibres in the form of an electrical signal. This signal is then transmitted to the muscular fibres, which contract. The principle of electrostimulation accurately reproduces the process observed during a voluntary contraction. The stimulator sends an electrical current impulse to the

nerve fibres, exciting them.

his excitation is then transmitted to the muscular fibres causing a basic mechanical response (= muscular twitch). The latter constitutes the basic requirement for muscular contraction. This muscular response is completely identical to muscular work controlled by the brain. In other words, the muscle cannot distinguish whether the command comes from the brain or from the stimulator. The parameters of the Compex programmes (number of impulses per second, contraction time, rest time, total programme time) subject the muscles to different types of work, according to muscular fibres. In fact, different types of muscular fibres may be distinguished according to their respective contraction speed: slow, intermediate and fast fibres. Fast fibres will obviously predominate in a sprinter, while a marathon runner will have more slow fibres. With a good knowledge of human physiology and a perfect mastery of the stimulation parameters of the various programmes, muscular work can be directed very precisely towards the desired goal (muscular reinforcement, increased blood flow, firming up, etc.).

2. Stimulation of the sensitive nerves

The electrical impulses can also excite the sensitive nerve fibres to obtain an analgesic or pain-relieving effect. The stimulation of the tactile sensitive nerve fibres blocks the transmission of pain by the nervous system. The stimulation of another type of sensitive fibres creates an increase in the production of endorphins and, therefore, a reduction of pain. With pain relief programmes, electrostimulation can be used to treat localised sharp or chronic pains as well as muscular pains. Caution: do not use the pain relief programmes for a long period without medical advice.

Benefits of electrostimulation

Electro-stimulation is a very effective way to make your muscles work:

- with significant improvement of different muscular qualities
- without cardio-vascular or mental fatigue
- with limited stress on the joints and tendons. Electrostimulation thus allows a greater quantity of work by the muscles compared with voluntary activity.

To be effective, this work must involve the greatest possible number of muscular fibres. The number of fibres working depends on the stimulation energy. The maximum tolerable energy should therefore be

used. The user controls this aspect of stimulation. The higher the stimulation energy, the greater the number of muscular fibres that are working and, therefore, the more significant the progress achieved. To maximise results, Compex recommends that you complement your electrostimulation sessions with other efforts, such as:

- regular exercise
- proper and healthy nutrition
- a balanced lifestyle

IV USAGE GUIDELINES

The use guidelines presented in this section should be considered as general rules. For all programmes, it is recommended that you read carefully the usage information and advice presented in chapter VII "Programmes and specific applications".

1. Electrode positions

For optimal results, use the electrode positions recommended by Compex. To do this, refer to the pictures and pictograms shown on the inside cover of the manual.

Each stimulation cable has two poles:

A positive pole (+) = red connection

A negative pole (-) = black connection

A different electrode must be connected to each pole.

Note: It is possible and normal to have an electrode arrangement that leaves one electrode connection free from a cable.

Depending on the characteristics of the current, efficacy can be optimised in certain programmes by placing the electrode connected to the positive pole (red connection) "strategically".

When working with a muscle stimulation programme (programme involving muscle contractions), it is important to place the positive electrode on the motor point of the muscle.

It is crucial to choose the right size electrodes (large or small) and correctly position these on the muscle group you want to stimulate to ensure the efficacy of the treatment. Therefore, always use the size of electrodes shown in the pictures. Unless you have other specific medical instructions, always follow the placement directions in the pictures.

Where necessary, look for the best possible position by slowly moving the positive electrode over the muscle until you find the point that will produce the best contraction or the most comfort for you.

Compex disclaims all responsibility for consequences arising from electrodes placed in other positions.

2. Stimulation positions

This position will vary depending on the position of the electrodes, the muscle group you wish to stimulate, and the programme you are using.

For programmes involving powerful muscular contractions, the muscle should always be stimulated in an isometric fashion. You must therefore fix the extremities of your limbs securely.

In this way, you provide maximum resistance to the movement and prevent any shortening of the muscle during the contraction, which could create cramp pains and serious stiffness after the session. For example, when stimulating the quadriceps, the user should be in a seated position with the ankles fixed with straps to prevent extension of the knees.

For other types of programmes (for example, the **Health** and **Active recovery** programmes), which do not involve powerful muscular contractions, position yourself as comfortably as possible.

To determine the stimulation position to be adopted depending on the chosen placement of the electrodes and programme, refer to chapter VII "Programmes and specific applications".

3. Adjusting stimulation energies

In a stimulated muscle, the number of recruited fibres depends on the stimulation energy. For programmes involving powerful muscular contractions, you must therefore use maximum stimulation energies (up to 694), always at the limit that you can endure, in order to recruit the maximum number of fibres.

4. Progression in the levels

In general, it is not advisable to go through the different levels quickly with the intention of reaching level 3 as fast as possible. In fact, the different levels correspond to progress with electrostimulation.

The simplest and most usual procedure is to start with level 1 and raise the level when changing to a new stimulation cycle.

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session

V THE **mi** TECHNOLOGY

mi stands for *muscle intelligence*TM (all items referring to this technology are preceded by the symbol **mi**).

This technology takes into account the specific features of each of our muscles and thus offers stimulation adapted to their characteristics.

It's simple, as the data is transmitted to the stimulator automatically!

It's personalised, as each of our muscles is unique!

The Energy stimulator bears the label **mi-READY, as it can use certain functions of the **mi** technology if a stimulation cable fitted with the **mi-SENSOR** system (sold separately) is connected.**

1. Practical rules of use

To access the **mi** functions, make sure that the stimulation cable fitted with the **mi-SENSOR** system (sold separately) is connected before you turn on the stimulator.

Avoid connecting the stimulation cable fitted with the **mi-SENSOR** system when the *Complex* is switched on.

To work correctly, the **mi-SENSOR** system must not be blocked or subjected to pressure.

During the stimulation session, the stimulation cable fitted with the **mi-SENSOR** system must always be connected to an electrode.

Mi-SENSOR (sold separately)

- This is a little sensor that links the stimulator to the electrodes. The **Mi-SENSOR** is the key component that lets you measure certain physiological characteristics of a muscle, analyse them and adapt the stimulation parameters accordingly. This tailoring of each programme to the changing condition of the muscle clearly enhances your comfort and the therapeutic efficacy.

Mi-SCAN

- This function adapts the electrostimulation session to the physiology of each user. Just before starting the work session, **Mi-SCAN** probes the chosen muscle group and automatically adjusts the stimulator parameters to the excitability of this area of the body. It is a truly personalised measurement.

N.B.: The **Mi-SCAN** function is operational as soon as a cable fitted with the **Mi-SENSOR** system (sold separately) is connected to the stimulator.


- This function is implemented at the start of the programme by a short sequence during which measurements are taken (a horizontal bar scans the figure at the left of the screen).
- Throughout the entire measurement test, you must remain completely still and be perfectly relaxed. The **Mi-SENSOR** system is very sensitive: the smallest contraction or the slightest movement could disturb the measurement test.
- Note that certain people might feel an unpleasant tingling sensation during the test.
- When the test is over, the symbol  appears and the programme can start.

VI HOW IT WORKS

You are strongly advised to carefully read the counter-indications and safety measures described at the start of this manual (chapter 1: "Warnings") prior to using your stimulator.

1. Description of the stimulator

A On/Off button

B "i" button used to:
increase stimulation energies in several channels
simultaneously

C Sockets for the 4 stimulation cables

D Stimulation cables

Channel 1 = blue Channel 2 = green
Channel 3 = yellow Channel 4 = red

E +/- buttons for the 4 stimulation channels

F Socket for the charger

G Rechargeable battery compartment

H Belt clip socket

2. Connections

Connecting the cables

The stimulator cables plug into the 4 sockets on the front of the device.

Four cables can be connected simultaneously to the four channels of the device. Both the sockets and the cables are colour-coded to simplify use and facilitate identification of the different channels.

blue = channel 1 green = channel 2
yellow = channel 3 red = channel 4

The Energy stimulator bears the label Mi -READY, in other words it can use certain Mi technology functions if a stimulation cable fitted with the Mi -SENSOR system is connected. This cable, sold separately, gives access to the Mi -SCAN function. It is connected to any one of the four stimulation channel sockets located on the front of the device and two electrodes are "clipped" on to it in the same way as for a conventional cable.

Connecting the charger

The Complex has considerable operating autonomy, as it uses rechargeable batteries.

To recharge them, use the charger supplied with your device and connect it to the front of the device, then plug the charger in a socket.

Remove any stimulation cables connected to the stimulator before recharging it.

We strongly recommend you fully charge the battery before using it for the first time as this will improve its autonomy and life span.

3. Preliminary settings

Before using the unit for the first time, you should select the working language of the device which is displayed on the options screen. Follow the instructions below. Afterwards, for the greatest comfort, the Complex offers you a number of setting options (operating language selection, display contrast setting, adjusting backlighting and volume setting). To change any of these settings, bring up the options screen by holding the On/Off button on the left of the stimulator for a few seconds.


fig.1

1b Use the channel 1 +/- button to choose the language you wish to use.

1c Use the channel 2 +/- button to adjust the contrast of the screen.

1d Use the channel 3 +/- button to adjust the volume.

1e Use the channel 4 +/- button to adjust the backlighting.

On: backlighting always on.

Off: backlighting always off.

Auto: backlighting activated whenever a button is pressed.

1a Press the On/Off button to validate the selected parameters. Your stimulator saves your options. It is now ready for use with the settings you selected.

4. Selecting a motivation

To switch on your stimulator, briefly press the On/Off button on the left of the Compex. A musical signal can be heard and a screen showing the different programme categories is displayed.

Before choosing a programme, you should select the desired category.


fig.1

2a Press the On/Off button to switch off the unit.

2b d Press the channel 1 and 3 +/- buttons to select a different motivation.

2e Press the channel 4 +/- button to confirm your choice and access the programme selection screen.

5. Selecting a programme

To choose a programme, it is particularly useful to consult chapter VII "Programmes and specific applications".

fig.3

3a Press the On/Off button to return to the previous screen.

3b Use the channel 1 +/- button to choose a programme.

3e Press the channel 4 +/- button to validate your selection and, depending on whether the programme bears either the symbol **START** or **▶▶▶▶▶**

a) **START** = the stimulation session starts immediately.

b) **▶▶▶▶▶** = a parameter setting screen appears.

6. Personalising a programme

The programme personalisation screen is not available for all programmes!


fig.4

4a Press the On/Off button to return to the previous screen.

4b Some programmes require manual selection of the muscle group you want to stimulate. This muscle group is shown in black on a small figurine above channel 1. Use the channel 1 +/- button to select your chosen group. The seven muscle groups proposed are shown successively in black on the small figurine:

Abdomen and lower back

Buttocks

Thighs

Legs and feet

Forearms and hands

Shoulders and arms

Thorax and back

Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

N.B.: If a cable fitted with the **mi-SENSOR** system (sold separately) is connected to the device, the muscle group is selected automatically.

4c Press the channel 2 +/- button to stop the warming-up sequence (closing the small animated rising convection symbol above the radiator).

4d Use the channel 3 +/- button to choose the programme difficulty level.

4e Use the channel 4 +/- button to confirm your choices and launch the programme.

7. During the stimulation session

Preliminary **mi-SCAN** test

*The **mi-SCAN** test is launched only if a cable fitted with the **mi-SENSOR** system (sold separately) is already connected to the stimulator when this is turned on!*

It is essential that you remain completely still and relaxed during the entire test.

If the **mi-SENSOR** system cable is connected, the test starts immediately after programme selection and personalisation.


fig.5

5a Press the On/Off button to interrupt the test at any moment. To relaunch the test from the start, press any of the 4 +/- buttons.


5e The **mi** logo above channel 4 indicates that the **mi-SENSOR** cable is connected to the stimulator. The figures in channel 4 will automatically change during the test as the sensor records different measurements. The 4 +/- buttons will be disabled for the duration of the test.

Once the test is complete, you will automatically be taken to the default programme start screen and prompted to increase the stimulation energies.

Adjusting stimulation energies

When you start a programme, you will be prompted to increase the stimulation energies. This is critical to the success of any treatment. To find out what energy

level needs to be reached for each programme, refer to the applications menu.

6a Press the On/Off button to place the unit in Pause mode.


6b c d e The Complex "bleeps" and the symbols of the four channels flash, changing from + to 000: the four channels are at 000 energy. You must increase the stimulation energy so that the stimulation can start. To do this, press the + buttons for the relevant channels until the desired setting is reached. If you want to increase the energy level of all four channels simultaneously, press the "i" button, located below the On/Off button. Press the "i" button twice to increase the levels in the first 3

channels, and 3 times to increase the levels in the first 2 channels. When you press the "i" button, the associated channels are highlighted in white on a black background.

Programme progression

Stimulation actually starts when the stimulation energy has been increased.

The examples reproduced below are sufficient to understand the general rules.

7a Press the On/Off button to interrupt the programme momentarily. To restart it, simply press the channel 4 +/- button. The session will resume at 80% of the energy levels that were in use prior to the interruption.

7 b c d e The different energies reached during the contraction phase are shown by a series of black bar graphs. Active rest phase energies are shown by hatched bar graphs.

N.B.: Active rest phase stimulation energies are automatically set at 50% of contraction energies. These can be modified during the rest phase. Once modified, they will be totally independent of the contraction energies.

N.B.: If your Compex emits a beeping sound and the symbols  under the active channels begin to flash, the stimulator is suggesting you increase the level of the stimulation energies. If you are working at the maximum tolerance level, simply ignore this message.

End of programme

At the end of each session, a small flag will be displayed on the screen and a short melody will be played. To switch off the stimulator, press the On/Off button.


fig.8

8. Consumption and recharging

Never recharge the stimulator without first disconnecting the stimulation cables.

Always use the charger supplied by Compex to recharge the batteries.

The Compex runs on rechargeable batteries. The battery's life depends on the programmes and stimulation energy used.

We strongly recommend you fully charge the battery before using your Compex for the first time as this will improve its autonomy and life span. If you do not use your device for a long period of time, please recharge the battery regularly.

Electricity consumption

The charge state of the battery is indicated by a small battery icon on the screen.


fig.9

If the battery icon contains just two lines, this means that power is running low.

Stop the session and recharge the unit.

If the symbol **START** normally displayed above the channel 4 +/- button is not visible and if the battery icon is flashing, this means that the battery is completely discharged. It is no longer possible to use the device. Recharge immediately.

Recharging

Remove all stimulation cables from the stimulator before recharging it. Connect the charger to the mains and then connect the stimulator to the charger. The charge menu shown below will automatically appear on the screen.


fig.10

The duration of the charging operation is shown on the screen (a complete charge may take 10 hours).

The battery icon is animated while the battery is recharging. When fully charged, the icon will be full and the total time taken to recharge the battery will flash on the screen. Simply disconnect the charger: the Compex will turn off automatically.

9. Problems and solutions

Electrode fault


fig.11

The Compex "bleeps" and alternatively displays the symbol of a pair of electrodes and an arrow pointing to the channel where a problem has been detected. In the example above, the stimulator has detected an error in channel 1.

- Check that electrodes are connected to this channel.
- Check whether the electrodes are old, worn, and/or the contact is poor: try using new electrodes.
- Try using the stimulation cable on a different channel. If the cable is still showing a fault, replace it.

Stimulation is not producing the usual sensation

- Check that all the settings are correct and ensure the electrodes are positioned properly.
- Change the positioning of the electrodes slightly.

The stimulation effect causes discomfort

- The electrodes are beginning to lose adhesion and no longer provide good contact on the skin.
- The electrodes are worn and need to be replaced.
- Change the positioning of the electrodes slightly.

The stimulator is not working

fig.12

- If an error screen appears while you are using the device, note the error number (in the example in fig. 12, error 1/0/0) and contact the nearest authorised customer support service.

VII PROGRAMMES AND SPECIFIC APPLICATIONS


The following applications are given by way of example. They will give you a better understanding of how electrostimulation sessions can be combined with voluntary activity. They will help you to determine the best procedure to follow according to your needs (programme selection, muscle group, time, electrode placement, body position). You can also choose a muscle group different to the one indicated, depending on the way you exercise or the area of the body that you want to stimulate.

The specific applications will provide you with information on electrode placement and the stimulation position to adopt. This information is given as a number (electrode positioning) and a letter (body position). These refer to the pictures showing placement of the electrodes and the pictograms showing the body positions on the inside cover of the manual.

If the Compex **is not fitted** with the **Mi-SENSOR** cable, the muscle group to be selected in the device is shown on the small figurine. If this does not appear (-), there is no need to choose a muscle group: selection is automatic.

If the Compex **is fitted** with the **Mi-SENSOR** cable (sold separately), the muscle group is selected automatically. There is therefore no need to look at the small figurine.

Example

For the 4 basic workout programmes: **Endurance, Resistance, Strength, Explosive strength**, we advise you to consult the training planner on our website www.compex.info. An interactive question and answer system takes you to a personalised training plan.

1. Your motivation: Sport

You are active and take part as often as possible in a physical or sporting activity. Even if you are not a specialist in a particular sport or a keen competitor, you are passionate about sport: you run, you swim, you cycle, sometimes you roller skate or just walk: in short, you take every opportunity to stay in peak physical condition! In this scenario, the **Sport** programmes will help you maintain your energy and vitality to benefit fully from all your exercise and sports, so as to stay on top form all year round.

Programmes

Endurance

Specially designed to improve oxygen absorption by the stimulated muscles, this programme allows you to exercise over a long period. It is the ideal top-up if you run, cycle, roller skate, etc.
Stimulation energy: highest tolerable energy (0-694)

Resistance

By subjecting your muscles to intense and prolonged effort, this programme leads to improved resistance to fatigue, during sustained workouts lasting a long time. It also develops certain muscles where good condition is important for everyday living (abdominal, dorsal, etc.).
Stimulation energy: highest tolerable energy (0-694)

Strength

This programme is very interesting for increasing your strength and muscle contraction speed during exercise or sport demanding very intense effort over a short time.
Stimulation energy: highest tolerable energy (0-694)

Explosive strength

Explosive strength, a key factor in some sporting activities of extremely brief duration, can be developed with this programme. It significantly increases the speed at which you will reach a high degree of strength and improves the efficiency of explosive gestures (bursts, shooting, sprinting, etc.).
Stimulation energy: highest tolerable energy (0-694)

Active recovery

Owing to the increased blood flow and the accelerated elimination of waste by the muscle contraction which it induces, this programme helps your muscle recovery. It has a relaxing and cramp-countering effect after a particularly protracted muscle workout or demanding exercise.
Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches

Muscle building

This programme will help you not only to obtain firm muscles which stand out, but also to develop in a balanced way your general musculature (strength, volume, tone). This way, you will get an attractive figure more quickly and more easily.
Stimulation energy: highest tolerable energy (0-694)

Body building

If you are seeking to increase your muscle mass, this programme is designed for you. In addition to a hard workout, this will help you to achieve the results you're looking for.
Stimulation energy: highest tolerable energy (0-694)

Combined workout

By subjecting your muscles to different workout sequences, this programme maintains your overall physical condition or reactivates your muscles after a period of rest and before more intense and targeted exercise or sport.

Stimulation energy: highest tolerable energy (0-694)

Demo

Refer to the quick start guide "Test your Complex in 5 minutes".

Specific applications

Your goal as a cyclist is to improve your endurance

To improve your endurance, you must develop to the maximum the oxygen supply to muscles stressed by sustained exercise. The **Endurance** programme will produce a significant increase in the consumption of oxygen by the muscles. The association of this programme with the **Oxygenation** programme, which leads to development of the intramuscular blood capillary network, is particularly interesting and will again improve your endurance.

Duration of the cycle: 8 weeks, 2-3 times a week

Example for 1 week

Mon.: Rest

Tue.: 1 x **Endurance** 8G 

Wed.: Cycling workout 1.5 hours, then 1 x **Oxygenation** 8G 

Thu.: 1 x **Endurance** 8G 

Fri.: Rest

Sat.: Cycling workout 60', then 1 x **Endurance** 8G 

Sun.: Cycling outing 2.5 hours, then 1 x **Oxygenation** 8G 

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

You want to improve the muscular qualities of your abdominal belt

The abdominal belt muscles play a determining role for efficiency of movement in most exercise and in sport.

However, improving the muscular qualities of the abdomen is all too frequently neglected owing to the tedious nature of voluntary exercises necessary. The **Resistance** programme provides a solution for this problem by imposing an exercise regime that is perfectly tailored to the physiological specifications of the muscular fibres that support the abdominal muscles.

Duration of the cycle: 6 weeks, 3 times a week

Programme: **Resistance** 10I  or 44I  (the entire abdominal belt), 11I  (rectus abdominus)

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

You want to prepare for a team sport before the season starts

During the preparatory pre-season period for team sports, specific muscle preparation is important. Stimulation of the quadriceps with the **Strength** programme will lead to an increase of the speed of contraction and muscular strength. The practical benefits will be obvious: improvement in starting and movement speed, jumping, kicking, etc. An **Active Recovery** session, taken after the most intensive workouts, accelerates the muscular recovery rate and reduces cumulative fatigue during the period of the season when the amount of training is high.

Duration of the cycle: 6-8 weeks, 3 times a week

Example for 1 week

Mon.: 1 x **Strength** 8G 

Tue.: Collective training, then 1 x **Active Recovery** 8G 

Wed.: 1 x **Strength** 8G 

Thu.: Collective training, then 1 x **Active Recovery** 8G 

Fri.: 1 x **Strength** 8G 

Sat.: Rest

Sun.: Collective training or friendly match, then 1 x **Active Recovery** 8G 

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

You want to improve the explosive strength of your thighs (long jump or high jump, sprinting, etc.)

To develop explosive muscular strength, voluntary training relies on tiring muscle training sessions that often include the risk of injury, since heavy weights are necessary. Incorporating the **Explosive Strength** programme in training eases muscular training sessions, obtaining both better results and more time for technical work.

Duration of the cycle: 6-8 weeks, 4 times a week

Example for 1 week

Mon.: 1 x **Explosive strength** 8G 

Tue.: Voluntary training on stadium

Wed.: 1 x **Explosive strength** 8G 

Thu.: Voluntary training with technical work on the trampoline

Fri.: 1 x **Explosive strength** 8G 

Sat.: Rest

Sun.: Voluntary training on stadium, followed by 1 x **Explosive strength** 8G 

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

You want to recover better and faster after exercise

A session with the **Active Recovery** programme, which accelerates drainage of toxins, will help you to relax your muscles and to restore your muscular balance faster, whatever the type of exercise taken. This programme is normally used between 15 minutes and 3 hours after exercise. If the result is insufficient, the session may be repeated immediately a second time. A session can also be repeated on the following day, in the event of persistent stiffness and cramp. An example of how to use this programme is shown above in the application "You want to prepare yourself for a team sport before the season starts".

You want to achieve harmonious shoulder musculature with moderate increase in muscle volume

Most voluntary exercise, such as jogging or cycling, puts little stress on the shoulder muscles. This is why it is particularly interesting to compensate this under use with electrostimulation sessions. The **Muscle Building** programme imposes a lot of specific work on the upper body muscles and thus ensures the harmonious development of the shoulders, with firm and well-defined muscles.

Duration of the cycle: 5 weeks, 4 times a week

Example for 1 week

Mon.: 45'-1 hour of voluntary exercise (jogging, swimming, cycling, fitness activity, etc.), then 1 x **Muscle Building** 17H 

Tue.: Rest

Wed.: 1 x **Muscle Building** 18C 

Thu.: 45'-1 hour of voluntary exercise (jogging, swimming, cycling, fitness activity, etc.)

Fri.: 1 x **Muscle Building** 17H 

Sat.: Rest

Sun.: 1 x **Muscle Building** 18C 

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

You want to develop your pectoral muscles

To obtain well-developed pectoral muscles, you must subject them to intensive work. In voluntary training, this type of work involves exercises with heavy loads. But this kind of work can be traumatic for the joints and the tendons involved.

Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

The **Body building** programme is the perfect solution for effectively increasing the muscular volume of the pectoral muscles, without the stress related to voluntary exercises.

Duration of the cycle: 5 weeks, 3 times a week

Example for 1 week

Mon.: Rest

Tue.: Voluntary training concentrating on the pectoral muscles, then 1 x **Body building** 43C , followed by 1 x **Oxygenation** 43D 

Wed.: Rest

Thu.: Voluntary training concentrating on the pectoral muscles, then 1 x **Body building** 43C , followed by 1 x **Oxygenation** 43D 

Fri.: Rest

Sat.: Voluntary training concentrating on the pectoral muscles, then 1 x **Body building** 43C , followed by 1 x **Oxygenation** 43D 

Sun.: Rest

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

You want to adopt a varied muscular training approach

Thanks to its eight sequences, which follow each other automatically, the **Combined Workout** programme imposes various kinds of stimulation on the muscles during the same session.

It cannot give priority to a particular muscular performance, but it achieves a correct level in all muscular work modes (endurance, resistance, strength, explosive strength).

Duration of the cycle: 6 weeks, 2-3 times a week

Example for 1 week

Mon.: 1 x **Combined Workout** 21D 

Tue.: Usual workout

Wed.: 1 x **Combined Workout** 21D 

Thu.: Usual workout

Fri.: 1 x **Combined Workout** 21D 
(optional)

Sat.: Usual workout or other exercise in the open air

Sun.: Rest

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

2. Your motivation: Health

Minor localised aches and pains and uncomfortable physical sensations are often the result of unavoidable and repetitive daily situations: incorrect posture sitting in front of a computer, prolonged standing, hours and hours spent in a car, a telephone regularly held between the shoulder and your cheek, etc. If stress and fatigue are added to all this, muscular tensions, back pains and other blood flow problems frequently appear. Using different stimulation principles, the **Health** programmes will cause your pains and unpleasant sensations to disappear, improving your well-being and your comfort in daily life.

Programmes

Pain management TENS

All painful phenomena, in particular pains due to sports-related pathological conditions, can be effectively alleviated using this programme that blocks the transmission of pain by the nervous system.

Stimulation energy: gradually increase the stimulation energy until a pronounced tickling sensation can be felt under the electrodes

Chronic pain

Thanks to its analgesic action, this programme provides rapid relief in case of chronic muscular pains (back, lower back pains, etc.). It reactivates the circulation, oxygenates and relaxes stiff muscles.

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches

Decontracture

This programme provides an excellent means for reducing sharp and localised muscular tensions (for example, lumbago or torticollis) and obtaining a relaxing and particularly pleasant effect.

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches

Oxygenation

Through a very strong increase of the blood flow and the improvement of the muscular oxygenation, this programme help to prevent muscular fatigue and maintain a long-lasting effort for a longer time. This programme also improves the recovery rate.

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches

Heavy legs

This programme is intended to eliminate the sensation of heaviness and discomfort in the legs which is felt in particular circumstances (standing for prolonged periods, excessively high ambient temperature, hormonal imbalance related to the menstrual cycle, etc.). It quickly brings a sensation of lightness and suppleness.

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches

Cramp prevention

By restoring the balance of the blood flow, this specific programme prevents the onset of cramps during the night or after physical effort.

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches

Relaxing massage

This programme is indicated for reducing unpleasant muscular tensions and creating a relaxing sensation of well-being, in particular after a tiring day. It induces a state of very pleasant optimal relaxation. Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches

*The **Pain management TENS, Chronic pain and Decontracture** programmes must not be used for a long period without medical advice.*

Specific applications

You have hand and wrist pains (e.g.: shooting pain caused by prolonged use of a computer mouse)

Remaining for a long time in certain inappropriate working positions and the micro-traumas caused by the repetitive use of a computer mouse can lead to a wide range of painful phenomena. Hand and wrist pain is very common and can be treated very effectively using the **Pain management TENS** programme.

Duration of the cycle: 1 week, min. once daily, then adapt according to any change in the pain. If required, the **Pain management TENS** programme can be repeated several times during the same day

Programme: **Pain management TENS** 39◇ -

You have muscular tensions in the back of the neck, middle of the back or lower back region

Pains in the back area or the back of the neck are due to a stiffness of the muscles, with accumulation of toxins. The pain can then become continuous or appear after only a few minutes spent working in an unfavourable position (for example, working in front of a computer screen). These chronic pains can be treated effectively with the **Chronic pain** programme which, in addition to combating pain, reactivates the circulation and relaxes the stiff muscles.

Duration of the cycle: 4 weeks, twice daily, with a 10 minute break between the 2 sessions

Programme: **Chronic pain** 15◇ ♠ (back of neck), 13◇ ♠ (middle of the back) ou 12◇ ♠ (lower back)

You have recent and acute muscular pains in the back of the neck or the lower back region (torticollis, lumbago)

Due to a sudden or inappropriate movement, an incorrect sleeping position or cold weather, a sharp pain may be felt in the back of the neck or lower back area. The person then presents with an acute and intense contracture of the muscles in this area and sharp pain. In such a situation, you must always consult a doctor to receive appropriate treatment. As well as this, the specific **Decontracture** programme can help to relax the muscles and eliminate pain.

Duration of the cycle: 1 week, daily

Programme: **Decontracture** 15◇ ♠ (torticollis) or 33◇ ♠ (lumbago)

You want to prepare for a seasonal endurance activity (e.g.: hiking, bicycle touring, cross-country skiing)

Physical activities that last a long time require strong muscles with good capillary circulation, so that the muscular fibres can be properly oxygenated.

If you don't exercise, or not regularly enough, this lack of muscular quality makes exercise uncomfortable and limits the speed of recuperation. It is also responsible for many unpleasant side effects (feeling of heaviness, stiffness, cramps or swellings). To restore your muscles' endurance and ensure full activation of the capillary system, the **Oxygenation** programme offers a specific and very effective mode of stimulation.

Duration of the cycle: 6-8 weeks before beginning the hike, 3 times a week, while alternating the muscular groups one day in two

Programme: **Oxygenation** 80  and **Oxygenation** 25 

You want to prevent the feeling of heaviness and discomfort in your legs

The feeling of heaviness in the legs is due to a temporary deficiency in the venous return. The insufficient supply of oxygen to the tissue, and in particular the muscles, that results from this situation, leads to this uncomfortable sensation in the legs. Use the **Heavy legs** programme to accelerate the venous return and obtain a significant relaxing effect on painful muscles.

Duration of the cycle: to be used on the calves, as soon as an occasional sensation of heaviness in the legs occurs

Programme: **Heavy legs** 25J -

You want to prevent calf cramp

Many people suffer from cramp in the calve muscles, which can appear spontaneously during rest at night or as a result of prolonged muscular effort.

The phenomenon of cramp may in part be caused by an imbalance in blood flow in the muscles. To improve the blood flow and prevent the onset of cramp, the **Compex** offers a specific stimulation programme that will allow you to obtain positive results.

Duration of the cycle: 5 weeks, daily, at the end of the day or in the evening

Programme: **Cramp prevention** 25 

You want to prevent uncomfortable muscular tension in the back of the neck

Remaining seated for long periods, associated with repetitive arm movements (as is frequently the case in front of a computer screen), can provoke increased tension in muscles in the back of the neck. Any other stress condition can also lead to a state of excessive muscular tension, which is often responsible for painful or uncomfortable sensations. The in-depth effects produced by the **Relaxing massage** programme ease such painful sensations, with a particularly relaxing effect.

Duration of the cycle: to be used on the muscles in the back of the neck, whenever an occasional sensation of painful muscular tension is felt; to be repeated, if necessary, in cases of particularly strong muscular tension

Programme: **Relaxing massage** 15 

3. Your motivation: Aesthetic

You are concerned about your well-being and your appearance, and want to have a well-balanced body. However, a sedentary lifestyle, an inadequately balanced diet and tight schedules may all be factors with a negative impact on the harmony of your figure. The targeted effects of the **Aesthetic** programmes will help you in your efforts to maintain a healthy lifestyle and benefit from the positive effects of intense muscular activity: well-toned muscles, a shapely figure and a firm body. Used regularly, these programmes will help you to regain or keep your figure.

Programmes

Toning

Using this programme will help you to tone your muscles and prepare them before the more intensive work of firming. You will quickly recover a toned, firm and shapely body.

Stimulation energy: highest tolerable energy (0-694)

Firming

This programme is very useful if you want to increase firmness, shape your muscles and restore the elastic quality of the supporting tissues. It also compensates for insufficient muscular activity, by making certain flabby areas work harder.

Stimulation energy: highest tolerable energy (0-694)

Shaping

This programme is your partner if you want to define and sculpt your body, when your muscles are already firm. Your figure will be visibly reshaped.

Stimulation energy: highest tolerable energy (0-694)

Abs

Combining the beneficial effects of the **Firming** and **Shaping** programmes, matched to the specificities of the muscles of the abdominal belt, this programme will enable you to reduce your waistline, while improving the muscular support of the abdominal wall. Stimulation energy: highest tolerable energy (0-694)

Buttocks

This programme subjects the buttock muscles to the most appropriate work for improving your figure. Its hyper-tensor effect is the ideal solution to restore or improve the tone of the buttocks region. Stimulation energy: highest tolerable energy (0-694)

Specific applications

You want to firm your thighs and sculpt your body

The **Toning**, **Firming** and **Shaping** programmes impose perfectly adapted and progressive work on the muscles.

This very intense muscular activity (hundreds of seconds of sustained contractions) first tones and firms your muscles, then redefines their contours to sculpt your body. In this example, the person wants to work on her thighs, but it is of course possible to stimulate other muscles as well.

Moreover, this stimulation procedure can also be applied to a number of muscle groups during the same cycle.

Duration of the cycle: 12 weeks

Weeks 1-4: 4 x **Toning** 8G \hat{A} /week

Weeks 5-8: 4 x **Firming** 8G \hat{A} /week

Weeks 9-12: 4 x **Shaping** 8G \hat{A} /week

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 **Shaping** session per week at the last level reached.

You want to reduce your waist and tone your abdominal belt

The **Abs** programme is specially designed to improve the support and firmness of the muscles of the abdominal belt. The different workout schemes, combined with a toning massage during the rest phases, guarantee beneficial results for your figure.

Duration of the cycle: 3 weeks, 4 times a week

Programme: **Abs** 49I -

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

Your target is to improve the firmness and shape of your buttocks

The considerable amount of progressive and specific work of the muscles in the buttocks region produced by the **Buttocks** programme effectively increases the firmness of these muscles. The shape of the buttocks improves gradually, with beneficial effects on the whole figure.

Duration of the cycle: 3 weeks, 4 times a week

Programme: **Buttocks** 9A,B,L -

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

VIII EMC TABLE (Electromagnetic Compatibility)

The Compex Energy needs special EMC precautions and must be installed and started according to the EMC information supplied in this manual.

All RF wireless transmission systems can affect the Compex Energy.

The use of accessories, sensors and cables other than those recommended by the manufacturer may result in stronger emissions or reduce the immunity of the Compex Energy.

The Compex Energy should not be used beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the Compex Energy works properly in the chosen configuration.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER CONCERNING ELECTROMAGNETIC EMISSIONS

The Compex Energy is intended for use in the electromagnetic environment specified below.

The customer or user of the Compex Energy should ensure that it is used in this environment.

Emissions test	Compliance	Electromagnetic environment - guide
CISPR 11 RF emissions	Group 1	The Compex Energy uses RF energy only for its internal operation. Consequently, its RF emissions are very low and are unlikely to interfere with any adjacent electrical device.
CISPR 11 RF emissions	Class B	The Compex Energy is suitable for use in any establishment, including a private dwelling and a place connected directly to the low voltage mains supply which powers residential buildings.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ emission oscillations IEC 61000-3-3	Not applicable	

RECOMMENDED SPACING BETWEEN A PORTABLE AND MOBILE COMMUNICATION APPLIANCE AND THE COMPEX ENERGY

The Compex Energy is designed for use in an electromagnetic environment in which radiated RF waves are controlled. The buyer or user of the Compex Energy can contribute to preventing electromagnetic interference by maintaining a minimum distance between RF portable and mobile communication appliances (transmitters) and the et Compex Energy according to the table of recommendations below and according to the maximum output power of the telecommunication appliance.

Maximum transmitter output power W	Spacing according to the frequency of the transmitter mCISPR 11		
	From 150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	From 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	From 800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

In the case of transmitters whose maximum output power is not shown in the table above, the recommended spacing of d metres (m) can be calculated using the appropriate equation for the transmitter frequency, where P is the maximum output power of the transmitter in watts [W] as set by the transmitter manufacturer.

NOTE 1: At 80 MHz and at 800 MHz, the spacing for high frequency amplitude is applied.

NOTE 2: These guidelines may not be appropriate for some situations. Electromagnetic wave propagation is modified by absorption and reflection due to buildings, objects and persons.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER -

Compex Energy is designed for use in the electromagnetic environment stipulated below.

The buyer or user of the Compex Energy must ensure it is used in this recommended environment.

Immunity test	Test level IEC 60601	Observance level
Electrostatic discharge (DES) CEI 61000-4-2	± 6 kV at the contact ± 8 kV in air	± 6 kV at the contact ± 8 kV in air
Rapid electrical transients in bursts of noise CEI 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable System battery-powered
Surge currents CEI 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Not applicable System battery-powered
Voltage trough, short power cuts and voltage variations on power supply input lines CEI 61000-4-11	< 5% VT (trough > 95% of UT) for 0.5 cycle < 40% VT (trough > 60% of VT) for 5 cycles < 70% VT (trough > 30% of VT) for 25 cycles < 5% VT (trough > 95% of VT) for 5 seconds	Not applicable System battery-powered
Magnetic field at mains frequency (50/60 Hz) CEI 61000-4-8	3 A/m	

NOTE: VT is the AC supply voltage before application of the test level.

- ELECTROMAGNETIC IMMUNITY

Compex Energy is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Energy must ensure it is used in this recommended environment.

Electromagnetic environment - guide

Floors must be wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity must be maintained at a minimum of 30%.

The quality of the power supply should be that of a typical commercial or hospital environment.

The quality of the power supply should be that of a typical commercial or hospital environment.

The quality of the power supply should be that of a typical commercial or hospital environment.
If the Compex Energy user requires continuous operation during mains power cuts, we recommend that the Compex Energy is powered by a UPS or a battery.

Magnetic fields at the mains frequency should be at the level of a representative site located in a typical commercial or hospital environment.

NOTE: VT is the AC supply voltage before application of the test level.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER -

Compex Energy is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Energy must ensure it is used in this recommended environment.

Immunity test	Test level IEC 60601	Observance level
RF immunity IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Signal lines Not applicable for power supplies of battery-powered systems
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz 10 V/m 26 MHz to 1 GHz	3 A/m 10 V/m

NOTE 1: From 80 MHz to 800 MHz, the spacing for high frequency amplitude is applied.

NOTE 2: These guidelines may not be appropriate for some situations.

- a** The field intensity from fixed transmitters, such as radio telephone base stations (cellular/wireless) and a mobile radio, amateur radios, AM and FM radio emissions and TV emissions cannot be predicted with any accuracy. It may therefore be necessary to consider an analysis of the electromagnetic environment of the site to calculate the electromagnetic environment coming from RF fixed transmitters. If the field intensity measured in the environment where the Compex Energy is located exceeds the appropriate RF observance level above, the Compex Energy should be monitored to ensure it is operating properly.

NOTE: VT is the AC supply voltage before application of the test level.

- ELECTROMAGNETIC IMMUNITY

Compex Energy is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Energy must ensure it is used in this recommended environment.

Electromagnetic environment - recommendations

Portable and mobile RF communication devices must only be used relative to the Compex Energy and its wiring at a distance which is not less than the spacing recommended and calculated using the appropriate equation for the transmitter's frequency. Recommended spacing

$$d = 1.2 \sqrt{P}$$

$$d = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$$

$$d = 2.3 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$$

Where P is the maximum output power of the transmitter in watts (W) set by the manufacturer's specifications and where d is the recommended spacing in metres (m). The field intensity of RF fixed transmitters, as determined by an electromagnetic survey^a must be less than the observance level to be found in each frequency range^b.

Interference may occur close to any appliance identified by the following symbol:

Electromagnetic wave propagation is modified by absorption and reflection due to buildings, objects and persons.

In the event of abnormal operation, new measures may then be imposed, such as realignment or movement of the Compex Energy.

b Above the frequency amplitude from 150 kHz to 80 MHz, the field intensity must be 3 V/m.

NOTE: VT is the AC supply voltage before application of the test level.

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© 12/2011 DJO, LLC
Compex® is a registered trademark of DJO, LLC.