
# GLOBUS

ITALIAN EXCELLENCE

-  **ELETTROSTIMOLATORI** - Manuale di utilizzo
-  **ELECTROSTIMULATORS** - User Manual
-  **ELECTROESTIMULADORES** - Manual de usuario
-  **ÉLECTROSTIMULATEURS** - Manuel de l'utilisateur

# GLOBUS

ITALIAN EXCELLENCE

	MANUALE DI UTILIZZO	pag. 3
	USER MANUAL	pag. 37
	MANUAL DE USUARIO	pag. 71
	MANUEL DE L'UTILISATEUR	pag. 105

# ITALIANO

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# ENGLISH

DEAR CUSTOMER, thank you for choosing a GLOBUS product. We remain at your entire disposal for any assistance or advice you may need.

The device has been manufactured in compliance with applicable technical regulations and has been certified in compliance with Directive 93/42/EEC as amended by 2007/47 Directive for medical devices, by Kiwa Cermet Italia S.p.A. Body n. 0476, in order to ensure product safety.

## TECHNICAL FEATURES

### Device

Size: 13 X 8 X 2 cm

Weight: 220 g

Container: ABS

Protection level: IP 22

Storage and shipping temperature: from -10°C to 45°C

Max. relative humidity: 30% - 75%

Made in PRC

### Charger

Type: FLO Switching Power Adaptor

Model: DK7-065-0200-EU

Input: 100-240 Vac - 50-60 Hz - 0,06 A

Output: 5,8 Vdc - 200 mA

Polarity: 

Battery pack: Ni-MH AAA 4.8 V 800 mAh

### Use conditions

Temperature: from 0°C to 35°C

Maximum relative humidity: from 15% to 93%

Atmospheric pressure: from 700 hPa to 1060 hPa

The values represent the limits allowed if the product or the accessories are not in their original package.

### Technical features of EMS and TENS currents

Available outputs: Channels 1-2

Constant current: Yes

Intensity: 0-100 mA peak-peak per channel

Waveform: rectangular, biphasic, symmetrical, compensated

Work frequency: 1-150 Hz

Recovery frequency: 1-150 Hz

Pulse amplitude: 50-400  $\mu$ s

Working time: from 1 to 30 seconds

Recovery time: from 0 to 1 minute

Frequency modulation range: continuous variation from 1 to 150 Hz

Minimum modulation time: 3 seconds

Amplitude modulation range: continuous variation from 50 to 400  $\mu$ s

### Technical features of iontophoresis currents

Available outputs: Channel 1

Constant current: Yes


Min. intensity: 0 mA/1000 Ohm

Maximum intensity: 10 mA/1000 Ohm step 0.1 mA/1000 Ohm

Treatment time: 1 minute

Max. time: 99 minutes

## EQUIPMENT

A. Device

B. 2 cables for electrode connection (for EMS, TENS and IONTOPHORESIS treatments)

C. Bag with 4 self-adhesive, re-usable electrodes (50x50 mm)  
(We recommend the use of these electrodes for small surfaces such as upper limbs, calves, cervical area...)

D. Bag with 4 self-adhesive, re-usable electrodes (50x90 mm)  
(We recommend the use of these electrodes for large surfaces such as thighs, abdomen, gluteus...)

E. Battery charger (see technical features)

F. Manual and warranty

G. Carrying case

The electrostimulation device is supplied with all the necessary cables and electrodes; therefore, when you open the package, check that the basic equipment is complete. If some elements are missing, contact your authorized retailer immediately.

Check that the device and the electrodes are intact.

The device can be combined with some optional accessories (see the features on our website [globuscorporation.com](http://globuscorporation.com)). To purchase these accessories, please contact your dealer.

Accessories that are not included in the equipment (to be purchased separately)

- Motor point pen
- Kit with 8 elastic bands for legs and thighs
- Kit with 4 elastic bands for thighs
- Face electrodes
- Splitter cables kit
- Gel
- Rectangular electrodes for iontophoresis (60x85 mm)
- Fast band
- Fast pad

- Probe series - Anal probes (periprobe anal probe - analia)
- Probe series - Vaginal probes (periprobe vaginal probe - novatys)

## **INTENDED USE**

The product has an estimated usable life of 5 years. It is advisable to return the device to the producer and/or authorized center every 2 years to perform maintenance and security checks. The number of treatments depends on the battery charge.

The estimated duration of the battery is of 6 months. Thereafter, it is advisable to replace it. The electrostimulators are designed to be used in the following operating environments:

- Domestic environment;
- Rehabilitation in general;
- Pain treatments in general;
- For aesthetic and sport purposes;

The use of this equipment is allowed to the patient (duly informed about the methods applied by the device).

## **CONNECTIONS**

### **Cable connection**

It is possible to use one or two cables, depending on how many channels are used (one or two). To connect the cables, insert the cable connector in the specific slot on the top of the unit. The slots are placed exactly above the corresponding (left or right) channel. Connect the cable to the channel you want to use.

For iontophoresis, use only channel 1.

### **Application of the electrodes**

Take the electrodes from the original packaging; all new electrodes have a sealed packaging. Ensure that the device is off.


To start, connect the two plugs of the cables to the electrodes, then remove the electrodes from their place and put them on the skin. To place the electrodes correctly, see the images at the end of this manual. After use, place the electrodes back in their specific place.

**WARNING:** do not unplug the electrodes if the device is operating.

## LABELLING AND SYMBOLS

	Warning.
	This symbol on your device indicates that it complies with the directives for medical devices (93/42/EEC 47/2007EEC). The number of the notified body is 0476.
	It indicates that this is a II class device.
	It indicates that this device has BF type parts.
	WEEE symbol (Waste of Electrical and Electronic Equipment). Recycling symbol. The WEEE symbol used for this product indicates that it cannot be treated as a household waste. Proper disposal of the product will contribute to protecting the environment. For more information on the recycling of this product, please contact the concerned office of your local body, the household waste management company or the store where the product was purchased.
	It indicates that the product has been produced in compliance with the directive 2011/65/EEC.
	It indicates the optimal temperature for the storage and transportation of the product.
	It informs the operator that it must read the manual before using the device.
IP 22	It indicates water protection.
	It indicates a compulsory behaviour.
	It indicates the pressure of the environment in which the device and the accessories are transported and stored.
	It indicates the humidity of the environment where the device and the accessories are used and stored.

	It indicates the manufacturer.
	It refers to the expiry date of the product.
	It refers to the production lot.
	It refers to the manufacturing date.
	Use indoor.

### Power supply

 **F/D<sup>®</sup> Switching Power Adaptor**  
**Model: DK7-065-0200-EU**  
**Input: 100-240Vac 50-60Hz 0.06A**  
**Output: 5.8Vdc 200mA**  
**Made in PRC** 


## Electrodes


- Clean and degrease the skin.
- Do not apply the electrode on wounds or injured skin.
- Connect the cable connector to the electrode connector.
- Remove the electrode.
- Apply on the skin.
- Start the program.
- At the end, turn off and put the electrode back in the package.
- Electrodes are for personal use.
- Do not remove the electrode by grabbing the connector.
- Electrodes should not touch each other.
- Do not apply the electrodes on the temples, the neck and in a transthoracic position.
- Do not leave the electrodes in the car.

## Display

## Keyboard

## WARNING AND CONTRAINDICATIONS

### Mandatory behavior

For maximum safety, the device must be used following the rules and the limitations of the user manual.

Treatments should not be performed on skin lesions.

If the package, the cable or the connector of the power supply show signs of wear or damage, replace them instantly.

The device should be connected to the mains with its power supply. Before doing this, make sure that power system specifications comply with the directives in force in your country. Make sure that the power supply is in a comfortable position and will be easily removed.

The manufacturer declines any responsibility for any use that differs from what is indicated and prescribed in this manual.

The full or partial reproduction of the texts and/or pictures contained in this manual by any electronic or mechanical means without prior written authorization of the manufacturer is forbidden.

### Warnings before the use

Do not use the device with other electronic devices simultaneously, especially if they maintain vital functions. See the enclosed tables for a correct use of the electromedical device. If the device is used nearby or upon other devices, make sure it works properly.

- It is recommended to read carefully the entire operating manual before using the device; keep carefully this operating manual.

- The device is capable of delivering current values exceeding 10mA.
  - Before each use always check the integrity of the device. This is a fundamental requirement for carrying out the therapy; do not use the device if the buttons or the cables are defective or malfunctioning. It should be used only by people over 18 and, in any event, who are able to understand and take action.
  - It should not be used for purposes other than transcutaneous neurostimulation.
  - It must be used following the indications and under the strict control of the physician or qualified physiotherapist.
  - The device must be used with the transcutaneous neurostimulation electrodes suitable for this use.
  - The device must be kept out of the reach of children.
  - With its current, it can disturb ECG monitoring devices.
  - It must not be used in a transthoracic mode as it could cause cardiac arrhythmia by superimposing its frequency to that of the heart. (Do not perform the treatment on the chest and the back simultaneously).
  - If there is any health problem, it must be used only after consulting a physician.
  - A simultaneous connection of a patient to a high frequency electrosurgical device can cause burns near the electrodes of the stimulator and therefore the stimulator may be damaged.
  - Once you have turned the device on, make sure the display shows the software version and the device model: it means that the device is working and ready to be used;
- If it does not, or the display does not show all the segments, turn it off and on again. If the problem persists, contact our help center and do not use the device.
- The sudden shutdown shortly after the starting indicates a low battery level. Recharge as reported in the section "HOW TO CHARGE THE BATTERIES".

### **Warnings during the use**

While using the electrostimulator, some indications should be followed:

- In the case of damaged cables, they must be replaced with original parts and not used anymore.
- Only use Globus marked electrodes.
- Great attention has to be paid when current densities exceed  $2\text{mA}/\text{cm}^2$  (effective value) for each electrode.

The device must be kept out of the reach of pets, as they could damage the device and contaminate the electrodes and other accessories with parasites.

- The cables of the electrostimulator should not be wrapped around people's neck to avoid any risk of strangulation and suffocation.
- Mobile and fixed radio communication devices could influence the functioning of the electromedical device: see the tables attached to this manual.

Special precautions for incontinence treatments.

- Patients with extra-urethral incontinence should not be subjected to treatment with the stimulator.

- Patients suffering from excessive incontinence due to evacuation disorders should not be treated with the electrostimulator.

- Patients with severe urinary retention in the upper urinary tract should not be treated with the stimulator.

- Patients with complete peripheral denervation of the pelvic floor should not be treated with the stimulator.

- Patients suffering from a total/subtotal prolapsed uterus/vagina must be stimulated with extreme caution.

- Patients with urinary infections should be treated for these symptoms before starting the stimulation with the stimulator.

- Before removing or touching the probe, it is necessary to turn the stimulator off or adjust the intensity of both channels to 0.0 mA.

- The treatment is a personalized medical prescription: do not lend the stimulator to other people.

## **SIDE EFFECTS AND CONTRAINDICATIONS**

Isolated cases of skin irritation have been observed in people with high skin sensitivity.

In case of allergic reaction to electrode gel, suspend the treatment and contact a specialist.

If during the treatment signs of tachycardia and extrasystole appear, suspend the treatment and contact your physician.

### **Contraindications**

Use is not recommended in the following cases:

- Stimulation of the anterior neck (carotid sinus).
- Patients with pacemakers.
- Patients with tumor diseases (consult the oncologist).
- Stimulation of the brain region.
- Pains with unknown etiology.
- Sores and dermatological diseases.
- Acute traumas.
- Stimulation on recent scars.
- Pregnancy.

- It is strictly forbidden to use the electrostimulator on the ocular area.
- Near body areas with metallic implants or infra-tissue metals (e.g. prostheses, osteosynthetic devices, coils, screws, plates), when using monophasic currents such as interferential and continuous currents (ionophoresis).

It is also recommended to use the device with caution in case of capillary fragility of the patient, as excessive stimulation could increase capillary ruptures.

## **MAINTENANCE AND CLEANING**

### **Maintenance and cleaning of the device**

- In case of real or alleged malfunctioning, do not tamper with the device and do not try to repair it by yourself.
- Do not intervene on the device and do not open it. Only specialized and authorized centers can repair it.
- Avoid violent impacts that may damage the device and lead to malfunctioning, even if it cannot be detected immediately.
- Use the electrostimulator in a dry and open environment. (Do not wrap it with other objects)
- Clean the electrostimulator only with disinfectant with sodium hypochlorite or quaternary ammonium salt (percentage: 0.2 - 0.3%) diluted with distilled water. After cleaning and disinfecting the device, dry it with a clean cloth.
- Always use the device with clean hands.
- It is recommended to use the device in a clean environment to avoid contamination with dust and dirt.
- It is recommended to use the device in a ventilated, well-aired space.
- The parts should be cleaned/disinfected after every use, unless otherwise indicated.

## **Accessories**

### **Use and storage of the electrodes and the cables**

After using multi-use single-patient and/or single-use electrodes, they must be stored using their plastic film and placed in their plastic bag.

Make sure that the electrodes do not touch each other and do not overlie one over the other.

Once the package has been opened, the electrodes can be used for 25-30 applications.

The electrodes must be always used with clean hands and they must be replaced if they are not perfectly in contact with the skin.

If non self-adhesive electrodes are used, it is advisable to clean their surface with proper cleansers that respect the requirements described in the manual.

The electrodes must be stored in their bag in an environment that respects the requirements described in the manual.

At the end of the treatment, unplug the cables from the connectors and clean them carefully with proper cleansers that respect the requirements described in the manual.

After cleaning and drying them, fold them up and place them in the plastic bags supplied with the cables.

### **How to charge the batteries**

The device is equipped with a nickel-metal hydride rechargeable battery pack (4.8V 800mAh).

When the display shows the low battery symbol, you should recharge it immediately. After the symbol appears, the device has a limited operating time and may not be able to execute the program completely. After turning the device off and disconnecting the electrodes, connect the device to the supplied charger inserting the connector in the specific inlet.

It is recommended to leave the device on charge for 8-12 hours to guarantee a longer duration on the battery pack.

The device can operate for about 15-20 hours, after the charge, depending on the executed treatments. It is recommended to charge the device for the above-mentioned time, in which it is possible to carry out about 70 treatments.

A possible warming of the battery area during the charge is to be considered as normal.

Never use a charger that is different from the one supplied with the device.

Please contact our customer care to replace the battery pack.


## **Disposal of the device**

Do not throw the device or parts of it into the fire, but dispose of the product in the specialized centers and respecting the regulations in force in your country. When the product has to be disposed of, the user can give it back to the retailer when purchasing a new device.

Correct waste separation or following the above-mentioned procedure contribute to avoid possible negative effects on environment and health and promote the reuse and/or recycle of device materials. The illegal disposal of the product entails the application of administrative fines in accordance with applicable regulations.

# INSTRUCTIONS FOR USE

## Programming

Press Fn for 2 seconds in the start screen to program the device with 6 modes (PED1 - PED6).

Press OK to pass from a parameter to the next one.

It is possible to change the duration of all programs from 5' to 100' with steps of 5'. Default duration 20'.

### PED1

It is possible to set pulse frequency and duration.

From 1 to 150 Hz with steps of 1 Hz.

From 50 to 250  $\mu$ s with steps of 10  $\mu$ s.

Default values: 100 Hz / 200  $\mu$ s.

### PED2

It will be possible to set pulse duration with a fixed frequency of 100 Hz.

From 50 to 250  $\mu$ s with steps of 10  $\mu$ s.

Default values: 200  $\mu$ s.

### PED3

It is possible to set frequency values.

From 1 to 150 Hz with steps of 1 Hz.

Default values: 100 Hz

Modulated pulse duration from 100  $\mu$ s to 250  $\mu$ s.

### PED4


It is possible to set pulse frequency and duration.

From 1 to 150 Hz with steps of 1 Hz.

From 50 to 320  $\mu$ s with steps of 10  $\mu$ s.

It is possible to change the following values:  $t_{UP} = 1$  s;  $t_{HOLD} = 5$  s;  $t_{DOWN} = 1$  s;  $t_{OFF} = 5$  s.

Default values: 30 Hz / 250  $\mu$ s.

### PED5

It is possible to set frequency values.

From 1 to 150 Hz with steps of 1 Hz.

Default values: 30 Hz

Modulated pulse duration from 100  $\mu$ s to 320  $\mu$ s.

## PED6

It is possible to set thold and toff values.

From 1 to 30 s with steps of 1 s.


It is possible to set the frequency value.

From 1 to 150Hz with steps of 1 Hz.

Pulse duration fixed at 300  $\mu$ s.

Default values: tHOLD = 2 s; tOFF = 10 s; 30 Hz.

tUP and tDOWN cannot be changed and have a duration of 1s.

## ELECTROSTIMULATION

### The advantages of using electrostimulation

Electrostimulation is not meant to replace physical activity, but should be considered as a complementary treatment.

Depending on the aims, it will be used by:

- Athletes, to integrate the normal training and improve performances.
- People who want to intervene on some imperfections.
- People who suffer from certain pathologies.
- People who suffered a trauma or during rehabilitation.

Electrostimulation is a technique that, by means of electric pulses that act on muscle motor points (motoneurons) or nerve endings (TENS pulses), causes a muscular contraction that is similar to voluntary contraction.

There are two methods of use:

- Muscle stimulation (ideal for strength development and aesthetic treatments).
- Nerve ending stimulation (ideal for pain treatments).

### Muscle Types

Muscles can be divided into three different types: striated or voluntary muscle, cardiac muscle and smooth or involuntary muscle.

Voluntary muscles include the skeletal muscle that moves the skeleton and the skin muscle that acts on facial expression.

The cardiac muscle and the smooth muscle are not controlled voluntarily.

The majority of human muscles are striated or voluntary, with approximately 200 muscles on each side of the body. (About 400 on the whole).

Skeletal muscles are the target of EMS.

## **The physiology of muscular contraction**

The skeletal muscle performs its functions through contraction. Muscle contraction causes joint movement and therefore skeleton movement. Muscle contraction takes place as follows: when a person decides to make a movement, the electric potential in the brain motor center is changed and an electric signal is sent to the muscle that will have to contract.

The electric signal is sent as a potential (voltage) change and an electric pulse passes through the motor nerve to the muscle that has to be stimulated.

After the pulse is received, Ca ions are released causing actin-myosin rapprochement. As a result, the distance between Z bands decreases causing muscle contraction

The necessary energy is provided by sugar and fat reserves of the human body. In other words, electric stimulation is not a direct source of energy, but acts as a tool that triggers muscle contraction.

The same mechanism is activated when muscle contraction is produced by EMS. In other words, they act as a natural pulse sent by the motor nervous system.

The muscle usually relaxes and goes back to its original state at the end of the contraction.

### **Isotonic and isometric contraction**

Isotonic contraction occurs when, in a body movement, muscles produce a state of constant tension. On the other hand, when the articular ends of a body part are blocked, muscles produce a tension that is called isometric contraction.

In the case of electrostimulation, isometric contraction is preferred as it can be controlled more effectively.

### **The distribution of the different fiber types in the muscle**

The relationship between the two main categories (type I and type II) can vary in a considerable way. There are muscular groups that are typically made up of type I fibers, like the soleus, and muscles that only have type II fibers, like the orbicular muscle; however, the majority of human body muscles is composed of a combination of the two types.

Studies on the distribution of fibers in the muscle have highlighted the close relationship between the (tonic or phasic) motoneuron and the functional features of the fibers it innervates; moreover, they have proved that a specific motor action (particularly in sports) can lead to a functional adaptation of fibers and change their metabolic features.

Motor unit type	Contraction type	Contraction frequency
Tonic ST	slow contraction I	0 - 50 Hz
Phasic FT	fast contraction II	50 - 70 Hz
Phasic FTb	fast contraction II b	80 - 120 Hz

To bring a tissue from rest to excitation with an induced electric stimulation (electrostimulator pulse), some conditions are necessary:

- The current amount
- The stimulation duration, which must be suited for the body part to be stimulated

On this basis, for tissue excitation, the relationship between the duration and the intensity of the stimulation is particularly important, not only the peak value of the intensity

This relationship changes for all the different muscles.

### **Recommended intensity**

The intensity value required to trigger contraction depends on patients, electrode placement, adipose tissue, perspiration, possible hairs on the treatment surface etc. Therefore, the same current intensity may give different sensations to different people, in different days or body sides. It is advisable to regulate the intensity during the same session to contrast accommodation, in order to obtain the same contraction.

The current intensity for the different phases is suggested with an approximate value and can be modified based on individual sensations.

- Moderate: the muscle is never tired, even in long treatments. The contraction is absolutely bearable and agreeable. First level of the intensity graph.
- Intermediate: the muscle is visibly contracted but does not provoke any joint movement. Second level in the intensity graph.
- High: the muscle undergoes a considerable contraction. Muscle contraction would extend or bend the limb if it were not blocked. Third level in the intensity graph.
- Maximal: the muscle undergoes a maximal contraction. This is a demanding process that should be performed only after several applications.

Moderate	from 10mA to 20mA
Intermediate	from 20mA to 30mA
High	above 30mA
Maximal	On the verge of the tolerance limit, always under the threshold of pain.

In treatment description, recommended intensity levels are indicated.  
NOTE Recommended current levels are just an indication.

## **Tens**

Transcutaneous Electrical Nerve Stimulation (TENS) is a selective stimulation of the large fibers of peripheral nerves, which contributes to closing the gate entrance to pain pulses and increases the release of endorphinic substances, reducing pain intensity. Therefore, TENS is particularly indicated to treat the severe and chronic pain caused by the main musculoskeletal disorders.

TENS currents reduce pain thanks to the following factors:

- a. Gate control theory
- b. Endorphin secretion
- c. Different sedative effects depending on frequency

## **Gate theory**

If the electric signals leading pain information to the brain are blocked, the perception of pain is eliminated. For instance, if we hit our head against an object, the first thing we do is to massage the traumatized area. In this way, we stimulate the receptors of touch and pressure. TENS in continuous mode and frequency modulation can be used to generate signals that are similar to those of touch and pressure. With sufficient intensity, their priority is so high that it prevails on pain signals. When the priority is obtained, the gate of sensory signals is opened and the pain gate is closed, impeding the passage of these signals to the brain.

## **Endorphin secretion**

When a nervous signal proceeds from the pain area to the brain, it spreads through a chain of joined connections called synapses. The synapse can be seen as the space between the end of a nerve and the start of another. When an electric signal reaches the end of a nerve, it produces substances called neurotransmitters that pass through the synapse and activate the start of the following nerve. The process is repeated until the signal reaches the brain. The opioids involved in pain reduction have the task of sliding into the synapse space and impeding the propagation of neurotransmitters. In this way, a chemical block of pain signals is obtained. Endorphins are opioids that are naturally produced by the body to tackle pain and they can act in both the marrow and the brain as effective analgesics. Tens can increase the natural production of endorphins and therefore, they decrease the perception of pain.

## **Incontinence programs**

Incontinence programs can be executed using specific anal and vaginal probes, which can be purchased separately. The probes must be covered by CE MDD certification, in accordance with the directive 93/42/EEC, and must be bipolar and equipped with a 2 mm female connector for 2 mm male cables.

### **Warnings**

Since these applications have a medical intended use, they must be carried out under medical supervision.

### **Use**

To use the probes correctly, follow the instructions provided by the manufacturer or the physician.

### **Maintenance**

For cleaning, sterilization and disinfection refer to manufacturer instructions.

### **Suggestions**

In case of deterioration, replace the probe and do not use it again.

## **Iontophoresis**

Iontophoresis is a form of electrotherapy through which pharmacological substances are transmitted inside the tissues thanks to a unidirectional continuous electric current.

Iontophoresis is based on the ionic dissociation capacity of some medicinal substances, which have very low molecular weight, after they are dissolved in water.

It is crucial to know if the active part of the medicine, after being dissociated in a ionic form, has positive or negative charge, with the aim of placing it correctly according to the direction of the electric flux.

The ions of the medicinal substance are transmitted inside the organism through cutaneous areas that oppose a low resistance to the current, reaching the cellular membranes that are thereby modified electrically. The home user can use iontophoresis treatments only after consulting a specialist who will prescribe treatment medicines and indications.

Current intensity should be regulated so as to be barely perceptible. We suggest starting the treatment from minimum intensity and increasing it gradually.

Therapy medicament **MUST NEVER BE APPLIED DIRECTLY ON THE SKIN**, but always on the absorbing surface of the electrode that corresponds to the polarity of the medicament, while the absorbing surface of the other electrode will have to be dampened with slightly salted water in order to facilitate current circulation.

## SPORT PROGRAM LIST

Name	Duo Tens	Elite S2	Genesy S2
Demo		P 1	P 1
Motor point pen		P 2	P 2
Lower limb warm-up		P 3	P 3
Lower limb pre-competition warm-up		P 4	P 4
Lower limb maximum strength		P 5	P 5
Lower limb endurance strength		P 6	P 6
Lower limb explosive strength		P 7	P 7
Lower limb reactivity		P 8	P 8
Lower limb capillarization		P 9	P 9
Lower limb active recovery		P 10	P 10
Lower limb muscle cool down		P 11	P 11
Lower limb endurance		P 12	
Lower limb decontracting		P 13	P 12
Decontracting	P 1		
Abdominal maximum strength		P 14	
Maximum strength	P 2		
Abdominal endurance strength		P 15	
Abdominal active recovery		P 16	
Abdominal cool down		P 17	
Cool down	P 3		
Abdominal capillarization		P 18	
Abdominal decontracting		P 19	
Upper limb warm-up		P 20	
Upper limb pre-competition warm-up		P 21	
Upper limb maximum strength		P 22	
Upper limb endurance strength		P 23	
Upper limb explosive strength		P 24	
Upper limb reactivity		P 25	
Upper limb capillarization		P 26	

<b>Name</b>	<b>Duo Tens</b>	<b>Elite S2</b>	<b>Genesy S2</b>
Upper limb active recovery		P 27	
Upper limb cool down		P 28	
Upper limb endurance		P 29	
Upper limb decontracting		P 30	

## **BEAUTY PROGRAM LIST**

<b>Name</b>	<b>Duo Tens</b>	<b>Elite S2</b>	<b>Genesy S2</b>
Leg anti-cellulite lipolysis		P 31	P 13
Hip-gluteus anti-cellulite lipolysis		P 32	
Abdominal anti-cellulite lipolysis		P 33	
Upper limb anti-cellulite lipolysis		P 34	
Leg drainage		P 35	P 14
Hip-gluteus drainage		P 36	
Abdominal area drainage		P 37	
Upper limb drainage		P 38	
Leg firming		P 39	P 15
Hip-gluteus firming		P 40	
Abdominal area firming		P 41	
Upper limb firming		P 42	
Leg toning		P 43	P 16
Dorsal toning		P 44	
Hip-gluteus toning		P 45	
Abdominal area toning		P 46	
Upper limb toning		P 47	
Leg sculpting		P 48	P 17
Hip-gluteus sculpting		P 49	
Abdominal area sculpting		P 50	
Pectoral sculpting		P 51	
Sculpting	P 4		
Leg vascularization		P 52	

<b>Name</b>	<b>Duo Tens</b>	<b>Elite S2</b>	<b>Genesy S2</b>
Upper limb vascularization		P 53	
Leg lipostress		P 54	
Gluteus lipostress		P 55	
Hip lipostress		P 56	
Abdominal area lipostress		P 57	
Upper limb lipostress		P 58	
Leg mass building		P 59	
Upper limb mass building		P 60	
Abdom. post-pregnancy drainage		P 61	P 18
Abdom. post-pregnancy lipolysis		P 62	P 19
Abdom post-pregnancy toning		P 63	P 20
Man abdominal area definition		P 64	
Man pectoral definition		P 65	
Woman leg mass building		P 66	
Woman up. limb mass building		P 67	
Swollen upper limbs		P 68	
Swollen legs		P 69	
Leg toning massage		P 70	
Face beauty 1		P 71	
Face beauty 2		P 72	
Breast toning 1		P 73	
Breast toning 2		P 74	

## TENS PROGRAM LIST

<b>Name</b>	<b>Duo Tens</b>	<b>Elite S2</b>	<b>Genesy S2</b>
Conventional antalgic Tens *	P 5	P 75	P 21
Endorphinic Tens *	P 6	P 76	P 22
Rotator cuff tendinitis *	P 7	P 77	P 23
Muscle pain *	P 8	P 78	P 24
Knee osteoarthritis *	P 9	P 79	P 25
Post-surgery pain *	P 10	P 80	P 26
Muscle injury *	P 11	P 81	P 27
Cervical pain *	P 12	P 82	P 28
Trapezius pain *	P 13	P 83	P 29
Sciatica *	P 14	P 84	P 30
Lumbago *	P 15	P 85	P 31
Epicondylitis *	P 16	P 86	P 32
Scapulo-humeral peri arthritis *	P 17	P 87	P 33
Carpal tunnel *	P 18	P 88	P 34
Osteoarthritis *	P 19	P 89	P 35
Menstrual pain *	P 20	P 90	P 36

## SPORT-PREVENTION PROGRAM LIST

<b>Name</b>	<b>Duo Tens</b>	<b>Elite S2</b>	<b>Genesy S2</b>
Ankle prevention		P 91	P 37
Knee prevention		P 92	P 38
Quadriceps muscle prevention		P 93	P 39
Lower limb reathletization			P 40
Shoulder reathletization		P 94	P 41
Shoulder prevention		P 95	P 42
Elbow prevention		P 96	P 43
Upper limb reathletization			P 44
Ankle muscle tone recovery			P 45

## REHABILITATION PROGRAM LIST

Name	Duo Tens	Elite S2	Genesy S2
Vastus medialis reinforcement *			P 46
Swollen ankles *			P 47
Hemiplegia upper limbs *			P 48
Hemiplegia lower limbs *			P 49
Recovery after ACL surgery *			P 50
Functional recovery lower limbs *			P 51
Quadriceps atrophy *		P 97	P 52
Shoulder subluxation prevention *		P 98	P 53
Muscle spasms *			P 54
Multiple sclerosis - Leg flexor extensor mm. *			P 55
Multiple sclerosis - Ankle flexor extensor mm. *		P 99	P 56

## INCONTINENCE PROGRAM LIST

Name	Duo Tens	Elite S2	Genesy S2
Mixed incontinence *		P 100	P 57
Stress incontinence *			P 58
Urge incontinence *			P 59

## IONOPHORESIS

Name	Duo Tens	Elite S2	Genesy S2
Iontophoresis *			P 60

\* Because of the presence of clinical programs, this product is a medical device. Therefore it is certified by the Kiwa Cermet Italia S.p.A. Body n. 0476 according to the 93/42/EEC EU directive for medical devices. The certification covers clinical applications.

## GENERAL NOTES ON ELECTRODE POSITIONING

The correct electrode positioning and the correct choice of their size are critical aspects for the electrostimulation to be effective.

To choose the size of the electrodes and as for their positioning it is necessary to refer to the images at the end of this manual. This information is also available on our website [globuscorporation.com](http://globuscorporation.com).

**NOTE** For all the programs that cause an important muscle contraction (such as, for example, strength, hypertrophy, toning and firming programs ...) it is important to position the electrode on the muscle motor point, that is the most sensitive point to stimulation.

If the electrode is not positioned exactly on the motor point, the contraction could be small and/or annoying. In this case it is necessary for the positive electrode to be moved a few millimeters up to feel an effective and comfortable contraction.

### The body position during the stimulation

The body position during the electrostimulation session depends on the body part involved and on the program type that is being carried out. During the treatment execution with high intensities, we suggest blocking the limbs in order to work in isometry. For instance, if you want to treat the quadriceps with a strength program, we suggest carrying out the treatment while sitting with the foot blocked, in order to avoid an involuntary leg extension during the contraction phase.

For all the programs that do not imply high intensity (massages, decontracting, drainage programs) the body position is not important, as long as it is comfortable.

### Electrode positioning for Tens programs

In the following pages of this manual you can find the images with the correct electrode positioning for tens treatments. If the localization of your pain is not included in the images represented, you can position the electrodes by forming a "square" on the painful area. Here's an example.

## TREATMENT DESCRIPTIONS

Following is a short description of the programs provided by the electrostimulator.

For further information about the programs you can download from our website [globuscorporation.com](http://globuscorporation.com), a complete guide, where you can find all the indications and suggestions to carry out the treatment in the correct way.

### SPORT

**WARM-UP:** increases muscular temperature for better training preparation. Increases the blood flow to and metabolism of the muscle. Recommended before training sessions or competitions.

Intensity: low or intermediate.

**MAXIMUM STRENGTH:** improves maximal force capability and increases muscular mass. Allows training with maximum load and reduced trauma risk as compared to classical training. Recommended for activities requiring great physical strength. Also useful during the early preparation period as it trains the muscles to receive significant stress.

Intensity: maximum, at pain threshold limit. During active Rest, adjust the intensity to obtain comfortable contractions.

**ENDURANCE STRENGTH:** improves the ability to display a high level of force for an extended period of time and to resist toxin accumulation, thus delaying muscular fatigue. The program involves a long series of lengthy contractions, alternated with brief active rest. Recommended for sports requiring intense muscular work for an extended period of time. Intensity: intermediate or maximum. Active Rest should induce muscular vibrations.

**EXPLOSIVE STRENGTH:** the Explosive Force program increases the ability to rapidly achieve maximum levels of force. The objective is to recruit the maximum number of muscle fibers in the shortest time possible. It is recommended to begin with a Maximal Force program to prepare the muscle for the Explosive Force program. Recommended in all sports requiring maximal force under natural load conditions, such as sprints, ball sports, and alpine skiing. Intensity: maximum muscular contraction at pain threshold limit. Active Rest should induce muscular vibrations.

**DECONTRACTING:** enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programs can be used whenever muscular relaxation is needed. Intensity: low to intermediate.

**CAPILLARIZATION:** increases the blood flow to the muscle and improves muscle resistance and recovery. An increased arterial capacity keeps the primary and secondary capillary network active, thus improving the tissue oxygenation system.

This delays the fatigue during a very intense work. Capillarization programs are especially recommended for sports that require endurance strength and muscle endurance. Intensity: intermediate.

**POST-COMPETITION/TRAINING ACTIVE RECOVERY:** improves recovery from muscular fatigue after training or competitions and enables easier training resumption. The stimulation varies from muscular vibrations to light contractions. This allows greater circulation and reduces muscular hardening it also alleviates pain thanks to the release of endorphins. This program is valid for all sports as it improves and speeds up the recovery of maximum muscular capacity. This program is recommended in the hours immediately after training or competitions.

Intensity: low to intermediate.

### **BEAUTY**

**DRAINAGE:** the coordinated stimulation obtained by using all channels of the electrostimulator produces an in-depth muscular massage effect. The massage activates venous and lymphatic systems by means of a mechanical effect. The purpose is to improve blood circulation and enable the toxin elimination. Intensity: low to intermediate.

**LIPOLYSIS:** the low frequency of this program positively affects the metabolism and increases the elimination of toxins and excess fluid. This is a method for reshaping critical areas, modifying the metabolism of lipocytes. Improves blood circulation. It helps to eliminate the orange peel skin. Intensity: low to intermediate.

**VASCULARIZATION:** enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programs can be used whenever muscular relaxation is needed.

Intensity: intermediate.

**FIRMING:** this program is suitable for sedentary individuals with flaccid muscles that require reactivation of muscular tone. It prepares the muscles for more intense stimulations. It can be performed every day on large muscle groups such as thighs, abdominal and dorsal muscles. Intensity: intermediate to high.

**TONING:** the toning program allows to develop and maintain muscular tone. This type of treatment may be combined with aerobic activity and is suitable for users of any level. The program is especially directed at all those who, mindful of their own physical condition, wish to integrate their own regular training. The program can be used up to three times a week. Intensity: intermediate to high.

**SCULPTING:** the sculpting program is the last specifically aesthetic step and it is ideal for those who want to improve the definition of their already toned muscles. The work is very intense and so the day after the treatment muscles could feel a bit stiff. In order to prevent this, it is possible to make stretching exercise before and after the session or add a cool-down program.

The program is suited for those who want to define muscles in the areas in which traditional exercise movements for muscle development are not sufficient. Intensity: from intermediate to high. However, this effect may be avoided by performing stretching exercises before and after the session or a winding- down program. The program is suitable for those who wish to achieve muscular definition in those areas where traditional muscular development exercises are insufficient.

Intensity: intermediate to high.

**MASS BUILDING:** specific program for those who want to increase muscle volume and mass. Intensity: maximum muscular contraction, on the verge of pain threshold. It is recommended to execute a muscle warm-up program before carrying out any firming, toning, sculpting and mass-building program.

## TENS

Uses the mechanisms of the nervous system to relieve both acute and chronic pain. Pulses are spread through nerves in order to stop pain propagation. The analgesic effect is obtained while the program is executed and it can last for a few hours. Its duration varies. TENS programs also improve circulation. Place the electrodes on the aching area. It can be used any time pain relieve is needed. Warning! TENS programs relieve pain but do not eliminate its cause. If pain persists, consult the physician!

**TENS (Antalgic-Endorphinic):** relieves pain and relaxes the muscle effectively. The TENS program uses the mechanisms of the nervous system to relieve both acute and chronic pain. Pulses are spread through nerves in order to stop pain propagation. Muscle relaxation increases circulation and helps muscle relaxation. It can be used any time pain relieve and muscle relaxation are needed. Warning! TENS programs relieve pain but do not eliminate its cause. If pain persists, consult the physician! Intensity: ANTALGIC TENS: non-painful tingling sensation. It is not supposed to induce any muscle contraction. ENDORPHINIC TENS: visible muscle vibrations.

**TENS (Cervical pains, Lumbago, Periarthritis, Epicondylitis, Sciatica, ...):** uses the mechanisms of the nervous system to relieve both acute and chronic pain. Pulses are spread through nerves in order to stop pain propagation. This program alternates stimulations with slight contractions that stimulate endorphin release. Warning! TENS programs relieve pain but do not eliminate its cause. If pain persists, consult the physician! Intensity: non-painful tingling sensation.

## **REHAB**

### **VASTUS MEDIALIS REINFORCEMENT AND RECOVERY AFTER ACL SURGERY:**

after a trauma or another interruption of the physical activity, muscles often lose strength and mass. It is important to reestablish normal muscle capability as soon as possible. The program gives the possibility of training specific muscles. It can be used any time it is necessary to restart the activity. Intensity: from moderate to intermediate.

**PREVENTION PROGRAMS:** stimulation aims at facilitating neuromuscular work, helping the mobilization of the injured joint. It prepares the limb to restart its normal activity. We recommend using it in the first phases of muscle trophism recovery. Intensity: Intermediate.

**REATHLETIZATION:** increase in muscle fiber diameter and capacity. Use is suggested in the first phases of recovery after an injury/trauma, as soon as the muscle shows a little volume or tone recovery. Intensity: Intermediate.

**HEMIPLEGIA:** if the patient suffers from a perception problem with attention reduction (hemineglect), the training of this ability can be profitably integrated with the treatment. The treatment should not be only passive, but should also be used actively during the training. The program can be used to facilitate the reeducation of motor faculties and as a method to reduce spasticity.

Intensity: depending on the muscle that has to be stimulated.

**INCONTINENCE:** specific programs for the stimulation of perineal muscles. It is recommended in case of urological problems and incontinence. The programs must be used with legal probes.

## WARRANTY

The device is guaranteed to the first user for a period of twenty-four (24) months from the date of purchase against defects in materials or of the manufacturing (if the user uses the device for professional purposes, is twelve, 12, months), provided that it is used properly and maintained under normal operating conditions.

Warranty coverage is limited in the following cases:

- Six (6) months for supplied accessories such as batteries, chargers, power supplies, cables, etc.;
- No warranty for extendable accessories and materials such as, for example, electrodes, etc..

The warranty is valid and enforceable in the country where the product was purchased. In the event that the product is purchased in any country of the European Community, the warranty is still valid in all its countries.

In order to take advantage of the warranty service, the user must comply with the following warranty clauses:

1. The products, and all accessories, to be repaired must be sent by and at the expenses of the customer in their original packages.
2. The product warranty is subject to the production of a fiscal document (fiscal receipt, receipted bill or sale invoice), attesting the product purchase date.
3. The repair work shall have no effect on the original warranty expiry date and shall neither renew nor extend it.
4. If, during the repair work, no defects are found, the costs related to inspection times shall in any case be charged.
5. The warranty becomes void where the fault has been caused by: impacts, falls, erroneous or improper use of the product, use of non-original power supply unit or charger, accidental events, alteration, replacement/detachment of the warranty seals and/or tampering with the product. The warranty does not cover damages caused during transportation when unsuitable packages are used.
6. The warranty does not cover the inability to use the product, other incidental or consequent costs or other expenses incurred by the purchaser.

The features and dimensions reported in this manual are not binding.

## FREQUENTLY ASKED QUESTIONS

**What happened if, when I increase intensity above 1 mA, it goes immediately back to 0 mA?**

Make sure that cables are connected, that electrodes are placed on the area to be treated and that they are not too worn out, as this could affect their conductivity.

**What electrodes should be used?**

It is recommended to use adhesive electrodes, as they improve stimulation quality and are more practical. If used with care, i.e. on clean skin, they can be used for up to 25-30 applications. The electrodes must always be replaced if they do not stay perfectly in contact with the skin.

**Where should I place the electrodes?**

This manual reports the images of electrode positioning for all body parts. (Respecting the indicated polarity is not necessary) It is sufficient to follow these instructions.

To verify the correct placement of the electrodes, use the special Motor Point Pen or follow this empirical method: place the electrodes as indicated in the pictures; start the stimulation; with your hand, slide the electrode in different directions on the muscle area. You will notice an increase or a decrease of the stimulation depending on the electrode position. Once you locate the point of highest stimulation, decrease channel intensity to zero (0.0 mA), place the electrode again and increase intensity gradually.

**Do Y cables permit to use more electrodes on the same channel?**

This allows you to work, for example, on quadriceps vastus medialis and vastus lateralis with the same channel; you can use both split cables to treat two limbs simultaneously, involving 4 muscles. Do not use them for medical applications.

**Does power decrease using Y cables?**

The power intensity of each channel does not vary. However, when Y cables are used to split one channel, the current is distributed on a wider muscle area, therefore the contraction will be weaker. Increase the intensity to obtain the same contraction.

**Can electrostimulation hurt?**

Muscle damage is very unlikely. An important principle that should be followed is to increase intensity gradually, observing muscle behaviour, and avoid to keep the limb completely outstretched. In case of doubt, please contact a specialist.

**Is it possible to use the electrostimulator during menstrual cycle?**

Interferences like anticipation, delay, increase or reduction of the cycle are extremely subjective and variable. In any case, it is recommended to avoid treatments in the abdominal area during menstrual cycle and immediately before or after it.

### **Is it possible to use the electrostimulator during lactation?**

Until now, the use of electrostimulation during lactation has not shown any side effects. Nonetheless, we suggest not to treat the thoracic region.

### **Are dermatological diseases (e.g. psoriasis, urticaria) contraindications for electrostimulation?**

Yes, do not treat skin areas affected by serious dermatological diseases.

### **When are the first results visible?**

The aesthetic results of electrostimulation are always subjective.

However, it is possible to say that, in the case of tonification, 3-4 regular and constant weekly sessions can already bring good results after 15 days, while electrolipolysis and electrodrainage require 40 days. It is possible to obtain better and quicker results if treatments are combined with good physical activity and a correct life style.

### **How many sessions can be performed on a weekly basis?**

As regards fitness and beauty, the number of sessions depends on the type of treatment. 3-4 sessions per week on alternate days are suggested for toning. Daily treatments are permitted for lipolysis and drainage programs.

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