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**MANUALE D'USO**

**USER MANUAL**

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

**MANUAL DE USO**

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**MY•STIM**

*elite*

0197

## EGREGIO CLIENTE

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# SUMMARY

SUMMARY	29
PROGRAM LIST Mod. ELITE	30
PROGRAM LIST Mod. MY STIM	36
ACCESSORIES	40
CONTRAINDICATIONS	41
CARE INSTRUCTIONS	42
- Taking care and cleaning the stimulator	
- Maintenance and use of electrodes and cables	
WARNINGS	43
- Precautionary measures	
- Use destination	
TECHNICAL FEATURES	44
BATTERY RECHARGEMENT	45
- How to charge the battery	
STIMULATOR AND DISPLAY DESCRIPTION	46
KEYBOARD DESCRIPTION	47
USING THE STIMULATOR	48
STIMULATION INTENSITY	51
- Open Circuit	
- Stimulation intensity	
TREATMENT DESCRIPTIONI SPORT, FITNESS AND BEAUTY, HEALTH AND REHAB	52
USE OF THE MOTOR POINTS FINDER	56
COMMON QUESTIONS	57
WARRANTY TERMS	58
ELECTRODE PLACEMENT	149

## Mod. ELITE - PROGRAM LIST SPORT MAN

Muscle	Program number	Program name	Electrode placement	Phase time (min.) Intensity
	Sport M. <b>1</b>	Capillarisation	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport M. <b>2</b>	Warm Up	As indicated in the figures at the end of the manual	10' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport M. <b>3</b>	Recovery post-comp./training	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport M. <b>4</b>	Maximum force lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>5</b>	Resistent force lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>6</b>	Explosive force lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>7</b>	Reactivity lower loimbs	I1-I2-I3-I4-I5-I6	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>8</b>	Relaxation lower limbs	I1-I2-I3-I4-I5-I6	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport M. <b>9</b>	Maximum force abdominal/dorsal	T1-T2-T3-T4-T5-T6	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>10</b>	Resistent force abdominals/dorsals	T1-T2-T3-T4-T5-T6	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>11</b>	Relaxation abdominals/dorsals	T1-T2-T3-T4-T5-T6	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport M. <b>12</b>	Maximum force upper limb/pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>13</b>	Resistent force upper limbs/pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>14</b>	Explosive force upper limbs/pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport M. <b>15</b>	Relaxation upper limbs pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** =Intermediate - **Max** = Maximum

## Mod. ELITE - PROGRAM LIST SPORT WOMAN

Muscle	Program number	Program name	Electrode placement	Phase time (min.) Intensity
	Sport W. <b>16</b>	Capillarisation	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport W. <b>17</b>	Warm up	As indicated in the figures at the end of the manual	10' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport W. <b>18</b>	Recovery post-comp./training	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport W. <b>19</b>	Maximum force lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>20</b>	Resistent force lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>21</b>	Explosive force lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>22</b>	Reactivity lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>23</b>	Relaxation lower limbs	I1-I2-I3-I4-I5-I6	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport W. <b>24</b>	Maximum force abdominals / dorsals	T1-T2-T3-T4-T5-T6	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>25</b>	Resistent force abdominals / dorsals	T1-T2-T3-T4-T5-T6	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>26</b>	Relaxation abdominals / dorsals	T1-T2-T3-T4-T5-T6	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Sport W. <b>27</b>	Maximum force lower limbs/pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>28</b>	Resistent force lower limbs/pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 20' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>29</b>	Explosive force lower limbs/pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' Int Max Int <input type="checkbox"/>
	Sport W. <b>30</b>	Relaxation lower limbs/pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** =Intermediate - **Max** = Maximum

## Mod. ELITE - PROGRAM LIST HEALTH and REHAB

Muscle	Program number	Program name	Electrode placement	Phase time (min.) Intensity
	Rehab <b>31</b>	DEMO	As indicated in the figures at the end of the manual	2' 2' 2' 2' Mod Mod Mod Mod
	Rehab <b>32</b>	Search motor points	As indicated in the figures at the end of the manual	5' Mod <input type="text"/> <input type="text"/> <input type="text"/>
	Rehab <b>33</b>	Antalgic Tens	As indicated in the figures at the end of the manual	30' Mod <input type="text"/> <input type="text"/> <input type="text"/>
	Rehab <b>34</b>	Endorphinic Tens	As indicated in the figures at the end of the manual	20' 20' Mod Mod <input type="text"/> <input type="text"/>
	Rehab <b>35</b>	Contractures	As indicated in the figures at the end of the manual	30' Mod <input type="text"/> <input type="text"/> <input type="text"/>
	Rehab <b>36</b>	Haematomas	As indicated in the figures at the end of the manual	20' Mod <input type="text"/> <input type="text"/> <input type="text"/>
	Rehab <b>37</b>	Sciatica	R7	30' Mod <input type="text"/> <input type="text"/> <input type="text"/>
	Rehab <b>38</b>	Atrophy recovery lower limbs	I1-I2-I4	5' 15' Int Ele <input type="text"/> <input type="text"/>
	Rehab <b>39</b>	Lumbago	T3	20' 20' Mod Mod <input type="text"/> <input type="text"/>
	Rehab <b>40</b>	Atrophy recovery upper limbs	S1-S2-S3-S4-S5-S6-S7	5' 15' Int Int <input type="text"/> <input type="text"/>
	Rehab <b>41</b>	Cervicals	R1	30' Mod <input type="text"/> <input type="text"/> <input type="text"/>
	Rehab <b>42</b>	Peryarthrititis	R2	15' 10' Mod Mod <input type="text"/> <input type="text"/>
	Rehab <b>43</b>	Epicondylitis	R3-R4	15' 10' Mod Mod <input type="text"/> <input type="text"/>

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** =Intermediate - **Max** = Maximum

# Mod. ELITE - PROGRAM LIST FITNESS MAN

Muscle	Program number	Program name	Electrode placement	Phase time (min.) Intensity
	Fitness M. <b>44</b>	Relaxant massage <b>Bio Pulse</b>	As indicated in the figures at the end of the manual	25' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>45</b>	Toning Massage	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>46</b>	vascularisation lower limbs	I1-I2-I3-I4-I5-I6	25' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>47</b>	Dtainage lower limbs	I1-I2-I3-I4-I5-I6	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>48</b>	Lypolisis lower limbs	I1-I2-I3-I4-I5-I6	30' 20' Int Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>49</b>	Firming lower limbs	I1-I2-I3-I4-I5-I6	20' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>50</b>	Toning lower limbs	I1-I2-I3-I4-I5-I6	15' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>51</b>	Sculpting / Firming lower limbs	I1-I2-I3-I4-I5-I6	20' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>52</b>	Mass lower limbs	I1-I2-I3-I4-I5-I6	15' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
		Fitness M. <b>53</b>	Drainage abdominal / dorsals	T1-T2-T3-T4-T5-T6
Fitness M. <b>54</b>		Lypolisis abdominals / dorsals	T1-T2-T3-T4-T5-T6	30' 20' Int Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>55</b>		Firming abdominals / dorsals	T1-T2-T3-T4-T5-T6	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>56</b>		Toning abdominals / dorsals	T1-T2-T3-T4-T5-T6	15' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>57</b>		Sculpting abdominals / dorsals	T1-T2-T3-T4-T5-T6	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>58</b>	Drainage Upp. limbs pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>59</b>	Lypolisis Upp. limbs pectorals	S1-S2-S3-S4-S5-S6-S7-S8	30' 20' Int Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>60</b>	Firming Upp. limbs pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>61</b>	Toning Upp. limbs pectorals	S1-S2-S3-S4-S5-S6-S7-S8	15' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>62</b>	Sculpting Upp. limbs pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness M. <b>63</b>	Mass Upp. limbs pectorals	S1-S2-S3-S4-S5-S6-S7-S8	15' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>


Key words: Mod = Moderated - Elev = Elevated - Int =Intermediate - Max = Maximum

## Mod. ELITE - PROGRAM LIST FITNESS WOMAN

Muscle	Program number	Program name	Electrode placement	Phase time (min.) Intensity
	Fitness W. <b>64</b>	Relaxant Massage <b>Bio Pulse</b>	As indicated in the figures at the end of the manual	25' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>65</b>	Toning massage	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>66</b>	vascularisation	11-12-13-14-15	25' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>67</b>	Drainage thighs and hips 1	11-12-13-14-15	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>68</b>	Drenaggio thighs and hips <b>Bio Pulse</b>	11-12-13-14-15	20' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>69</b>	Lypolisi thighs and hips 1	11-12-13-14-15	30' 20' Int Int <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>70</b>	Lipolisi thighs and waist 2 <b>Bio Pulse</b>	11-12-13-14-15	20' 20' Int Int <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>71</b>	Firming thighs and hips 1	11-12-13-14-15	20' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>72</b>	Firming thighs and hips 1 <b>Bio Pulse</b>	11-12-13-14-15	20' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>73</b>	Toning thighs and hips	11-12-13-14-15	15' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>74</b>	Sculpting thighs and hips	11-12-13-14-15	20' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>75</b>	Drainage gluteus 1	16	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>76</b>	Drainage gluteus 2 <b>Bio Pulse</b>	16	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>77</b>	Lypolisi gluteus 1	16	30' 20' Int Int <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>78</b>	Lypolisi gluteus 2 <b>Bio Pulse</b>	16	20' 20' Int Int <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>79</b>	Firming gluteus 1	16	20' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
	Fitness W. <b>80</b>	Firming gluteus 2 <b>Bio Pulse</b>	16	20' Ele <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** = Intermediate - **Max** = Maximum

## Mod. ELITE - PROGRAM LIST FITNESS WOMAN

Muscle	Program number	Program name	Electrode placement	Phase time (min.) Intensity
	Fitness W. <b>81</b>	Toning gluteus 1	I6	15' 
	Fitness W. <b>82</b>	Toning gluteus 2	I6	15' 
	Fitness W. <b>83</b>	Sculpting gluteus 1	I6	20' 
	Fitness W. <b>84</b>	Sculpting gluteus 2	I6	20' 
	Fitness W. <b>85</b>	Drainage abdominals 1	T5-T6	20' 
	Fitness W. <b>86</b>	Drainage abdominals 2 <b>Bio Pulse</b>	T5-T6	20' 
	Fitness W. <b>87</b>	Lypolisis abdominals 1	T5-T6	30' 
	Fitness W. <b>88</b>	Lypolisis abdominals 2 <b>Bio Pulse</b>	T5-T6	20' 
	Fitness W. <b>89</b>	Firming abdominals 1	T5-T6	20' 
	Fitness W. <b>90</b>	Firming abdominals 2 <b>Bio Pulse</b>	T5-T6	20' 
	Fitness W. <b>91</b>	Toning abdominals	T5-T6	15' 
	Fitness W. <b>92</b>	Sculpting abdominals	T5-T6	20' 
	Fitness W. <b>93</b>	Firming breast	S9	15' 
	Fitness W. <b>94</b>	Sculpting breast	S9	10' 
	Fitness W. <b>95</b>	Drainage upper limbs	S1-S2-S3-S4-S5-S6-S7	20' 
	Fitness W. <b>96</b>	Lypolisis upper limbs	T5-T6	30' 
	Fitness W. <b>97</b>	Capillarisation face	V1	15' 
	Fitness W. <b>98</b>	Face lifting effect	V1	15' 

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** = Intermediate - **Max** = Maximum

# Mod. MY STIM PROGRAM LIST SPORT

N° Program	Program number	Electrode placement	Phase time (min.) Intensity
<b>1</b>	DEMO	As indicated in the figures at the end of the manual	2' 2' 2' 2' Mod Mod Mod Mod
<b>2</b>	Search motor points	As indicated in the figures at the end of the manual	2' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>3</b>	Capillarisation	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>4</b>	Warm Up	As indicated in the figures at the end of the manual	10' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>5</b>	Tension relieving	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>6</b>	Post competition/ workout recovery	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>7</b>	Maximal power - lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' <input type="checkbox"/> Int Max Int
<b>8</b>	Resistance - lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' <input type="checkbox"/> Int Max Int
<b>9</b>	Explosive power - lower limbs	I1-I2-I3-I4-I5-I6	5' 12' 5' <input type="checkbox"/> Int Max Int
<b>10</b>	Maximal power - torso and upper limbs	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' <input type="checkbox"/> Int Max Int
<b>11</b>	Resistance - torso and upper limbs	S1-S2-S3-S4-S5-S6-S7-S8	5' 12' 5' <input type="checkbox"/> Int Max Int
<b>12</b>	Explosive power - torso and upper limbs	S1-S2-S3-S4-S5-S6-S7-S8	5' 12' 5' <input type="checkbox"/> Int Max Int

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** = Intermediate - **Max** = Maximum

## Mod. MY STIM - PROGRAM LIST FITNESS

N° Program	Program number	Electrode placement	Phase time (min.) Intensity
<b>13</b>	Relaxing massage	As indicated in the figures at the end of the manual	25' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>14</b>	Toning massage	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>15</b>	Vascularisation - lower limbs	I1-I2-I3-I4-I5-I6	25' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>16</b>	Drainage - lower limbs	I1-I2-I3-I4-I5-I6	30' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>17</b>	Lipolysis - lower limbs	I1-I2-I3-I4-I5-I6	40' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>18</b>	Firming - lower limbs	I1-I2-I3-I4-I5-I6	20' 10' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>19</b>	Toning - lower limbs	I1-I2-I3-I4-I5-I6	20' 10' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>20</b>	Modelling - sculpting - lower limbs	I1-I2-I3-I4-I5-I6	20' 20' Max Int <input type="checkbox"/> <input type="checkbox"/>
<b>21</b>	Increasing mass - lower limbs	I1-I2-I3-I4-I5-I6	15' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>22</b>	Drainage - abdominal belt	T1-T2-T3-T4-T5-T6	30' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>23</b>	Lipolysis - abdominal belt	T1-T2-T3-T4-T5-T6	40' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>24</b>	Firming - abdominal belt	T1-T2-T3-T4-T5-T6	20' 10' Int Int <input type="checkbox"/> <input type="checkbox"/>
<b>25</b>	Toning - abdominal belt	T1-T2-T3-T4-T5-T6	20' 10' Int Int <input type="checkbox"/> <input type="checkbox"/>
<b>26</b>	Modelling - abdominal belt	T1-T2-T3-T4-T5-T6	20' 20' Int Int <input type="checkbox"/> <input type="checkbox"/>
<b>27</b>	Drainage - upper limbs and pecs	S1-S2-S3-S4-S5-S6-S7-S8	30' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>28</b>	Lipolysis - upper limbs and pecs	S1-S2-S3-S4-S5-S6-S7-S8	40' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>29</b>	Firming - upper limbs and pecs	S1-S2-S3-S4-S5-S6-S7-S8	20' 10' Int Int <input type="checkbox"/> <input type="checkbox"/>
<b>30</b>	Toning - upper limbs and pecs	S1-S2-S3-S4-S5-S6-S7-S8	20' 10' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>31</b>	Modelling - upper limbs and pecs	S1-S2-S3-S4-S5-S6-S7-S8	20' 20' Max Int <input type="checkbox"/> <input type="checkbox"/>
<b>32</b>	Increasing mass - upper limbs and pecs	S1-S2-S3-S4-S5-S6-S7-S8	15' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** =Intermediate - **Max** = Maximum

## Mod. MY STIM - PROGRAM LIST *FITNESS*

N° Program	Program number	Electrode placement	Phase time (min.) Intensity
<b>33</b>	Drainage - thighs and hips	I1-I2-I3-I4-I5	30' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>34</b>	Lipolysis - thighs and hips	I1-I2-I3-I4-I5	40' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>35</b>	Firming - thighs and hips	I1-I2-I3-I4-I5	20' 10' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>36</b>	Toning - thighs and hips	I1-I2-I3-I4-I5	20' 10' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>37</b>	Drainage - gluteus	I6	30' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>38</b>	Lipolysis - gluteus	I6	40' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>39</b>	Firming - gluteus	I6	20' 10' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>40</b>	Toning - gluteus	I6	20' 10' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>41</b>	Modelling - gluteus	I6	20' 20' Ele Int <input type="checkbox"/> <input type="checkbox"/>
<b>42</b>	Firming - breast	S9	15' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>43</b>	Modelling - breast	S9	10' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>44</b>	Capillarisation - face	V1	15' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>45</b>	Face lifting effect	V1	15' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** =Intermediate - **Max** = Maximum

## Mod. MY STIM - PROGRAM LIST *HEALTH and REHAB*

N° Program	Program number	Electrode placement	Phase time (min.) Intensity
<b>46</b>	Analgesic effect by tens for acute pain	As indicated in the figures at the end of the manual	30' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>47</b>	Analgesic effect by tens for chronic pain	As indicated in the figures at the end of the manual	30' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>48</b>	Endorphin effect by tens	As indicated in the figures at the end of the manual	30' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>49</b>	Hematoma	As indicated in the figures at the end of the manual	20' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>50</b>	Sciatica	R7	30' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>51</b>	Recovery from atrophy - lower limbs	I1-I2-I4	5' Int 15' Ele <input type="checkbox"/> <input type="checkbox"/>
<b>52</b>	Recovery from atrophy - upper limbs	S1-S2-S3-S4-S5-S6-S7-S8	5' Int 15' Ele <input type="checkbox"/> <input type="checkbox"/>
<b>53</b>	Backache	T3	20' Mod 20' Mod <input type="checkbox"/> <input type="checkbox"/>
<b>54</b>	Cervical pain	R1	30' Mod <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>55</b>	Periarthritis	R2	15' Mod 10' Mod <input type="checkbox"/> <input type="checkbox"/>

Key words: **Mod** = Moderated - **Elev** = Elevated - **Int** = Intermediate - **Max** = Maximum

# ACCESSORIES

Accessories:

- a. 4 electrode connection cables  
(coloured cables for EMS, TENS treatments)
- b. Case of 4 reusable adhesive electrodes (5 x 5 cm)
- c. Case of 4 reusable adhesive electrodes (5 x 9 cm)
- d. Bag for transport
- e. Battery charger (see technical features)
- f. Operating manual complete of electrode placement guide
- g. Unit

This stimulator is supplied complete of cables and electrodes. When you open the packing, verify the simulator equipment is complete. If some of the above- mentioned accessories are missing, you should immediately contact the shopkeeper. Also check the stimulator and electrodes integrity.

**Electrodes:** The electrodes are for single use they are self-adhesive and don't require gel. The electrode cable is female. They are branded CE complied with MDD 93/42/EEC requirements. All the information can be subjected to changes without notice.

## STANDARD SYMBOLS

It refers to the manufacturer

Warning

0197

This symbol on your device indicates that it complies with the directives on medical devices (93/42/CEE 47/2007CEE). The number of the notified unit is 0197.

It indicates that this is a II class device.

It indicates that the product has been produced respecting the directive 2002/95/CE.

It informs the operator of a compulsory conduct

It refers to the pressure of the storage and transport environment where the device and its accessories are used

It refers to the humidity of the storage and transport environment where the device and its accessories are used.

WEEE symbol (Waste of Electrical and Electronic Equipment). Recycling symbol. The WEEE symbol used for this product indicates that the device may not be disposed of as a household product. Properly dispose of the product to help protect the environment.

It indicates the optimal temperatures for the storage and transportation of the product.

# CONTRAINDICATIONS

## Contraindications

- Do not stimulate the neck region (carotid sinus).
  - Patients who have a demand type cardiac pacemaker.
  - Stimulation should be not used on cancer patients (previous permission received by the oncologist).
  - Do not stimulate the cerebral zone.
  - Do not use on undiagnosed pains until aetiology is established.
  - Do not use on sores and dermatological diseases.
  - Do not use on acute traumas.
  - Do not stimulate recent cicatrix.
  - Do not use during pregnancy.
  - Do not use during driving car.
  - It is absolutely forbidden to stimulate the ocular and genital zone
  - It is also forbidden the use of the stimulator in the water (it is not protected against the water).
- Do not use the stimulator in environment where combustible gas, inflammables, or oxygen are present.
- It is also properly to use with caution in case of capillary fragility because an excessive stimulation could cause the break of more quantity of capillaries. Anyway, it is always preferable to contact your doctor before using.

## Adverse possible effects

Skin irritation can occur in rare case. In case of allergic reactions by the use of gel, stop the treatment and contact the specialist doctor.

If during the stimulation treatment tachycardia or extrasistolia should appear, it is better to stop the treatment. The sensation will pass immediately.

# CARE INSTRUCTION

## Taking care and cleaning the stimulator

- In case of damage of the covering apparatus, this has to be replaced.
- In case of real or presumed spoilt, do not make any repairs by yourself. Only authorized centres can carry out repairs.
- Do not give a shock to the product by clashing, overturning or shaking. It may cause trouble or accident in a long run although no immediate trouble occurs.
- Store the electro stimulator at a place free from a wet areas and water. Avoid installing in places subject to humidity, high temperature, dust or places exposed to direct sunbeam or near to sources of high electromagnetic emission: The stimulator must be used free and not wrapped or covered.
- Do not use thinner, gasoline, kerosene, cleaning powder or chemicals; for cleaning the machine and the cables, it is sufficient a soft duster lightly humid.
- After a long period of no use, it is advisable to check the correct working of the stimulator and of the battery. (Contact authorized centres).
- Do not use the apparatus during the battery rechargement.
- The stimulator is connected to the electricity by means of the apposite feeder. Before the connection check that electric system is complied with the requirements in force in your country.
- It is advisable not to use the apparatus together with other electronic equipment above all those used for support of vital functions.

## Maintenance and use of electrodes and cables

Before placing the electrodes on your body, be sure to always wash and dry the skin in order to ensure good conductivity.

After the use of the reusable mono-user electrodes, store them in the apposite seat and insert in a plastic bag cleaned and closed.

Check that there is no contact or overlapping among the electrodes.

The electrodes can be used for 25/30 applications.

If the electrodes loose their conductivity and adhesion, they could be old and need replacement. In case of use of electrodes not adhesive, it is advisable to clean the surface with water.

Use only the electrodes recommended by the supplier.

Be sure that the cables and electrodes are cleaned before any application.

In case of cable deterioration, replace them with original spare parts.

# WARNINGS

## Precautionary measures

- Only person with mental competence must use the device.
- The electro stimulator must be used only for transcutaneous neuro-stimulation.
- The neurostimulation apparatus should only be used under the continued supervision of a physician.
- The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation.
- The user must keep the device out of reach of children.
- Electronic monitoring equipment (such as ECG monitors and ECG alarms) may not operate properly when stimulation is in use
- Stimulus delivered by this device may be sufficient to cause electrocution: Electrical current of this magnitude must not flow through the thorax because it may cause a cardiac arrhythmia. (Do not stimulate pectoral and dorsal muscles simultaneously)
- Do not use the apparatus in case of any health problem without consulting your doctor.
- Do not use the stimulator near short wave devices.

## Use destination

Read the operating manual carefully in order to avoid carrying out installation, use and maintenance non-adequate that could damage the apparatus, the operator and the user.

The electrostimulators are designed to be used in the following operating environments:

- domestic environment;
- clinics;
- physiotherapy centers;
- rehabilitation centers;
- general pain treatments;
- beauty and sport purposes.

The use of this device is permitted to the patient (appropriately informed about the use conditions of the device) and to medical staff.

# TECHNICAL FEATURES

## TECHNICAL FEATURES

Number of channels:	4
Constant current:	Yes
Intensity:	0-100 mA with load 1000 Ohm
Wave form:	rectangular, biphasic, symmetric, compensated
Work frequency:	0,3-150 Hz
Recovery frequency:	0,3-150 Hz
Amplitude:	40-450 $\mu$ s.
Working time:	from 1 to 30 sec.
Recovery time:	from 0 to 1 min.
Frequency Range:	continuous variation from 1 to 150 Hz modulation minimum time 3 sec.
Amplitude Range:	continuous variation from 40 to 450 $\mu$ s. minimum frequency 1 Hz maximum frequency 150 Hz. modulation minimum time 5 sec.
Ramp-up time:	from 0 to 15 sec.
Ramp-down time:	from 0 to 15 sec.
Power supply:	accumulator Ni-MH rechargeable 7,2 V 1,5 Ah
Adaptor:	Model: Dk7-088-0200 Input: 100-240V 50-60Hz 0,2A 8.8V - 200mA  Output: 120mA (1000 ohm)
Devices case:	ABS
Protection:	IP 20
Dimension:	160x99x35,4
Weight:	404 gr.
Storage and transportation:	from -25°C to 70°C
Max. relative humidity:	0-93%
Atmospheric pressure:	from 700 hPa to 1060 hPa
Use conditions	
Temperature of use:	from 0°C to 40°CMax.
relative humidity:	from 30% to 75%
Atmospheric pressure:	from 700 hPa to 1060 hPa

## **Device removal**

Don't throw the device out but recycle or get rid of them appropriately following the removal rules in force in your country.

# BATTERY RECHARGEMENT

## How to charge the battery

The stimulator is supplied together with a rechargeable battery. The batteries are Nickel Metal hydrate type with peak efficiency without memory effect. When the battery indicator lights up on the display it is advisable to recharge the batteries.

Shut off the stimulator and the electrodes; connect it to the supplied charger (Model MED), inserting the plug in the special outlet (see the side figure).

The LED on the charger is lit when the connection is correctly executed. The first recharge must be absolutely complete and it has not to be stopped for any reason.

A complete recharge takes about 10 hours; should the stimulator accidentally stay connected for longer than (72 hours) the battery will not be damaged.

Only use the charger supplied with the stimulator. Battery should only be replaced by an authorized centre.

## Attention

Do not use a charger different from that supplied with the stimulator.



During the charger the stimulator **MUST NOT** be connected to the electrodes.

If the charger is connected with the stimulator switched on, it will automatically switch off.

Should the cable or connector show signs of wear and tear you should replace it immediately.

# STIMULATOR AND DISPLAY DESCRIPTION


## Warnings and information for first use

Before using the stimulator read the following information:

- Read this operating manual with attention;
- Examine with care the stimulator and its parts for individuating eventual external damages;
- if at the switching on all the segments should appear on the display, shut off the stimulator and start again. In case the segments would be present again, contact the assistance centre without using the machine.
- The stimulator is ready to use when on the display the software version (left low number).
- if the battery is not sufficiently charged, the machine will switch off immediately. It is therefore advisable to recharge the battery following the instruction mentioned in paragraph "Battery Rechargement".

# KEYBOARD DESCRIPTION

**NOTE :** 3" = By pressing the switch for three seconds you activate the funcion.

<b>ON/OFF/OK</b>	Confirms the selection of a program in progress and activates the pause 3" = Switching on/switching off.
<b>LEFT</b>	3" = during the program returns to the previous phase.
<b>UP</b>	During the program it increases the four channels intensity simultaneously.
<b>DOWN</b>	During the program it decreases the four channels intensity simultaneously.
<b>RIGHT</b>	3" During the program it allows to skip to the subsequent phase.
<b>TIME</b>	During the program displays the remaining total time or the remaining time of the ongoing phase.
<b>fn</b>	Pressed alone it allows the change of the duration phase.
<b>Intensity</b>	Increases/decreases the intensity of corresponding channels.

# USING STIMULATION

## PRACTICAL USE

For a correct use, proceed as follow:

- check that the stimulator is not connected to the electric system for the battery rechargement;
- connect the cables to the special output;
- attach the electrodes to the cables already connected to the stimulator;
- place the electrodes on the skin;

## Start up and navigation

To turn on press the button On/Off (OK) for three seconds until a bip is heard

Low left of the display, The software version and the model will be displayed for three seconds

### **How to select the program.**

By buttons Up e Down select the number of the program. Press OK to confirm. During the choice the following parameters are displayed:  
number of phases the program name (up) and the program duration (left up).

### **Starting program.**

After the program selection, use switches to increase the intensity of each channel or switch to increase the intensity on all channels simultaneously.

### **Display during the treatment.**

During the treatment execution the display will show you the following parameters: the number of program (top left), how many phases the program consists of, in which phase you are (the current square is filled), the remaining time of the ongoing phase (top right).

### **How to put the program in pause.**

You can put the program in Pause mode by pressing switch OK of Joypad. The intensity turns to 0. Start stimulation again by pressing OK.

The program continues where stimulation was shut off.

During the pause the lightening PAUSE description will appear.

### **How to stop the program**

To stop the program hold down the bottom OK for 3 seconds.

### **How to display the remaining time of a program.**

During the program, the remaining time of the ongoing phase is displayed. If you want to see the remaining total time press switch TIME.

It will be displayed for two seconds and all the phase indicators are filled.

### **How to change the phase time.**

To change the phase time, press switch Function fn; the time indicator of the phase in progress, starts blinking; press switches Up and Down of Joypad to change it.

The new entered time will be automatically confirmed after two seconds of visualisation.

### **Jumping phase.**

While stimulating, you have the option to jump to the next or previous sequence. Hold down the push button right or left for 3 seconds.

### **Increasing and decreasing the intensity**

To increase or decrease the intensity on each channel, press the corresponding switches Up or down.

If you want to increase or decrease all the channels simultaneously, press switches Up or Down of Joypad.

# STIMULATION INTENSITY

## Open circuit

Should the current circuit be open (due to a cable breakage or electrodes not correctly placed) and the user increases the intensity more than 10 mA, the intensity automatically shuts off to 0 mA.

It is therefore necessary before starting any treatments, to be sure that the cables are connected to the device and the electrodes have a good contact with your skin.


Also check if the electrodes are too old and need replacement.

## Stimulation intensity

The needed intensity value to obtain a special contraction, is completely subjective. It can depend on a series of factors such as the position of electrode, the adipose stratum, the perspiration etc. For these reasons the same current intensity value can give different sensations from a person to another, from a day to another or furthermore you could feel more contraction on the left limb rather than on the right one.

During the same working session, it will be necessary to set the intensity as to obtain the same contraction level. This is due to the so-called accommodation phenomenon

- Moderate: the muscle is not tired; the contraction is tolerable and comfortable. First level in the intensity graph.
- Intermediate: the muscle visibly contracts but it does not cause any articular movement. Second level in the intensity graph
- Elevated: the muscle sensibly contracts. The contraction would cause a articular movement for this reason the limb must be blocked. Third level in the intensity graph
- Maximum: the muscle is subjected to a maximum contraction. It is a challenging work that should be carried out after several applications.

PS. In the program guide also the advisable intensity values are indicated. We propose you an intensity value just as indication; you have to regulate it on the bases of your sensations.

# TREATMENTS

Here below you can find a brief description of the program included in the electro stimulator

## SPORT

### **Warm Up:**

Increases the muscle's temperature to make it more prepared for training. It also increases the blood supply and the metabolism in the muscle. Recommended before physical training or competition.

**Intensity:** moderate or intermediate.

### **Maximum Force:**

Improves the maximum force capacity and increases muscle mass. It allows to train with maximum loads with fewer risks of injuries compared to classic heavy strength training. Recommended for all sports where you need maximum strength.

**Intensity:** Maximum muscle contraction without pain. Active rest should give comfortable muscle vibrations.

### **Resistance Force:**

Improves the capacity to develop a high level of force over a prolonged period of time and increases the ability to withstand toxin accumulation and therefore delays the fatigue of the muscle. To achieve a very high intensity of work the objective is to induce a large number of long strong contractions alternating with short active rests. This is recommended for sports needing intense muscle work over a prolonged period of time.

**Intensity:** intermediate or maximum. The Active Rest should give muscle vibrations.

### **Explosive Force:**

This program increases the capacity to rapidly achieve a maximum force expression. The objective is to use the largest amount of muscles in the shortest possible time. It is advisable to start with a maximum force program in order to prepare the muscle for the Explosive Force. Recommended for all sports characterized by maximum force expressions with natural load e.g. sprint, ball sports, alpine skiing. **Intensity:** maximum, muscle contraction without pain. The Active Rest should give muscle vibrations.

### **Reactivity:**

To improve the reactivity it is necessary to increase the quick contraction of the fast fibres and the receptive activity of the myo-tactical nervous fibres. It is advisable to combine this treatment with pliometric and proprioceptive activities in order to complete the reactivity training. The muscle contractions are brief, very strong and at a very high speed. For all sports needing quick acceleration and reaction capacity e.g. ball sports, combat sports, sprint.

**Intensity:** maximum, muscle contraction without pain. The Active Rest should give muscle vibrations.

### **Muscle Relaxation:**

Enables a very effective muscle relaxation thanks to the comfortable vibrations that increase the circulation. These programs can be used at any occasion when muscle relaxation is needed.

**Intensity:** from moderate to intermediate.

**Capillarisation:**

Increases the blood supply to the muscle and improves the resistance quality and recovery. This circulation increasing, allows to maintain active the primary and secondary capillary system that improves the tissue oxygenation. This gives a delayed fatigue during very intense work. The capillarisation programs are especially recommended for sports where resistance strength and endurance are important.

**Intensity:** intermediate.

**Recovery post competition/training:**

It increases the ability to regain muscular capacity after training or competition. The stimulation varies from muscular vibrations to light contractions that allows a circulation improvement and a minor muscle hardening. This program also relieves the pain due to the endorphin production. This program is recommended for all sports because it helps to recover rapidly the maximum capacity of the muscle. It is advisable to be used after 8–24 hours of training or competition.

**Intensity:** from moderate to intermediate

## **FITNESS AND BEAUTY**

**Drainage:**

The combined stimulation achieved by the four channels of the stimulator produces a muscular pumping effect, starting at the extremity of the limb and spreading all the way to its root. The pumping action, mainly activates the venous but also the lymphatic systems by mechanical effect. Its objective is to improve the circulation and to promote the elimination of toxins.

**Intensity:** from moderate to intermediate.

**Lypolysis:**

The low frequency of this program has a positive effect on the metabolism; it increases the capacity to eliminate toxins and unwanted fluids. This method enables to model the undesirable areas modifying the metabolism of lipocytes. It improves circulation and increases the capacity to eliminate the orange skin appearance. **Intensity:** from moderate to intermediate.

**Vascularisation:**

It enables a very effective decontraction of the muscle due to the comfortable vibrations that increase the circulation and help the muscle to relax. These programs can be used at any occasion when muscle relaxation is needed.

**Intensity:** intermediate.

**Firming:**

The program is suitable for sedentary people who have weak muscles that have lost some of their bulk. It prepares the muscle to more intense stimulations. It can be used every day on large muscle groups such as the thighs, abdominal and dorsal muscles.

**Intensity:** from intermediate to elevated.

**Toning:**

The toning program offers genuine maintenance work for muscular tone. This type of training session is of moderate intensity, making it suitable for everyone without causing muscle soreness. The program is specially aimed at those who, attentive to their physical condition, are looking for a complement to their regular training. It can be undertaken three times a week.

**Intensity:** from intermediate to elevated

### **Sculpting:**

The sculpting program is the last phase in a complete aesthetic program. The program is aimed at those who desire to obtain muscular sculpting in place where traditional muscle-building exercises are insufficient. The sculpting program produces intense work and for this reason it is normal for some muscle soreness to be felt during the hours following sessions of this nature. You can prevent the soreness by executing some stretching exercises before and after the session or you can associate a relaxation program.

**Intensity:** from intermediate to elevated

### **Muscular mass increasing:**

Special program to increase the volume and the muscular mass.

**Intensity:** maximum muscle contraction without pain.

It is advisable to carry out a warm up program before starting any Firming, Toning, Sculpting, Body Sculpting and Mass program.

## **HEALTH**

Tens utilizes the nervous system's own pain relief mechanism to ease both acute and chronic pain. When stimulating, impulses are sent through the nerves to block the pain impulses. The analgesic effect occurs during the program execution but it can also last for hours. The duration of this effect varies. Tens also increases the circulation. Place the electrode on the painful area and start the program. You can carry out this program on any occasion when pain relief is needed.

Attention! TENS relieves the pain but it does not treat the cause. If the pain persists, seek medical advice.

### **TENS (Antalgic-Endorphinic)**

It offers effective relaxation of the muscle, as well as pain relief. The TENS utilizes the nervous system's own pain relief mechanism to ease both acute and chronic pain. When stimulating, impulses are sent through the nerves to block the pain impulses.

You can carry out this program on any occasion when pain relief and muscle relaxation are needed.

Attention! TENS relieves the pain but it does not treat the cause. If the pain persists, seek medical advice.

**Intensity:** ANTALGIC TENS Pins and needle sensation without pain. Should not cause muscle contraction. ENDORPHINIC TENS: visible muscle vibrations.

**TENS:** (Cervical, Lumbago, Periarthritis, Epicondylite, Tunnel Carpal, Sciatica, Knee pain, Achilles' Tendon):

The TENS utilizes the nervous system's own pain relief mechanism to ease both acute and chronic pain. When stimulating, impulses are sent through the nerves to block the pain impulses. In this program we use the stimulation alternative with light contractions that stimulate the release of endorphins.

Attention! TENS relieves the pain but it does not treat the cause. If the pain persists, seek medical advice.

**Intensity:** pins and needle sensations without pain.

## **REHAB**

### **Atrophy Recovery**

Resuming physical activity after a period of interruption. After injuries or other interruption in physical activity, the muscles often lose strength and muscle mass. It is important to restore normal muscle capacity as soon as possible. The Atrophy Recovery program, gives the possibility to train specific muscle that are affected. Can be used on any occasions when resuming activity is needed.

**Intensity:** from moderate to intermediate.

### **Muscle Relaxation**

This program enables a very effective relaxation of the muscle due to comfortable muscle vibrations that increase the circulation and help the muscle to relax. Tensed muscles are a common problem e.g. after an injury or after overworking of muscle. This kind of program can be used in any occasion when muscle relaxation is needed.

**Intensity:** moderate.

### **Haematoma:**

This program increases the absorption of the muscle haematoma caused by, for example, a contusion. It also improves the circulation in the involved arm. This program should not be used within 72 hours of the muscle injuries. Always seek medical advice for further information.

**Intensity:** from intermediate to moderate.

This is a medical device for the presence of clinic stimulation program. It is therefore certified by the Notified Body Cermat n° 0197 complied with 93/42/CEE requirements. The certification covers clinic applications

## USE OF THE MOTOR POINTS FINDER\*

User manual:

- 1) Connect the cable with the red pin to the finder and the black pin to an electrode.
- 2) Place the electrode on the interested muscle following the location shown with a (-) symbol (negative) in the corresponding picture of the electrode placement section in the electro stimulator user manual.
- 3) The Motor point is in the zone shown with a (+) symbol (positive) in the corresponding picture of the electrode placement section in the electro stimulator user manual.
- 4) We suggest applying conductive gel to ease the search of the motor point.
- 5) After connecting the cable to the electro stimulator, select the program "Search Motor points" and progressively increase the intensity (5-10mA).
- 6) Slightly press the finder on the skin and move it around the research area. The motor point should be where the electric intensity feeling is the highest.
- 7) After reducing the intensity, disconnect the finder and connect another electrode which must be placed at the exact location where the finder found the motor point.

The Motor point finder is used to locate the exact position where to place the electrodes to ensure maximum muscle contraction.

\*The Motor point finder is a GLOBUS accessory not included and available at any Globus dealer.

# COMMON QUESTIONS

## ***Why does the intensity return immediately to 0r when increasing its value?***

The unit has a control of current supply. In case the circuit is open the intensity is automatically shut off. It is therefore necessary, before starting any treatment, to be sure that the cables are well connected and the electrodes completely placed on the skin.

## ***What kind of electrodes must be used for electro stimulation?***

We suggest you to use adhesive electrodes that are more practical and improve the stimulation quality. If used with care, they will last for about 25 - 30 applications. The electrodes need to be replaced when a bad connection with your skin occurs.

## ***Where to place the electrodes?***

On the back of this manual you will find a comprehensive electrode placement guide. If you want to verify the correct placement of electrodes you can also use the special motor point pen or follow this method: place the electrode as indicated by the picture, start the stimulation and with your hand move the electrode without removing it and let the skin sliding above the muscle. You will notice a change in contraction according to the different positions. Place the electrode where the contraction is greater. Always remember to decrease the intensity to zero before moving the electrodes.

## ***What to do when the discharged battery indicator appears?***

You can continue your stimulation until it starts diminishing; then turn off the unit, disconnect the electrode cables and connect the charger to the special plug. Connect the charger to an external socket 220 V.

## ***The electro stimulation can hurt me?***

It is very difficult that electro stimulation can cause muscular problem. A very important principle to follow during stimulation is to increase the intensity gradually and to avoid pain.

## ***Is it possible to use the electric stimulator during the menstrual cycle?***

Possible interferences, like anticipations, delays, accentuations or reductions of the cycle can be subjective and variable. It is however advisable to avoid treatments in the abdominal zone during the cycle.

## ***Is it possible to use the electric stimulator during lactation?***

Until now, it has not found collateral effects in relation to the electric stimulator use during lactation. It is however advisable not to treat the thoracic zone.

## ***Does the presence of dermatological diseases like psoriasis, urticaria, etc. prevent the electric stimulator use?***

Yes, it is advisable not to treat epidermic zones, which present important dermatological diseases.

## ***When is it possible to see the first results?***

The aesthetic results obtained are always subjective. It is however possible to state that you can achieve good results about tonicity after 15 days with 3 or 4 sessions a week; for lypolysis and drainage treatment at least 40 days are necessary. Better and speeder results are obtained if treatments are combined with a good physical activity and a correct life style

## ***How many sessions I can weekly perform?***

As regards the physical training we suggest you to follow the weekly program indications reported on the Globus Personal Trainer. The weekly frequency of Fitness and aesthetic program can be varied according to the kind of treatment. Toning treatments for example should be performed 3 - 4 times a week lypolysis and drainage instead can be used every day.

# Warranty terms

The device is guaranteed to the first user for a period of twenty-four (24) months from the date of purchase against defects in materials or of the manufacturing, twelve (12) months if the user uses the device for professional purposes, provided that it is used properly and maintained under normal operating conditions.

Warranty coverage is limited in the following cases:

- Six (6) months for supplied accessories subject to wear such as batteries, chargers, power supply units, cables, G-trode handpiece.
- Ninety (90) days for the media containing software such as, for example, CD-ROMs, memory cards, etc...
- No warranty for extendable accessories and materials such as, for example, electrodes, etc..

The warranty is valid and enforceable in the country where the product was purchased. In the event that the product is purchased in any country of the European Community, the warranty is still valid in all its countries.

In order to take advantage of the warranty service, the user must comply with the following warranty clauses:

1. The products, and all accessories, to be repaired must be sent by and at the expenses of the customer in their original packages.
2. The product's warranty is subject to the production of a fiscal document (fiscal receipt, receipted bill or sale invoice), attesting the product's purchase date.
3. The repair work shall have no effect on the original warranty expiry date and shall neither renew nor extend it.
4. If, during the repair work, no defects are found, the costs related to inspection times shall in any case be charged.
5. The warranty becomes void where the fault has been caused by: impacts, falls, erroneous or improper use of the product, use of non-original power supply unit or charger, accidental events, alteration, replacement/detachment of the warranty seals and/or tampering with the product. The warranty does not cover damages caused during transportation when unsuitable packages are used.
6. The warranty does not cover the inability to use the product, other incidental or consequent costs or other expenses incurred by the purchaser.

NOTE: Before returning the device for repairs, we recommend to read carefully the user instructions contained in the manual and visit Globus website.

In case you need to return your product for assistance, contact your dealer or contact Globus Service. The manufacturer reserves the right to make changes without previous notice. The features and dimensions reported in this manual are not binding.

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**Prodotto da: SAVIA L.t.d. - 4nd Industrial District, Feng Huang Village (Shenzen) - Made in PRC - sales@tyece.com**

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**This product as been manufactured by: SAVIA L.t.d. - 4nd Industrial District, Feng Huang Village (Shenzen) - Made in PRC - sales@tyece.com**

This product has been manufactured according to the technical regulations in force and it has been subjected to certifications according to 93/42/CEE directive for medical devices in order to assure the product's security. The user must follow all prescriptions and restrictions shown in the operating manual in order to keep the maximum security level of the product. The manufacturer declines all responsibility for every different use from what indicated in this operating manual. All complete or partial reproductions of texts and/or pictures of this operating manuals made by any electronic or mechanic equipments are strictly forbidden without written permission of the manufacturer.

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