Tesmed+

Max 830

ELETTROSTIMOLATORE PROFESSIONALE
PROFESSIONAL ELECTRICAL MUSCLE STIMULATION
STIMULATEUR ÉLECTRIQUE PROFESSIONNEL
PROFESSIONELLE ELEKTROSTIMULATIONSGERÄT
ELECTROESTIMULADOR PROFESIONAL

BREVETTO INTERNAZIONALE
INTERNATIONAL PATENT
BREVET INTERNATIONAL
INTERNATIONALES PATENT
PATENTE INTERNACIONAL

MANUALE D'USO - INSTRUCTION MANUAL - MODE D'EMPLOI
GEBRAUCHSANWEISUNGEN - MANUAL DE INSTRUCCIONES

Indice

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| Prog. | Onda | Parametri | Indicazione Trattamenti |
|-------|---------------|---|--|
| 1 | Waims-System® | 1 → 2 → 3 → 4 → (pause 1) 1 → 2 → 3 → 4 → | Azione Linfodrenante**, Rilassante, Recupero attivo, Defaticante |
| 2 | Waims-System® | 1 → 2 → 3 → 4 → (pause 2) 1 → 2 → 3 → 4 → | |
| 3 | Waims-System® | 1 → 2 → 3 → 4 → (pause 3) 1 → 2 → 3 → 4 → | |
| 4 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 → (pause 1) 1 → 2 → 3 → 4 → 3 → 2 → 1 → | |
| 5 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 → (pause 2) 1 → 2 → 3 → 4 → 3 → 2 → 1 → | |
| 6 | Tapping 1 | Tpulse=100µs Freq.=1÷100Hz Time=1'÷60' | Massaggio muscolare, Riattivante, Cellulite |
| 7 | Tapping 2 | Tpulse=200µs Freq.=1÷100Hz Time=1'÷60' | |
| 8 | Tapping 3 | Tpulse=300µs Freq.=1÷100Hz Time=1'÷60' | |
| 9 | Tapping 4 | Tpulse=400µs Freq.=1÷100Hz Time=1'÷60' | |
| 10 | Press 1 | Tpulse=120µs Freq.=1÷150Hz Time=1'÷60' | Forza resistente, Forza esplosiva, Sviluppo forza |
| 11 | Press 2 | Tpulse=240µs Freq.=1÷150Hz Time=1'÷60' | |
| 12 | Press 3 | Tpulse=300µs Freq.=1÷150Hz Time=1'÷60' | |
| 13 | Kneading 1 | Tpulse=40µs Freq.=1÷150Hz Time=1'÷60' | Tonificazione, Modellamento, Cellulite, Rassodamento |
| 14 | Kneading 2 | Tpulse=160µs Freq.=1÷150Hz Time=1'÷60' | |
| 15 | Kneading 3 | Tpulse=200µs Freq.=1÷150Hz Time=1'÷60' | |
| 16 | Kneading 4 | Tpulse=340µs Freq.=1÷150Hz Time=1'÷60' | |
| 17 | Modulation 1 | Tpulse=80µs Freq.=1÷150Hz Time=1'÷60' | Drenante, Stretching |
| 18 | Modulation 2 | Tpulse=140µs Freq.=1÷150Hz Time=1'÷60' | |
| 19 | Burst 1 | Tpulse=180µs Freq.=1÷150Hz Time=1'÷60' | Aerobico, Ipertonia |
| 20 | Burst 2 | Tpulse=260µs Freq.=1÷150Hz Time=1'÷60' | |

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contents

| | | | | |
|---|---|---|---|----|
| > | > | > | symbols used | 26 |
| > | > | > | contraindications | 27 |
| > | > | > | warnings | 28 |
| | | | Precautions | 28 |
| | | | Possible side effects | 28 |
| > | > | > | electro stimulation | 29 |
| | | | Electrodes | 29 |
| | | | Tiredness | 29 |
| | | | Advantages of using electro stimulation | 29 |
| | | | Results | 30 |
| > | > | > | high technology | 30 |
| > | > | > | battery charger | 30 |
| | | | Replacing the battery | 31 |
| | | | Connections | 31 |
| | | | Disconnecting the electrodes | 31 |
| > | > | > | equipment description | 32 |
| | | | LCD display and keys | 32 |
| > | > | > | use | 32 |
| | | | Switching on | 32 |
| | | | Modes | 32 |
| | | | Programs 01 to 20 | 32 |
| | | | Programs 21 to 99 | 33 |
| | | | Intensity control | 33 |
| | | | Frequency control (for programs 01 to 20) | 33 |
| | | | Error warning | 33 |
| | | | "Keyboard block" function | 34 |
| | | | "ON block" function | 34 |
| > | > | > | maintenance | 34 |
| > | > | > | disposal | 35 |

contents

| | | | | |
|---|---|---|--|---------|
| > | > | > | normative references | 35 |
| | | | TESMED MAX 830 technical data | 36 |
| | | | Battery charger technical data | 36 |
| | | | Guide and manufacturer's statement: electromagnetic emission and immunity | 37 |
| > | > | > | how to choose the right treatment to suit individual needs | 39 |
| > | > | > | electrode polarization | 39 |
| > | > | > | maintenance | 39 |
| > | > | > | personalised programs | 40 |
| | | | (programs 01 to 20) | |
| > | > | > | Tesmed Max 830 programs | 41 |
| | | | (programs 21 to 99) | |
| | | | FITNESS treatments | 41 |
| | | | BEAUTY treatments | 42 |
| | | | SPORT treatments | 45 |
| | | | MASSAGE TENS treatments | 48 |
| > | > | > | positioning of electrodes | 127-134 |

symbols used

Important note, if not respected, there may be a risk for the user and/or the equipment may be damaged. We recommend that you read the information in this manual carefully and keep it in a safe place.

Low frequency parts applied.

Internally fed equipment (DC).

Parts that must be disposed of as special waste.

Read the operator's manual

Manufacturer

Double electrical isolation

IP

Protection against dust and liquid penetration

contraindications

ATTENTION:

Electro stimulation is not recommended in areas with any of the following:

- inflammation, open wounds, strained or torn muscles, dermatitis, bacterial infections, eritema or skin irritation;
- muscular hernias;
- fractured bones;
- prosthesis or osteosynthesis;
- varicose veins or vascular disease in general;
- neoplasia;
- mucous membranes.
- in the thoracic area to people with a cardiac arrhythmia

This medical device is prohibited in the following cases:

- pregnant women;
- growing children (risk of juvenile heart attack);
- patients wearing pace-makers (risk of interference with the functioning of the pace-maker);
- implanted defibrillators;
- alteration of the peripheral nerves, muscular disease and serious damage to the joints;
- patients with metal joint prosthesis and/or osteosynthesis and/or intratissue/endotissue metal;
- serious heart problems;
- patients suffering from epilepsy;
- patients with bone tumours;
- patients with serious high blood pressure;
- patients suffering from serious systematic, debilitating or dermatological disease;
- patients undergoing treatment with an HF device, as this could cause burns and damage to the stimulator itself;
- near any short wave or microwave therapy equipment, as this could cause instability in the stimulator output;
- unconscious patients or patients unable to communicate;
- people suffering from psychophysical problems that prohibit their exposure to external factors;
- in the presence of a mix of flammable anaesthetic with air, oxygen or nitrogen protoxide.

Never place the electrodes on the carotid artery and/or on a spastic muscle.

Use of the TESMED MAX 830 for therapy purposes should be supervised and recommended by a doctor. Interrupt application immediately in the event of discomfort and obtain medical advice.

warnings

- You must take into account the technical characteristics of each program and read the relevant paragraphs in the User Manual.
- The Tesmed Max 830 must not be used for purposes other than those indicated.
- You must begin each treatment with the current emission value set at zero (min channel power setting) and then gradually increase and/or reduce the power to suit your tollerability.
- Max power does not necessarily mean better performance or results.
- We recommend using original electrodes only to guarantee top performance, as these have been developed specifically for the Tesmed Max 830. Non original electrodes could burn the skin.
- The Tesmed Max 830 must be kept out of the reach of children.
- Electrocardiogram monitoring equipment may be disturbed by the current used for the electro stimulator.
- We do not recommend using the Tesmed Max 830 on children under the age of 14 or the elderly over 75.
- The electrodes must not be placed less than 10 cm from the eyeballs.
- Do not use the equipment if there are special health problems without first obtaining medical advice from a doctor.
- Should the equipment malfunction, do not try to open or repair the equipment, but send the product in its original packing to **Feldi S.r.l.**, describing the problems found on the inside of the package (see warranty conditions).
- The bipolar cables not connected to the body for the treatment must be disconnected from the device
- Device not suitable for simultaneous use with high frequency surgical devices
- The application of electrodes close to the chest may increase the risk of heart fever
- The electrodes should never be placed on the carotid breast, along the head, directly on the eyes, covering their mouths and / or on a spastic muscle.
- The Tesmed Max 830 electrostimulator must be powered only to a very low voltage supplied by the supplied power supply. If the power supply and the cable that connects the unit is damaged, it must be replaced by a power supply and cable of the same model supplied with the product.

Precautions

- Do not tamper with the equipment or try to open it.
- Keep the electro stimulator away from sources of heat, powder, sprayed water, humid atmospheres and avoid direct exposure to the sun.
- Avoid violent knocks that could damage the equipment.

Possible side effects

- Isolated cases of skin irritation may occur in people with highly sensitive skin. Should there be an allergic reaction to the electrode gel, stop the treatment and consult a specialist doctor. It is normal for the skin to be slightly red after the treatment in the areas where the electrodes were placed: this will disappear after roughly 30 minutes. We strongly advise against treating any areas of the body with serious skin problems.

electro stimulation

Electro stimulation is a technique that provokes muscular contractions very similar to voluntary contraction thanks to electrical impulses acting on the motor points of the muscles. As everyone knows, muscles contract when they receive an impulse coming from the nervous system and then return to their normal length once they are no longer stimulated thanks to their elasticity. A motor nerve fibre splits into several ramifications before entering the muscle and ends at the motor end plate. The nerve fibre and the group of muscular fibres innervated by this make up the motor unit. A muscular fibre responds to the electrical stimulus only when the intensity of the stimulus reaches a threshold value, which is specific to each fibre. The most obvious use of the electro stimulation of innervated muscles is in kinesiotherapy of atrophy, functional reanimation after a long periods of immobility, in vasodilation, pain reducing and in stimulating the treated tissues to combat cellulite (PASSIVE EXERCISE). When treating a patient with electro stimulation, the following are crucial: the correct positioning of the electrodes, the stimulation program and the wave type.

Electrodes

The positioning of the electrodes is crucial. It must provide accurate stimulus parameters for the zone being treated since the conductivity of body tissues varies according to their structure, aqueous and ionic make-up. The electrodes must be in perfect contact with a thoroughly clean skin (also depilated, if necessary). This ensures better penetrability and permeability of the skin when receiving the stimulus to suit the application to be carried out. **For better electrode conductivity we recommend they be replaced when they loose adherence. They can also be cleaned with running water or wiped with a damp cloth; leave to dry before using them again for a few more applications.** NB: With electro stimulation, the current density diminishes as the tissue consistency increases.

We recommend using Tesmed electrodes exclusively, as electrodes of different sizes and composition could cause burns.

Tiredness

The high percentage of simultaneous contraction/relaxing of the muscles due to electrical stimulus often results in a drop in the relaxation of the neurotransmitters, thus resulting in muscle fatigue. The scientific explanation is as follows: the stimulus tends to provoke the recruitment of large diameter motor neurones (these have a low threshold), thus exciting the more consistent and faster muscle fibres (FIBRES 2 or white fibres). These tire more easily, but they use less effort than FIBRES 1 or red fibres. **In other words, the higher the frequency of stimulus, the more tired the muscle becomes.** For this reason, it is important to include relaxation programs and get medical advice. Having said this however, Tesmed's sports programs already contain a relaxation phase. **During treatments using the electro stimulator, every increase in intensity must be done gradually. If there is EXCESSIVE contraction, you must reduce the intensity.** In fact, stimulation must not provoke muscular tetany (i.e. long continuous contractions). **There are no standard rules for the number of contractions and the duration and frequency of the sessions. These must be set to suit each person and the condition of the tissues to be treated.** When treating areas affected by CELLULITE, you should normally start the treatment by gradually raising the intensity until the patient starts to feel a slight tingling, but without creating excessive contraction of the muscles or an unpleasant sensation of localised heat.

To sum up: the vasomotor and metabolic effects induced by the passage of an electrical current considerably increase the consistency and quality of tissues. The best level of stimulus is when the patient perceives the current passing through the tissue; this should provoke a light initial contraction of the muscles (a slight trembling), but never too intense.

There are two important advantages here:

- 1) you avoid the risk of muscular over-contraction and the accumulation of lactic acid;
- 2) the frequency of muscular contractions may be much higher than with traditional gym, to get best results faster.

Advantages of using the electro stimulator

In beauty: improvement of capillary micro circulation with more oxygen supplied to tissues, thus producing better fat mobilisation; toning up and firming up of tissues; draining of excess liquids; an increase in metabolism and cellular activities; a reduction of localised fatty tissues and cellulite.

In sports: isolation of the muscular group being treated; saving of physical and psychological energy; recruitment of a larger number of muscular fibres than is the case with voluntary contraction; training of under-used muscular fibres during normal work-out sessions; shorter recovery periods; protection of the body's system of tendons and muscles; improved blood circulation. Muscular damage due to electro stimulation is most unlikely, as long as the intensity is increased gradually. Any interference in the menstrual cycle (e.g. early, late, heavy or light periods) are extremely subjective and variable. We recommend, therefore, to avoid treatments in the abdominal area during the actual period and the days immediately before and after this. To date, no side effects have been noted concerning the use of electro stimulation during lactation. However, we recommend not treating the upper body area during this time.

We strongly advise against treating any parts of the body with serious skin problems.

Results

All results in the aesthetics field are naturally subjective. We can, in any case, sustain that in the case of TONE UP, good results can be seen after just 15 days with 5-7 regular and constant weekly sessions; cellulite treatment will require 30-40 days. Better and quicker results may be obtained if treatments are added to suitable physical activity and a correct lifestyle.

high technology

Tesmed Max 830 is designed and built with cutting-edge electronic components, a high storage microcontroller with great memory, processing and storage capacity. This means that the Tesmed Max 830 is an extremely professional product with superior performance, yet compact size. The software and the hardware in the Tesmed Max 830 guarantee a special self-calibration function for the output parameters to suit the different characteristics of every human body.

The Tesmed Max 830 programs include a sequential pulsed treatment. Thanks to the new patented Waims-system®, the Tesmed Max 830 electro stimulator offers a valid alternative to manual massage (lymphatic drainage). This stimulation is an exclusive feature of the Tesmed Max 830 that, thanks to this draining, helps blood circulation remarkably, since the extracellular liquids and built-up fats are moved and can thus be removed by the body. This type of massage leaves your legs and arms feeling lighter and gives the skin an even look, at the same time considerably reducing the "orange peel" effect. All areas treated with the patented WAIMSSYSTEM® gain better tone and are firmer.

battery charger

The standby of the battery depends on the treatments used, the number of channels used and the working intensity.

The symbol  appears on the Tesmed Max 830 display when the battery is running low; if this symbol  starts flashing during a treatment, we recommend that you recharge the battery as soon as the treatment has ended. If the symbol  appears just a few minutes after switching on the

Tesmed Max 830 (e.g. during parameter setting), you must recharge the battery immediately before starting any treatment.

To recharge the Tesmed Max 830 battery, disconnect all the cables, plug in the battery charger and connect the low voltage connector to the electro stimulator (fig. 2 pag. 123).

We recommend recharging the battery immediately after purchase.

Complete recharging of the battery takes approximately 8-10 hours.

During the first phase of recharging (the first 3-4 hours), the symbol  flashes on the display. Then it is steady and so the equipment is ready for use. For max autonomy, we recommend continuing with the recharging for a total of 8-10 hours. The battery charger can always stay connected to the Tesmed Max 830 for a longer period without damaging the equipment. After complete recharging, you may find that the display is blank.

- IF THE DEVICE DURING THE FIRST LOAD DOES NOT READ THE "FULL CHARGE" OF THE BATTERY AND THERE IS STILL THE SYMBOL OF THE FLASHING BATTERY , THE DEVICE THEREFORE IS FULLY USEFUL AND THE READING OF THE FULL LOAD WILL REACH MAINLY WITH THE SUCCESSIVE LOADS.

Never recharge the equipment when the cables and electrodes are connected to a person and never recharge the battery by using a different battery charger: always use the one supplied with the equipment. Do not short-circuit the battery.

Replacing the battery

The Tesmed Max 830 works thanks to a rechargeable battery at the rear of the equipment. Its performance may drop after multiple recharging.

Should the battery's autonomy be far lower than normal, you must replace it with another battery. NB: use only Tesmed approved battery types.

- To access the battery compartment, simply remove the removeable belt hook and slide the cover as shown by the arrow (fig. 3 pag. 123).
- The battery is connected to the equipment by means of a connector; to replace the battery, extract it and then slide the connector out from the equipment (fig. 4 pag. 123).
- Take special care over the connection when inserting the new battery: the connector has a flat surface on one side that must be inserted facing upwards (fig. 5 pag. 123).

NB: Used batteries must be disposed of properly to protect the environment.

Connections

The electrical pulses generated by the Tesmed Max 830 are transmitted to the muscles by means of self-adhesive electrodes. Before switching the Tesmed Max 830 on and before using the intensity command, always apply the electrodes to the body part to be treated and then connect them via the cables to the Tesmed Max 830: insert the pins on the bipolar cables in every pair of electrodes used for the stimulation and then connect them to the equipment by inserting them in the relevant colour coded seat on the rear of the Tesmed Max 830 (fig. 6 pag. 123). If this procedure is not followed, there is the risk that the equipment will exceed a set threshold and so produce an error message. Thanks to an in-depth study and the TMS "Tesmed Multiplexer used System" technology used in the Tesmed Max 830, the equipment is designed to work correctly without having to apply a specific polarity to the electrodes, thus making it easier to apply, connect and use.

Disconnecting the electrodes

After treatment, disconnect the pins on the bipolar cables as shown figure 7 opag. 123. Take care not to pull the connection cables from the electrodes.

equipment description (fig. 8 pag. 124)**LCD display and keys**

| | | | |
|---|--------------------------------------|---|--|
| a | Tesmed logo | k | ON/OFF and select channels double function key |
| b | Program number | l | Intensity symbol |
| c | Program symbol | m | Keyboard block and ON block symbol |
| d | Low battery symbol | n | Timer symbol |
| e | Frequency symbol | o | Treatment time and remaining time |
| f | Selected output channels (4) | p | Clock symbol |
| g | Relative intensity for each 4 output | q | Output channels |
| h | Function keys (MODE) | r | Battery recharger socket |
| i | Decrease key (⏮) | s | Battery door |
| j | Increase key (⏭) | t | Removable belt hook |

use

The Tesmed Max 830 has **99** programs (01 to 99).

Programs **01** to **20** can be personalised and you can change the frequency and treatment time values. Programs **21** to **99** are pre-set with a fixed duration and cannot, therefore, be changed.

Before using the Tesmed Max 830, you must choose the treatment type you want, using the personalised or set treatment tables, find the relevant program and then apply the electrodes to the body part you want to treat and connect them to the equipment.

Switching on

To switch the Tesmed Max 830 on, press key  for approximately 2 seconds: the system runs a self-test procedure, during which the display shows all the available functions (fig. 9 pag 125); then appears on the display what shown in fig. 10 pag 125.

- "Display lighting" function:

The display automatically switches off after about 20 seconds when the equipment is not in use. The display is lit again when you press any key.

- "Self switching-off" function

The Tesmed Max 830 automatically switches off is not used for more 5 minutes in order to save the battery.

Modes

Press key  repeatedly to select the number of channels you want to use to suit the treatment and then key  to select the program (fig. 11 pag 125).

Programs 01 to 20

Programs 01 to 20 let you personalise the treatment time, once selected. Press key  and the relevant symbol **TIME** starts flashing on the display (fig. 12 pag 125).

Use keys   to adjust the treatment time: min 1 minute to max 60 minutes. Press key  again to control the intensity.

Programs 21 to 99

After having selected a program from 21 to 99, press key to go straight to the intensity control.

Intensity control

By "intensity" we mean the value of current emitted to the body part and thus the degree of stimulation causing the muscle to contract. Symbol appears on the display (fig. 13 pag 125): you can now adjust the intensity using keys and (1÷99). The intensity rises or falls on all the selected channels at the same time. When you set the intensity, the clock begins to flash indicating that the equipment is providing energy (fig. 14 pag 125). You can, however, access each channel during use to control the intensity independently: simply press key to select the channel whose intensity you want to change and then use keys and to increase or decrease the value. Fig. 15 pag 125 shows that the intensity values for the 4 channels can be different.

The display constantly updates the treatment time until completion, as well as showing all the set parameters, meaning that you can always control the treatment.

At the end of the treatment, the Tesmed Max 830 clears the intensity setting for all output channels. You can now switch the equipment off by pressing key for 2 seconds or set new parameters to begin a new treatment. Use keys and , to increase or decrease the value of any parameter (in steps of one unit). Keep pressed for faster increases/decreases.

NB: *To perceive a change in intensity in a single channel or all channels, changes are only made when the muscle is contracting. This avoids the risk of uncontrolled or sudden changes in intensity. During the treatment you can quickly clear the intensity setting by pressing key briefly. All the default parameter settings are restored and so you can start increasing/decreasing the intensity again to resume the treatment.*

Frequency control (for programs 01 to 20)

Programs 01 to 20 let you correct the pulse frequency value at will. The "frequency" is the number of pulses per second. When you change the frequency, you feel an increase or decrease in the speed of stimulation. These programs are identified by a special symbol on the display (fig. 16 pag 125). After starting the treatment, press key for 3 seconds to enter the frequency control mode. Press keys and repeatedly to increase or to decrease the frequency. The symbol flashes (at a varying speed) to indicate frequency control mode. Press key once to quit frequency control mode.

NB: *if in frequency control mode during a treatment and you want to clear the intensity, briefly press key twice. All the default parameter settings are restored and so you can start increasing/decreasing the intensity again to resume the treatment.*

Error warning

The Tesmed Max 830 automatically goes into error mode if an electrode comes loose. In this case, the equipment display will show "er" (flashing for about 5 seconds) against the channel with the current distribution problem (fig. 17 pag 125). All pulses are stopped immediately and "00" appears as the intensity value on the display. The error message can mean that:

- There are no electrodes connected to the channel in error mode.
- The electrodes are not properly adhered to the skin: reposition the electrodes or replace if too old.
- The cable used to connect the electrodes is faulty: check by inverting two cables.

The system automatically resets itself once the intensity value increases. Make sure the electrode is connected properly and then use key  to select the channel in error mode and adjust the intensity. If the problem persists, the equipment repeats the error warning: check the state of the cables and that the electrodes adhere properly to the body. The error warning may be delayed when using low intensity values, or even not detected, as it is influenced by the difference in each person's resistance.

"Keyboard block" function

Press key  and, at the same time, press key  and release. Symbol  appears on the display (fig. 18 pag 125). Follow the same procedure to unblock the keyboard.

This function, when enabled, disables all the keys and so avoids the risk of accidental changes to the set parameters during treatment.

"ON block" function

Press key , and, at the same time, press key , for 3 seconds. Symbol  and the Tesmed logo appear on the display (fig. 19 pag 125). Follow the same procedure to switch the equipment on.

We recommend using this function when storing the equipment in its case or bag, thus avoiding the risk of accidentally switching it on.

maintenance

The Tesmed Max 830 is a highly sophisticated product and must be handled with care. The maintenance tips below will help you ensure its long life.

- Do not allow it to get wet. Rain, humidity and liquids or condensation generally contain corrosive mineral substances that can damage the electronic circuits.
- Keep the Tesmed Max 830 out of the reach of children.
- Do not leave the Tesmed Max 830 in particularly hot rooms. Excessive temperatures can affect the life of the electronic circuits, damage the batteries and deform the plastic parts.
- Do not leave the Tesmed Max 830 in particularly cold rooms, as condensation could form when it heats up during use and so damage the electrical circuits.
- Do not open the Tesmed Max 830. Any tampering with the equipment by inexperienced people could damage it.
- Do not allow it to fall, be knocked or shaken, as the electronic circuits inside could be irreparably damaged.
- Do not use the battery for purposes other than those intended.
- Never use the battery charger or batteries if damaged.
- Do not short circuit the battery.
- Use only the battery charger provided. Contact Customer Service in the event of a malfunction.
- Never leave the battery in hot or cold places, such as a closed car in the height of Summer or during the Winter, as this will reduce its capacity and life.
- Never throw used batteries in the fire!
- Follow local regulations to dispose of used batteries, e.g. recycle them. Never throw them away with household waste.

- To clean the Tesmed Max 830 (with the equipment switched off), use a soft cloth and alcohol, but no solvents, as these could ruin the plastic parts.

Use only original spare-parts and products for maintenance.

Follow these maintenance times to keep the Tesmed Max 830 in top working condition:

- A)** Check the state of the electrode cables and check for signs of oxidation (monthly if used regularly; fortnightly if used intensely).
- B)** Replace the electricity conducting electrode as soon as these show signs of wear. The manufacturer's warranty is invalidated if the above maintenance is not done. The manufacturer guarantees a 5 year working life of the equipment provided it is used and maintained correctly. The equipment, batteries and its accessories must be disposed of properly (treat as special waste) to protect the environment.

disposal

The equipment, batteries and its accessories must be disposed of properly (treat as special waste) to protect the environment.

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TESMED MAX 830 TECHNICAL DATA

| | |
|--|--|
| Power supply | Nickel-metal-hydrate (NI-MH) rechargeable battery 4.8V = 750 mA/h |
| Casing | Polycarbonate |
| Pulse repetition frequency | 1Hz - 150Hz |
| Max energy per pulse | E= 3200 μ J |
| Max peak power / mean power | 10 W per pulse / 5 mW |
| Mean current at max frequency Mean current at min frequency | 2 mA 0,2 mA |
| Max peak current | 120mA |
| Pulse shape | Rectangular, biphasic (positive and negative), no continuous component |
| Pulse duration (Positive/Negative) | da 40 μ s a 400 μ s |
| Pre-set programs | 99 |
| Number of channels | 4 |
| Dimensions (mm) | 128 x 58 x 18 |
| Weight (gr) | 145 |
| Applied part: Low Voltage | BF |
| Conditions of use | Continuous use |
| All electrical data are provided for a charge of 1000 ohm per channel | |

Equipment not suitable for use in the presence of inflammable anaesthetic mixtures with air or oxygen or nitrous oxide.

BATTERY CHARGER TECHNICAL DATA

| | |
|---|--|
| Electrical input | 100-240V \sim , 50/60Hz, 180mA |
| Electrical output | 6V \sim , 60mA |
| Polarity, external \emptyset , internal \emptyset |  3mm, 1.1mm |
| Compliance with technical standards | EN 60601-1, 3a ed. |

Guide and manufacturer's statement - electromagnetic immunity

The Tesmed Max 830 electro stimulator is designed to work in the electromagnetic environment described below. The customer or user of the Tesmed Max 830 electro stimulator must guarantee that it used in such an environment.

| Immunity test | IEC 60601 test level | Conformity level | Electromagnetic environment - guide |
|---|---|--------------------------------|--|
| Electrostatic discharge (ESD) IEC 6100-4-2 | ± 8kV contact ± 15kV in air | ± 8kV contact ± 15kV in air | Floors must be wooden, concrete or tiled. If floors are covered with a synthetic material, the relative humidity should be at least 30%. |
| Transiting / fast electric trains IEC 6100-4-4 | ± 2kV for power lines ±1kV for input/output lines | Not applicable | The Tesmed Max 830 electro stimulator is a device with DC power supply from an incorporated battery |
| Overvoltage IEC 61000-4-5 | ±1kV differential mode ±2kV common mode | Not applicable | The Tesmed Max 830 electro stimulator is a device with DC power supply from an incorporated battery |
| Voltage drops, brief cuts and variations in voltage on supply input lines IEC 61000-4-11 | 5% U_T (>95% drop in U_T) for 0,5 cycles 40% U_T (60% drop in U_T) for 5 cycles 70% U_T (30% drop in U_T) for 25 cycles <5% U_T (>95% drop in U_T) for 5 seconds | Not applicable | The Tesmed Max 830 electro stimulator is a device with DC power supply from an incorporated battery |
| Magnetic field at mains frequency (50/60 Hz) IEC 61000-4-8 | 30 A/m | Not applicable | The Tesmed Max 830 electro stimulator is a device with DC power supply from an incorporated battery |

NOTE: U_T is the AC mains voltage rate prior to application of the test level.

Recommended separation distances between portable and mobile radio frequency equipment and the Tesmed Max 830 electro stimulator

The Tesmed Max 830 electro stimulator is designed to work in an electromagnetic environment where irradiated RF interference levels are kept under control. The customer or the user of the Tesmed Max 830 electro stimulator may help avoiding electromagnetic interference by guaranteeing a min distance between mobile and portable RF communication devices (transmitters) and the Tesmed Max 830 electro stimulator, as recommended below and to suit the max output power of the radio communication devices.

| Max output power of the specified transmitter (W) | Separation distance at the transmitter frequency (m) | | |
|---|--|-----------------|-----------------|
| | 150kHz - 80MHz | 80MHz - 800MHz | 800MHz - 2,5GHz |
| | $1,17 \sqrt{P}$ | $1,17 \sqrt{P}$ | $1,17 \sqrt{P}$ |
| 0,01 | 0,12 | 0,12 | 0,23 |
| 0,1 | 0,37 | 0,37 | 0,74 |
| 1 | 1,17 | 1,17 | 2,33 |
| 10 | 3,7 | 3,7 | 7,37 |
| 100 | 11,7 | 11,7 | 23,3 |

Transmitters with a max output power not indicated above: the recommended separation distance in metres (m) can be calculated using the applicable equation for the transmitter frequency, where P is the max nominal output power of the transmitter in Watts (W) as indicated by the manufacturer of the transmitter and d is the recommended separation distance in metres (m).

NOTE 1: At 80 MHz and 800 MHz, the highest frequency interval is applied.

NOTE 2: These guidelines may not apply to all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.

Guide and manufacturer's statement - electromagnetic immunity

The Tesmed Max 830 electro stimulator is designed to work in the electromagnetic environment described below. The customer or user of the Tesmed Max 830 electro stimulator must guarantee that it used in such an environment.

| Immunity test | IEC 60601 test level | Conformity level | Electromagnetic environment - guide |
|--------------------------------|----------------------------|------------------------|--|
| RF carried IEC 61000-4-6 | 3Vrms da 150kHz a 80MHz | Not applicable | Portable and mobile RF communication devices should not be used any closer to any part of the Tesmed max 830 electro stimulator, including the cables, than the recommended separation distance calculated using the applicable equation for the transmitter frequency. Recommended separation distance: $d = 1,17\sqrt{P}$ $d = 1,17\sqrt{P}$ 80MHz - 800MHz $d = 2,33\sqrt{P}$ 800MHz - 2,5GHz |
| RF irradiated IEC 61000-4-3 | 3V/m da 80MHz a 2,7GHz | 3V/m 80MHz - 2,7GHz | where P is the max nominal output power of the transmitter in Watts (W) as indicated by the manufacturer of the transmitter and d is the recommended separation distance in metres (m). The intensity of fixed RF transmitters (as established by an electromagnetic survey of the site) could be less than the conformity level in each frequency interval. There may be interference near devices marked:  |

NOTA 1: At 80 MHz and 800 MHz, the highest frequency interval is applied.

NOTA 2: These guidelines may not apply to all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.

A) a. Field intensity for fixed transmitters, such as radio phone stations (mobile phones and cordless phones) and walkie-talkies, radio ham equipment, AM and FM radio equipment and TV transmitters cannot be foreseen in theory or with precision. To assess an electromagnetic field causes by fixed RF transmitters, you need to carry out an electromagnetic survey of the site. If the field intensity found in the place where the Tesmed Max 830 electro stimulator is to be used exceed the applicable conformity level above, you will need to observe the normal working of the Tesmed Max 830 electro stimulator. If there are anomalies, extra measures will be needed, such as a different position or direction of the Tesmed Max 830 electro stimulator.

B) b. The field intensity in the frequency range 150 kHz to 80 MHz should be less than 3V/m.

Guide and manufacturer's statement - electromagnetic emissions

The Tesmed Max 830 electro stimulator is designed to work in the electromagnetic environment described below. The customer or user of the Tesmed Max 830 electro stimulator must guarantee that it used in such an environment.

| Emission test | Conformity level | Electromagnetic environment - guide |
|---|------------------|--|
| RF emissions CISPR 11 | Group 1 | The Tesmed Max 830 electro stimulator only uses RF energy for its internal functioning. Therefore, its RF emission is very low and is unlikely to cause any interference with any electronic equipment nearby. |
| RF emissions CISPR 11 | Class B | — |
| Harmonic emissions IEC 61000-3-2 | Not applicable | The Tesmed Max 830 electro stimulator is a device with DC power supply from an incorporated battery |
| Voltage fluctuation emissions / flicker IEC 61000-3-3 | Not applicable | The Tesmed Max 830 electro stimulator is a device with DC power supply from an incorporated battery |

how to choose the right treatment to suit individual needs

- You can choose the most suitable treatment from those in the table of treatments below, consulting the list of photos for details of how to position the electrodes on the muscles in question.
- The photo number is shown next to each treatment, in the "Photo" column. If more than one photo is indicated for the same treatment (e.g. **1, 37, 40-44**), select photo **1, 37, 40, 41, 42, 43** or **44** as required.
- We recommend that you do not treat a group of muscles with more than one session a day.
- We recommend that you do not exceed 2 hours of electro stimulation per day, even if involving different groups of muscles.
- **The treatment times shown in the tables is not binding: actual treatment times may vary considerably from one person to another.**
- If in doubt about which treatment to use, we advise you consult a doctor or a trainer.
- During the treatment, the output intensity is controlled automatically at each change in phase for greater safety. You can then increase/decrease the intensity at will.
- To ensure your safety, you cannot increase the intensity during the pause between one pulse and the next.

***NB:** The programs in the **Tesmed Max 830** are modulates to provide various stages* of massage (1 to max 5). Some programs have one or more Waims-system® stages. During treatment (and especially with the Waims-system®) we always recommend that you adapt the intensity of each channel, as your sensitivity will vary various from one group of muscles to another.*

* The number of stages in a program corresponds to the number of wave types used in it.

electrode polarization

Most commonly available electro stimulators adopt a common electrode standard: they either contain continuous current (DC) or generate modified waves (different positive/negative values, with more of one or the other). Hence the need to specify the polarity. After careful study and with the use of TMS "Tesmed Multiplexer used System" technology in the Tesmed Max 830 means that this equipment can work correctly without having to give a specific polarity to the the electrodes, thus making them easier for you to apply, connect and use.

maintenance

Follow the same program as for the treatment, but only 2-3 times a week and not on two consecutive days, unless otherwise indicated by your doctor or trainer

PERSONALISED PROGRAMS

(programs 01 - 20)

positioning of
electrodes

- You can change the time and frequency of the programs in the table below at your discretion.

| Prog. | Wave | Parameter | Treatment Indication |
|-------|---------------|---|---|
| 1 | Waims-System® | 1 → 2 → 3 → 4 (pause 1) 1 → 2 → 3 → 4 | Lymph draining, Relaxing, Active recovery, Tiredness |
| 2 | Waims-System® | 1 → 2 → 3 → 4 (pause 2) 1 → 2 → 3 → 4 | |
| 3 | Waims-System® | 1 → 2 → 3 → 4 (pause 3) 1 → 2 → 3 → 4 | |
| 4 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 1) 1 → 2 → 3 → 4 → 3 → 2 → 1 | |
| 5 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 2) 1 → 2 → 3 → 4 → 3 → 2 → 1 | |
| 6 | Tapping 1 | Tpulse=100µs Freq.=1÷100Hz Time=1'÷60' | Decontracting massage, Reactivating, Cellulite |
| 7 | Tapping 2 | Tpulse=200µs Freq.=1÷100Hz Time=1'÷60' | |
| 8 | Tapping 3 | Tpulse=300µs Freq.=1÷100Hz Time=1'÷60' | |
| 9 | Tapping 4 | Tpulse=400µs Freq.=1÷100Hz Time=1'÷60' | |
| 10 | Press 1 | Tpulse=120µs Freq.=1÷150Hz Time=1'÷60' | Stamina, Explosive force, Development force |
| 11 | Press 2 | Tpulse=240µs Freq.=1÷150Hz Time=1'÷60' | |
| 12 | Press 3 | Tpulse=300µs Freq.=1÷150Hz Time=1'÷60' | |
| 13 | Kneading 1 | Tpulse=40µs Freq.=1÷150Hz Time=1'÷60' | Tone-up, Modelling, Cellulite, Firm-up |
| 14 | Kneading 2 | Tpulse=160µs Freq.=1÷150Hz Time=1'÷60' | |
| 15 | Kneading 3 | Tpulse=200µs Freq.=1÷150Hz Time=1'÷60' | |
| 16 | Kneading 4 | Tpulse=340µs Freq.=1÷150Hz Time=1'÷60' | |
| 17 | Modulation 1 | Tpulse=80µs Freq.=1÷150Hz Time=1'÷60' | Draining, Stretching |
| 18 | Modulation 2 | Tpulse=140µs Freq.=1÷150Hz Time=1'÷60' | |
| 19 | Burst 1 | Tpulse=180µs Freq.=1÷150Hz Time=1'÷60' | Aerobics, Winding-down |
| 20 | Burst 2 | Tpulse=260µs Freq.=1÷150Hz Time=1'÷60' | |

TESMED PROGRAMS

(programs 21 - 99)

positioning of electrodes

- Programs 21 - 99 are used for treatment with fixed time and stages (not alterable).

| FITNESS | | Prog. | Photo | Time Prog. | n. stages | Days |
|---------|--------------------------------|-----------------|-------------------|------------|-----------|------|
| 1 | Increase arm and pectoral mass | 23 | 1,37,40÷44 | 23' | 3 | 30 |
| 2 | Increase leg mass | 59 | 51÷61,73÷76,78,79 | 21' | 3 | 30 |
| 3 | Boost pectoral mass | 59 | 37 | 21' | 3 | 30 |
| 4 | Boost triceps mass | 58 | 41 | 21' | 3 | 30 |
| 5 | Boost biceps and forearm mass | 58 | 44 | 21' | 3 | 30 |
| 6 | Boost abdomen mass man | 76 | 45÷48 | 29' | 5 | 30 |
| 7 | Boost quadriceps mass | 58 | 54 | 21' | 3 | 30 |
| 8 | Boost upper back mass | 58 | 68 | 21' | 3 | 30 |
| 9 | Stretching triceps | 65 | 64, 66 | 10' | 1 | 30 |
| 10 | Stretching shin muscles | 65 | 13,31,74,75 | 10' | 1 | 30 |
| 11 | Stretching quadriceps | 65 | 54 | 10' | 1 | 30 |
| 12 | Stretching wrist muscles | 65 | 31,34÷36,76÷79 | 10' | 1 | 30 |
| 13 | Lymphatic draining legs | WAIMS® 1 | 10,14,58,76 | 30' | 1 | 15 |
| 14 | Lymphatic draining arms | WAIMS® 2 | 3 | 30' | 1 | 15 |
| 15 | Aerobic | 62 | 1÷79 | 21' | 3 | 15 |
| 16 | Anaerobic | 68 | 1÷79 | 18' | 3 | 15 |
| 17 | Warm-up | 21 | 1÷79 | 10' | 2 | * |
| 18 | Taut muscles | 60 | 1÷79 | 21' | 3 | 30 |
| 19 | Inflamed muscles | 61 | 1÷79 | 21' | 3 | 30 |

* Follow your trainer or doctor's advice
/ No indications provided

| BEAUTY | | | Prog. | Photo | Time Prog. | n. stages | Days |
|--------|-----------------------------------|--------|-------|----------------|------------|-----------|------|
| 20 | Lymphatic draining legs woman | WAIMS® | 66 | 10,11,33 | 23' | 5 | 15 |
| 21 | Lymphatic draining arms woman | WAIMS® | 74 | 3 | 40' | 5 | 15 |
| 22 | Lymphatic draining legs man | WAIMS® | 87 | 58,76 | 40' | 5 | 15 |
| 23 | Lymphatic draining arms man | WAIMS® | 74 | 64 | 40' | 5 | 15 |
| 24 | Soft cellulite arms | | 92 | 20 | 47' | 2 | 30 |
| 25 | Hard cellulite arms | | 89 | 20 | 45' | 2 | 30 |
| 26 | Cellulite abdomen | | 95 | 4÷7 | 47' | 2 | 30 |
| 27 | Cellulite hips | WAIMS® | 93 | 8,9 | 54' | 3 | 30 |
| 28 | Cellulite inside thighs | WAIMS® | 97 | 12÷14 | 49' | 4 | 30 |
| 29 | Soft cellulite outside thighs | WAIMS® | 91 | 31 | 52' | 3 | 30 |
| 30 | Hard cellulite outside thighs | WAIMS® | 98 | 31 | 50' | 3 | 30 |
| 31 | Thick cellulite side thighs | WAIMS® | 90 | 31,33 | 54' | 3 | 30 |
| 32 | Hard cellulite glutei + thighs | WAIMS® | 99 | 29÷31 | 50' | 3 | 30 |
| 33 | Light cellulite glutei | WAIMS® | 94 | 28 | 52' | 3 | 30 |
| 34 | Thick cellulite glutei | WAIMS® | 96 | 28 | 55' | 3 | 30 |
| 35 | Hard cellulite glutei | WAIMS® | 88 | 28 | 50' | 3 | 30 |
| 36 | Muscle definition | | 57 | 1÷79 | 21' | 5 | 30 |
| 37 | Localised slimming arms | | 26 | 2,20 | 24' | 4 | 60 |
| 38 | Localised slimming abdomen | | 76 | 4÷7 | 29' | 5 | 60 |
| 39 | Localised slimming hips | | 77 | 8, 9 | 19' | 5 | 60 |
| 40 | Draining inside thighs woman | | 38 | 12÷14,32 | 20' | 2 | 15 |
| 41 | Localised slimming outside thighs | | 37 | 31 | 26' | 4 | 30 |
| 42 | Localised slimming glutei | | 37 | 27,28 | 26' | 4 | 30 |
| 43 | Localised slimming inside thighs | | 37 | 12÷14,32 | 26' | 4 | 30 |
| 44 | Draining arms | WAIMS® | 22 | 20 | 30' | 4 | 15 |
| 45 | Draining arms and pectorals | | 76 | 1,37,44 | 29' | 5 | 15 |
| 46 | Draining abdomen area | | 76 | 4÷7 | 29' | 5 | 15 |
| 47 | Draining thighs | | 38 | 13 | 20' | 2 | 15 |
| 48 | Draining thighs + glutei | | 82 | 26 | 20' | 1 | 15 |
| 49 | Draining legs | | 38 | 10,11,33,58,76 | 20' | 2 | 15 |
| 50 | Draining heavy legs woman | | 82 | 11 | 20' | 1 | 15 |
| 51 | Draining glutei + thighs | | 82 | 29÷31 | 20' | 1 | 15 |

* Follow your trainer or doctor's advice
/ No indications provided

| BEAUTY | | Prog. | Photo | Time Prog. | n. stages | Days |
|--------|------------------------------|------------------|---------------------------------|------------|-----------|------|
| 52 | Draining inside thighs man | 38 | 54,55,56 | 20' | 2 | 15 |
| 53 | Draining back of legs | 38 | 10,11,14,33 | 20' | 2 | 15 |
| 54 | Modelling arms | 29 | 2,20,40÷44, 64, 66, 67 | 28' | 5 | 30 |
| 55 | Modelling breasts | 78 | 1 | 22' | 4 | 30 |
| 56 | Modelling abdomen area | 76 | 4÷7 | 29' | 5 | 30 |
| 57 | Modeling and dorsal wellness | 48 | 15,16,19, 21÷24,62,68,70 | 28' | 5 | 30 |
| 58 | Modelling thighs + hips | WAIMS® 71 | 8 | 28' | 3 | 30 |
| 59 | Modelling legs | 73 | 10÷14, 32÷36 | 15' | 3 | 30 |
| 60 | Modelling glutei | 77 | 28 | 19' | 5 | 30 |
| 61 | Firm-up | 70 | 1÷79 | 21' | 3 | 30 |
| 62 | Firm-up arms + pectorals | 78 | 1,2,20,37,38, 40÷42,44 | 22' | 4 | 30 |
| 63 | Firm-up breasts | 78 | 1 | 22' | 4 | 30 |
| 64 | Firm-up slack abdomen woman | 76 | 4÷7 | 29' | 5 | 30 |
| 65 | Firm-up abdomen area | 76 | 4÷7 | 29' | 5 | 30 |
| 66 | Firm-up back muscles | 79 | 22,23,62,68,70 | 24' | 4 | 30 |
| 67 | Firm-up thighs | 77 | 13,32 | 19' | 5 | 30 |
| 68 | Firm-up inside thighs | 77 | 12÷14,32 | 19' | 5 | 30 |
| 69 | Firm-up outside thighs | 77 | 31 | 19' | 5 | 30 |
| 70 | Firm-up legs | 79 | 10÷14, 32÷36 | 24' | 4 | 30 |
| 71 | Firm-up hips + glutei | 77 | 26 | 19' | 5 | 30 |
| 72 | Firm-up glutei | 70 | 28 | 21' | 3 | 30 |
| 73 | Anti-stretchmarks breasts | 72 | 1 | 15' | 5 | / |
| 74 | Anti-stretchmarks abdomen | 89 | 4÷7 | 45' | 2 | / |
| 75 | Anti-stretchmarks arms | 72 | 2,20 | 15' | 5 | / |
| 76 | Anti-stretchmarks thighs | 84 | 31 | 15' | 5 | / |
| 77 | Anti-stretchmarks glutei | 84 | 28 | 15' | 5 | / |
| 78 | Anti-stretchmarks hips | 84 | 8,9 | 15' | 5 | / |
| 79 | Tone-up | 69 | 1÷79 | 21' | 3 | 30 |
| 80 | Tone-up arms + pectorals | 69 | 1,2,20,37,38, 40÷42,44 | 21' | 3 | 30 |

* Follow your trainer or doctor's advice
/ No indications provided

| BEAUTY | | Prog. | Photo | Time Prog. | n. stages | Days |
|--------|--|-------|-------------------------------|------------|-----------|------|
| 81 | Tone-up pectorals woman | 80 | 1 | 21' | 3 | 30 |
| 82 | Tone-up pectorals man | 69 | 37,38 | 21' | 3 | 30 |
| 83 | Tone-up arms man | 36 | 40,41,43,44,66 | 20' | 2 | 30 |
| 84 | Tone-up forearms man | 36 | 43 | 20' | 2 | 30 |
| 85 | Combined tone-up upper body woman WAIMS® | 47 | 21÷24 | 26' | 5 | 30 |
| 86 | Tone-up abdomen woman | 75 | 4÷7 | 20' | 4 | 30 |
| 87 | Tone-up abdomen post-natal | 67 | 4÷7 | 21' | 3 | 30 |
| 88 | Combined tone-up abdomen + pectorals man WAIMS® | 56 | 38 | 26' | 5 | 30 |
| 89 | Tone-up upper body muscles | 69 | 70 | 21' | 3 | 30 |
| 90 | Tone-up abdomen man WAIMS® | 56 | 45÷48 | 26' | 5 | 30 |
| 91 | Tone-up abdomen man (gutting) | 75 | 46 | 20' | 4 | 30 |
| 92 | Tone-up flexors man | 69 | 74 | 21' | 3 | 30 |
| 93 | Tone-up back muscles | 64 | 62,68 | 21' | 3 | 30 |
| 94 | Combined tone-up abdomen + flanks man | 75 | 46 | 20' | 4 | 30 |
| 95 | Combined tone-up upper back + shoulders man | 69 | 68 | 21' | 3 | 30 |
| 96 | Combined tone-up shoulders + upper arms man | 69 | 64 | 21' | 3 | 30 |
| 97 | Tone-up flanks WAIMS® | 86 | 9,50 | 26' | 4 | 30 |
| 98 | Tone-up hips post-natal | 80 | 8,9 | 21' | 3 | 30 |
| 99 | Tone-up hips man (love handles) | 76 | 49 | 29' | 5 | 30 |
| 100 | Tone-up adductors man | 69 | 59 | 21' | 3 | 30 |
| 101 | Tone-up inside thighs (anterior) | 80 | 12÷14, 53÷55,59 | 21' | 3 | 30 |
| 102 | Tone-up thighs + hips | 80 | 8 | 21' | 3 | 30 |
| 103 | Tone-up outside thighs | 80 | 31 | 21' | 3 | 30 |
| 104 | Tone-up inside thighs | 80 | 12÷14,32 | 21' | 3 | 30 |
| 105 | Tone-up glutei | 80 | 28 | 21' | 3 | 30 |
| 106 | Tone-up quadriceps man | 69 | 54 | 21' | 3 | 30 |
| 107 | Tone-up calves | 80 | 77 | 21' | 3 | 30 |
| 108 | Tone-up legs man | 69 | 51-61,73÷79 | 21' | 3 | 30 |
| 109 | Tone-up legs woman | 80 | 10÷14,30, 31,33÷36 | 21' | 3 | 30 |

* Follow your trainer or doctor's advice
/ No indications provided

| SPORT | | Prog. | Photo | Time Prog. | n. stages | Days |
|-------|---|-------|-------------------------|------------|-----------|------|
| 110 | Relaxing arms + pectorals | 81 | 1,37,38,40,43,44 | 20' | 1 | * |
| 111 | Relaxing abdominals | 81 | 45÷48 | 20' | 1 | * |
| 112 | Relaxing back muscles | 81 | 19,23,24,64,68 | 20' | 1 | * |
| 113 | Relaxing legs | 81 | 51÷60,73÷79 | 20' | 1 | * |
| 114 | Wind-down legs | 82 | 51÷60,73÷79 | 20' | 1 | * |
| 115 | Sport wind-down legs | 82 | 51÷60,73÷79 | 20' | 1 | * |
| 116 | Strength | 32 | 1÷79 | 23' | 3 | 30 |
| 117 | Explosive strength | 43 | 1÷79 | 23' | 3 | 30 |
| 118 | Explosive strength arms + pectorals | 43 | 37,40,43 | 23' | 3 | 30 |
| 119 | Explosive strength legs | 43 | 51÷60,73÷79 | 23' | 3 | 30 |
| 120 | Max strength arms + pectorals | 51 | 37,40,43 | 23' | 3 | 30 |
| 121 | Max strength abdominals | 42 | 45÷48 | 23' | 4 | 30 |
| 122 | Max strength back muscles | 51 | 62,68 | 23' | 3 | 30 |
| 123 | Max strength legs | 42 | 51÷60,73÷79 | 23' | 4 | 30 |
| 124 | Stamina | 25 | 1÷79 | 26' | 4 | 30 |
| 125 | Stamina arms + pectorals | 44 | 37,40,43 | 23' | 3 | 30 |
| 126 | Stamina abdominals | 53 | 45÷48 | 23' | 3 | 30 |
| 127 | Stamina back muscles | 44 | 62,68 | 23' | 3 | 30 |
| 128 | Stamina legs | 25 | 51÷60,73÷79 | 26' | 4 | 30 |
| 129 | Active recovery | 34 | 1÷79 | 20' | 2 | * |
| 130 | Recovery after races/training | 27 | 1÷79 | 20' | 2 | * |
| 131 | Aerobic stamina | 62 | 1÷79 | 21' | 3 | 30 |
| 132 | Warm-up | 21 | 1÷79 | 10' | 2 | * |
| 133 | Athletic tone-up legs | 69 | 51÷60,73÷79 | 21' | 3 | 30 |
| 134 | Athletic tone-up arms | 69 | 40-44,66,67 | 21' | 3 | 30 |
| 135 | Football training warm-up | 49 | 1÷79 | 10' | 2 | * |
| 136 | Football training max strength | 51 | 1÷79 | 23' | 3 | * |
| 137 | Football training explosive strength | 52 | 1÷79 | 23' | 3 | * |
| 138 | Football training stamina | 53 | 1÷79 | 23' | 3 | * |
| 139 | Football training active recovery | 54 | 1÷79 | 15' | 1 | * |
| 140 | Football training recovery after training | 55 | 1÷79 | 20' | 1 | * |

* Follow your trainer or doctor's advice
/ No indications provided

| SPORT | | | Prog. | Photo | Time Prog. | n. stages | Days |
|-------|---|--------------------------|-------|--|------------|-----------|------|
| 141 | Football training lymphatic draining | WAIMS[®] | 87 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 142 | Ski training warm-up | | 40 | 1÷79 | 10' | 2 | * |
| 143 | Ski training max strength | | 42 | 1÷79 | 23' | 4 | * |
| 144 | Ski training explosive strength | | 43 | 1÷79 | 23' | 3 | * |
| 145 | Ski training stamina | | 44 | 1÷79 | 23' | 3 | * |
| 146 | Ski training active recovery | | 45 | 1÷79 | 15' | 5 | * |
| 147 | Ski training recovery after training | | 46 | 1÷79 | 20' | 1 | * |
| 148 | Ski training lymphatic draining | WAIMS[®] | 87 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 149 | Cycling training warm-up | | 21 | 1÷79 | 10' | 2 | * |
| 150 | Cycling training max strength | | 23 | 1÷79 | 23' | 3 | * |
| 151 | Cycling training explosive strength | | 24 | 1÷79 | 23' | 3 | * |
| 152 | Cycling training stamina | | 25 | 1÷79 | 26' | 4 | * |
| 153 | Cycling training active recovery | | 63 | 1÷79 | 15' | 3 | * |
| 154 | Cycling training recovery after training | | 27 | 1÷79 | 20' | 2 | * |
| 155 | Cycling training lymphatic draining | WAIMS[®] | 87 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 156 | Swimming training warm-up | | 40 | 1÷79 | 10' | 2 | * |
| 157 | Swimming training max strength | | 42 | 1÷79 | 23' | 4 | * |
| 158 | Swimming training explosive strength | | 43 | 1÷79 | 23' | 3 | * |
| 159 | Swimming training stamina | | 25 | 1÷79 | 26' | 4 | * |
| 160 | Swimming training active recovery | | 83 | 1÷79 | 15' | 3 | * |
| 161 | Swimming training recovery after training | | 46 | 1÷79 | 20' | 1 | * |
| 162 | Swimming training lymphatic draining | WAIMS[®] | 74 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |

* Follow your trainer or doctor's advice
/ No indications provided

| SPORT | | Prog. | Photo | Time Prog. | n. stages | Days |
|-------|--|-------|---|------------|-----------|------|
| 163 | Running training warm-up | 21 | 1÷79 | 10' | 2 | * |
| 164 | Running training max strength | 23 | 1÷79 | 23' | 3 | * |
| 165 | Running training explosive strength | 24 | 1÷79 | 23' | 3 | * |
| 166 | Running training stamina | 33 | 1÷79 | 23' | 3 | * |
| 167 | Running training active recovery | 34 | 1÷79 | 20' | 2 | * |
| 168 | Running training recovery after training | 27 | 1÷79 | 20' | 2 | * |
| 169 | Running training lymphatic draining WAIMS® | 87 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 170 | Tennis training warm-up | 49 | 1÷79 | 10' | 2 | * |
| 171 | Tennis training max strength | 42 | 1÷79 | 23' | 4 | * |
| 172 | Tennis training explosive strength | 52 | 1÷79 | 23' | 3 | * |
| 173 | Tennis training stamina | 33 | 1÷79 | 23' | 3 | * |
| 174 | Tennis training active recovery | 83 | 1÷79 | 15' | 3 | * |
| 175 | Tennis training recovery after training | 27 | 1÷79 | 20' | 2 | * |
| 176 | Tennis training lymphatic draining WAIMS® | 74 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 177 | Volleyball training warm-up | 21 | 1÷79 | 10' | 2 | * |
| 178 | Volleyball training max strength | 51 | 1÷79 | 23' | 3 | * |
| 179 | Volleyball training explosive strength | 43 | 1÷79 | 23' | 3 | * |
| 180 | Volleyball training stamina | 25 | 1÷79 | 26' | 4 | * |
| 181 | Volleyball training active recovery | 83 | 1÷79 | 15' | 3 | * |
| 182 | Volleyball training recovery after training | 46 | 1÷79 | 20' | 1 | * |
| 183 | Volleyball training lymphatic draining WAIMS® | 74 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 184 | Basketball training warm-up | 21 | 1÷79 | 10' | 2 | * |
| 185 | Basketball training max strength | 51 | 1÷79 | 23' | 3 | * |
| 186 | Basketball training explosive strength | 52 | 1÷79 | 23' | 3 | * |
| 187 | Basketball training stamina | 25 | 1÷79 | 26' | 4 | * |

* Follow your trainer or doctor's advice
/ No indications provided

| SPORT | | Prog. | Photo | Time Prog. | n. stages | Days |
|-------|--|-------|--|------------|-----------|------|
| 188 | Basketball training active recovery | 45 | 1÷79 | 15' | 5 | * |
| 189 | Basketball training recovery after training | 27 | 1÷79 | 20' | 2 | * |
| 190 | Basketball training lymphatic draining WAIMS® | 87 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |

* Follow your trainer or doctor's advice
/ No indications provided

| MESSAGES TENS (soft intensity) | | Prog. | Photo | Time Prog. | n. stages | Days |
|-----------------------------------|--|-------|-------|------------|-----------|------|
| 191 | Relaxing massage neck + shoulderse | 28 | 15 | 15' | 1 | 15 |
| | | 48 | | 28' | 5 | 30 |
| | | 39 | | 25' | 5 | 30 |
| 192 | Massage upper + lower back woman WAIMS® | 48 | 22 | 28' | 5 | 30 |
| | | 85 | | 40' | 5 | 15 |
| | | 39 | | 25' | 5 | 30 |
| 193 | Massage upper + lower back man WAIMS® | 48 | 70 | 28' | 5 | 30 |
| | | 85 | | 40' | 5 | 15 |
| | | 39 | | 25' | 5 | 30 |
| 194 | Massage legs woman | 31 | 11 | 22' | 4 | 15 |
| 195 | Relaxing massage calves (cramp) | 81 | 77 | 20' | 1 | 15 |
| 196 | Anti-stress massage | 41 | 1÷79 | 20' | 1 | 15 |
| 197 | Wind-down massage | 50 | 11 | 10' | 2 | 15 |
| 198 | Energetic massage | 31 | 1÷79 | 22' | 4 | 15 |
| 199 | Relaxing massage | 30 | 1÷79 | 15' | 3 | 15 |
| 200 | Toning massage | 35 | 1÷79 | 21' | 3 | 15 |
| 201 | Massage TENS all body | 48 | 1÷79 | 28' | 5 | 30 |

* Follow your trainer or doctor's advice
/ No indications provided

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| Prog. | Onde | Paramètres | Indication Traitements |
|-------|---------------|---|---|
| 1 | Waims-System® | 1 → 2 → 3 → 4 (pause 1) 1 → 2 → 3 → 4 | Lymphodrainage, Relaxant, Récupération active, Défatiguant |
| 2 | Waims-System® | 1 → 2 → 3 → 4 (pause 2) 1 → 2 → 3 → 4 | |
| 3 | Waims-System® | 1 → 2 → 3 → 4 (pause 3) 1 → 2 → 3 → 4 | |
| 4 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 1) 1 → 2 → 3 → 4 → 3 → 2 → 1 | |
| 5 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 2) 1 → 2 → 3 → 4 → 3 → 2 → 1 | |
| 6 | Tapping 1 | Tpulse=100µs Freq.=1÷100Hz Time=1'÷60' | Décontractant, Réactivant, Cellulite |
| 7 | Tapping 2 | Tpulse=200µs Freq.=1÷100Hz Time=1'÷60' | |
| 8 | Tapping 3 | Tpulse=300µs Freq.=1÷100Hz Time=1'÷60' | |
| 9 | Tapping 4 | Tpulse=400µs Freq.=1÷100Hz Time=1'÷60' | Force résistante, Force explosive, Développement de la force |
| 10 | Press 1 | Tpulse=120µs Freq.=1÷150Hz Time=1'÷60' | |
| 11 | Press 2 | Tpulse=240µs Freq.=1÷150Hz Time=1'÷60' | |
| 12 | Press 3 | Tpulse=300µs Freq.=1÷150Hz Time=1'÷60' | Tonification, Modelage, Cellulite, Raffermissment |
| 13 | Kneading 1 | Tpulse=40µs Freq.=1÷150Hz Time=1'÷60' | |
| 14 | Kneading 2 | Tpulse=160µs Freq.=1÷150Hz Time=1'÷60' | |
| 15 | Kneading 3 | Tpulse=200µs Freq.=1÷150Hz Time=1'÷60' | |
| 16 | Kneading 4 | Tpulse=340µs Freq.=1÷150Hz Time=1'÷60' | Drainant, Stretching |
| 17 | Modulation 1 | Tpulse=80µs Freq.=1÷150Hz Time=1'÷60' | |
| 18 | Modulation 2 | Tpulse=140µs Freq.=1÷150Hz Time=1'÷60' | Aérobique, Hypertonie |
| 19 | Burst 1 | Tpulse=180µs Freq.=1÷150Hz Time=1'÷60' | |
| 20 | Burst 2 | Tpulse=260µs Freq.=1÷150Hz Time=1'÷60' | |

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| Prog. | Welle | Parameter | Angabe Behandlung |
|-------|---------------|---|---|
| 1 | Waims-System® | 1 → 2 → 3 → 4 (pause 1) 1 → 2 → 3 → 4 | Lymphknotendrainage Entspannung Aktive Kräfterückgewinnung Erholung |
| 2 | Waims-System® | 1 → 2 → 3 → 4 (pause 2) 1 → 2 → 3 → 4 | |
| 3 | Waims-System® | 1 → 2 → 3 → 4 (pause 3) 1 → 2 → 3 → 4 | |
| 4 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 1) 1 → 2 → 3 → 4 → 3 → 2 → 1 | |
| 5 | Waims-System® | 1 → 2 → 3 → 4 → 3 → 2 → 1 (pause 2) 1 → 2 → 3 → 4 → 3 → 2 → 1 | |
| 6 | Tapping 1 | Tpulse=100µs Freq.=1÷100Hz Time=1'÷60' | Lösen von Verspannungen Reaktivierung Zellulitis |
| 7 | Tapping 2 | Tpulse=200µs Freq.=1÷100Hz Time=1'÷60' | |
| 8 | Tapping 3 | Tpulse=300µs Freq.=1÷100Hz Time=1'÷60' | |
| 9 | Tapping 4 | Tpulse=400µs Freq.=1÷100Hz Time=1'÷60' | |
| 10 | Press 1 | Tpulse=120µs Freq.=1÷150Hz Time=1'÷60' | Widerstandskraft Explosive Kraft Kraftentwicklung |
| 11 | Press 2 | Tpulse=240µs Freq.=1÷150Hz Time=1'÷60' | |
| 12 | Press 3 | Tpulse=300µs Freq.=1÷150Hz Time=1'÷60' | |
| 13 | Kneading 1 | Tpulse=40µs Freq.=1÷150Hz Time=1'÷60' | Stärkung Modellierung Zellulitis, Kräftigung |
| 14 | Kneading 2 | Tpulse=160µs Freq.=1÷150Hz Time=1'÷60' | |
| 15 | Kneading 3 | Tpulse=200µs Freq.=1÷150Hz Time=1'÷60' | |
| 16 | Kneading 4 | Tpulse=340µs Freq.=1÷150Hz Time=1'÷60' | |
| 17 | Modulation 1 | Tpulse=80µs Freq.=1÷150Hz Time=1'÷60' | Drainage Stretching |
| 18 | Modulation 2 | Tpulse=140µs Freq.=1÷150Hz Time=1'÷60' | |
| 19 | Burst 1 | Tpulse=180µs Freq.=1÷150Hz Time=1'÷60' | Aerobik Hypertonie |
| 20 | Burst 2 | Tpulse=260µs Freq.=1÷150Hz Time=1'÷60' | |

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| SPORT | | Prog. | Foto | Dauer Prog. | Anzahl Phasen | Tage |
|-------|---|-------|--|-------------|---------------|------|
| 163 | Training Laufen Warm Up | 21 | 1÷79 | 10' | 2 | * |
| 164 | Training Laufen maximale Kraft | 23 | 1÷79 | 23' | 3 | * |
| 165 | Training Laufen explosive Kraft | 24 | 1÷79 | 23' | 3 | * |
| 166 | Training Laufen Widerstandskraft | 33 | 1÷79 | 23' | 3 | * |
| 167 | Training Laufen Aktive Erholung | 34 | 1÷79 | 20' | 2 | * |
| 168 | Training Laufen Erholung nach dem Training | 27 | 1÷79 | 20' | 2 | * |
| 169 | Training Laufen Lymphknotendrainage WAIMS® | 87 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 170 | Training Tennis Warm Up | 49 | 1÷79 | 10' | 2 | * |
| 171 | Training Tennis maximale Kraft | 42 | 1÷79 | 23' | 4 | * |
| 172 | Training Tennis explosive Kraft | 52 | 1÷79 | 23' | 3 | * |
| 173 | Training Tennis Widerstandskraft | 33 | 1÷79 | 23' | 3 | * |
| 174 | Training Tennis Aktive Erholung | 83 | 1÷79 | 15' | 3 | * |
| 175 | Training Tennis Erholung nach dem Training | 27 | 1÷79 | 20' | 2 | * |
| 176 | Training Tennis Lymphknotendrainage WAIMS® | 74 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 177 | Training Handball Warm Up | 21 | 1÷79 | 10' | 2 | * |
| 178 | Training Handball maximale Kraft | 51 | 1÷79 | 23' | 3 | * |
| 179 | Training Handball explosive Kraft | 43 | 1÷79 | 23' | 3 | * |
| 180 | Training Handball Widerstandskraft | 25 | 1÷79 | 26' | 4 | * |
| 181 | Training Handball Aktive Erholung | 83 | 1÷79 | 15' | 3 | * |
| 182 | Training Handball Erholung nach dem Training | 46 | 1÷79 | 20' | 1 | * |
| 183 | Training Handball Lymphknotendrainage WAIMS® | 74 | 3,10÷14,22, 25,33,34,51, 53,54,56,57, 58,70,75,76 | 40' | 5 | * |
| 184 | Training Basketball Warm Up | 21 | 1÷79 | 10' | 2 | * |
| 185 | Training Basketball maximale Kraft | 51 | 1÷79 | 23' | 3 | * |
| 186 | Training Basketball explosive Kraft | 52 | 1÷79 | 23' | 3 | * |
| 187 | Training Basketball Widerstandskraft | 25 | 1÷79 | 26' | 4 | * |

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Feldi S.r.l. - Via Iseo 46, 25030 - Erbusco (BS) - Italy